


CATALOGUE

2020-2021



Thank you for taking the time to go through our new catalogue for 2020-2021!

Since 1996, MVS In Motion has been your solid business partner as manufacturer and master distributor for a wide range of physiotherapy, rehabilitation and sports & fitness products. Those products reach the end-users through you, our valued distributors and dealers.

2020 has been a remarkable year due to the health crisis caused by the Covid-19 pandemic. We all have learned to adapt to this new reality, which will continue to influence the business environment we all operate in. You can rely on us to remain your solid business partner.

This catalogue will again bring some new products to our range: we are proud to introduce the revolutionary Wheeleo, a rolling cane, and are adding new Oximeter Junior, Shotgun Thermometer, Blood Pressure meter and Neurological Hammers to our Measurement range, just to list a few. Page 5 shows all the new products. We continue to offer Private Label products as well, contact us for inquiries.

An important milestone in our sector, the EU's new Medical Device Regulation (MDR) is delayed for one year and will now come into effect May 2021. We are in the process of ensuring that all our Medical Devices will be compliant by May 2021, in fact some are already meeting the new requirements today. In relation to MDR, we are also working hard to obtain ISO13485-certification in 2021.

We have had our challenges so far but we did not sit back, and worked hard to improve our business and serve you better. We will continue to do so in 2021!

Greetings on behalf of the entire team.

Sincerely,
Sandra Vrancaert
CEO/Owner



Sandra Vrancaert
CEO / Owner
sandra.vrancaert@mvs-in-motion.com



Tom Lambrechts
Business Development / Sales
tom.lambrechts@mvs-in-motion.com



Manuela Stobbelaar
Logistics
logistics@mvs-in-motion.com



Monique Breuer
Customer Service
customer.service@mvs-in-motion.com



Patrick Steenackers
Warehousing



Tom Benoy
Sales
tom.benoy@mvs-in-motion.com



Griet Van de Velde
Marketing
marketing@mvs-in-motion.com



Patrick Verlinden
Administration
info@mvs-in-motion.com



Andy Huari
Warehousing

CONTENT

01 - Resistive Exercise

6



03 - Upper & Lower Body Exercise

33



05 - Fitness

46



07 - Thermal Treatment

58



02 - Hand Therapy

23



4 - Wellness

41



06 - Strength

55



08 - Measurement

62



OUR OWN BRANDS



MoVeS products are aimed towards physiotherapy and rehabilitation, predominantly resistive exercise. Resistive exercise products have various levels of resistance, which also makes them perfect for fitness. MoVeS products are generally CE certified.



Our Mambo Max products are the go-to-products for wellness, fitness and strength training at an affordable price tag. Combine various products for a complete work-out, without leaving your home.

ONLINE

www.moves-you.com
www.moves-band.com

The virtual home of our MoVeS and MoVeS FIT product line. This B2C platform connects our distributors directly to the end users, with product information, a shop locator and numerous free downloadable exercises. The website has been upgraded and is mobile friendly.

www.mambomax.com

Completely revitalised! This is the virtual home of our Mambo Max product line. This B2C platform connects our products directly to the end users, with product information and numerous free downloadable exercises.

www.mvs-in-motion.com

MVS In Motion's B2B site, with all available products and digital catalogue. Ask for your login to download all high resolution pictures, product information, logos, videos and many more marketing tools.

Social Media













Here you will find our latest news and useful information. Like or follow our pages for a higher brand awareness to create more demand. Our pages are continuously monitored to answer your questions as soon as possible. We are active on LinkedIn and YouTube. Other social media channels can be added later.



By making exercises available on our websites, not only do we reduce our ecological footprint by cutting down on the number of unused prints, it also allows us to keep costs low and offer you better value-for-money.

LEGEND

In order to keep all the information in this catalogue as clear and easy accessible as possible, we use a system of icons that are shown next to the product where applicable.



NEW PRODUCTS

Band & Tubing Scissors



22

Wheeleo



33

Spiky Foam Roller



41

Knobble® II



45

Electric Compressor Pump



48

PearlPack
Knee/Elbow Wrap



59

Shotgun Thermometer



62

Fingertip Pulse Oximeter
Junior



63

Blood Pressure Monitor



63

Dejerine Hammer



73

Troemner Hammer



73



LEADER IN TECHNOLOGY

What is Snap-Stop®?

Snap-Stop® makes exercising with latex bands much safer. This unique feature will avoid sudden snapping, even if there's a cut or puncture. When working out with a resistive exercise band, the band has to endure extremely high tensions. A small cut can have a big consequence...

MoVeS and MoVeS FIT Bands are unique in their field. Thanks to our extensive research and development we have achieved an innovative production process: when used properly, the latex band will not suddenly snap, even when exposed to high forces. This does not apply to a very old band, or band that has been exposed outdoors for some time, as well as an abused band. Using a fresh and clean sheet of band is recommended.

Snap-Stop® for safer resistive exercising.



See for yourself

One-Stop-Shop for all your bands and tubes

MVS In Motion is your one-stop-shop for all resistive exercise products. With MoVeS we offer both original colour and FIT colour sequence. In addition, we offer a wide range of Latex-Free resistive exercise bands in original colours. Our tubes and loops are also available in the original and FIT colour sequence.



The latex bands, tubes and (super)loops are made of high quality latex from huge rubber tree plantations in Malaysia. All products have an extraordinary elasticity and linearity. The Snap-Stop® feature on our latex bands, dipped production process for our tubes and the multi-layered structure of our (super)loops makes them your ideal partner for your physiotherapy trainings or intensive work-outs. The different resistances allow you to use MoVeS for rehabilitation as well as for fitness.

MoVeS Latex-Free is made from TPE, to the same high quality standards as the regular MoVeS bands. Never before were the characteristics from Latex-Free bands so close to those from their latex counterpart.

MoVeS Band: 8 resistances (Tan - Gold) & 5 lengths (1,5m - 2,5m - 5,5m - 22,5m - 45,5m)

MoVeS FIT Band: 6 resistances (Peach - Grey) & 4 lengths (1,5m - 2,5m - 5,5m - 45,5m)

MoVeS Latex-Free Band: 5 resistances (Yellow - Black) & 3 lengths (1,5m - 5,5m - 45,5m)



Original Colours

Tan	Extra Light
Yellow	Light
Red	Medium
Green	Heavy
Blue	Extra Heavy
Black	Special Heavy
Silver	Super Heavy
Gold	Ultra Heavy

Same Resistance



FIT Colours

Peach	Light
Orange	Medium
Lime Green	Heavy
Blueberry	Extra Heavy
Plum	Special Heavy
Grey	Super Heavy





- Well-known colours used by physiotherapists all over the world
- Packaging with two tone layout brings more variety
- MoVeS brings premium quality to an affordable price range.




- "FIT" - colours are more trendy and youthful than original colours, reaching out to a wider customer base such as athletes and home users
- Sell MoVeS FIT next to MoVeS to widen customer base = EXTRA SALES



Resistance - MoVeS Band vs. Market Leader





Light
Medium
Heavy
Extra Heavy
Special Heavy



Market Leader:
Some non-linear behaviour
+ green and blue are almost
identical

MoVeS Band:
Best combination of
resistance / linearity

Other alternative:
Very high resistance



MoVeS Bands can be stretched up to 600% elongation! MoVeS Bands are therefore your perfect partner for a rehabilitation programme or a full body work-out. Even after an effective training with amazing results, MoVeS Bands allow you to use the same band again. When you pre-stretch your band, you will be able to achieve higher forces without having to worry about the 'elastic stop', the point where the band feels like a regular string. Most competitor's bands are not to be stretched more than 3 times their resting length, while a MoVeS Band achieves easily the double, maintaining its amazing elasticity.

These specifications apply to all MoVeS resistances. The high-quality material enables you to have a much wider array of resistances and possibilities compared to the competitors. At the start of your exercise, the band already gives some resistance, which allows you to achieve much better results than when the band is at its resting length.

Note:

- 50% elongation = 1,5 x lay-flat length
- 100% elongation = 2 x lay-flat length
- 150% elongation = 2,5 x lay-flat length



MoVeS Band



MoVeS Bands are premium quality, low powder (in fact, less than 2mg/dm², so we're officially powder-free) latex resistive exercise bands. The unique Snap-Stop® fabrication process makes this product extremely tear resistant. Using the MoVeS Bands for resistive exercises provides positive as well as negative forces on the muscles and improves the range of motion, the strength and the cooperation of the muscle groups. The 8 different resistance levels are determined by the thickness of the band. MoVeS bands are available in 5 lengths and come with a full colour manual. The band has a width of +/- 127mm.

Dispenser 25 bands of 1,5m in a nicely designed OTC dispenser box. Each band comes in a separate polybag.

- 1,5m 10 pre-cut bands of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header and barcode.
- 2,5m 10 pre-cut pieces of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header and barcode.
- 5,5m One small roll of our premium band.
- 22,5m One medium roll of our premium band.
- 45,5m One large roll of our premium band.



COLOUR	TOTAL KG / ELONGATION %				
	50%	100%	150%	200%	250%
Yellow	1,0	1,3	1,5	1,9	2,3
Red	1,5	2,0	2,5	3,2	3,7
Green	1,6	2,5	3,0	3,5	4,1
Blue	2,2	3,0	3,6	4,4	5,2
Black	2,5	3,4	4,0	4,9	6,1
Silver	4,1	5,8	7,2	7,9	9,5
Gold	5,8	7,9	10,0	12,0	14,0

Polybag (10-pack)
1,5m / 2,5m



Box with one Roll
5,5m / 22,5m / 45,5m



Dispenser 1,5m



PLASTIC
FREE






COLOUR	RESISTANCE	MOVES BAND					
		1,5m Dispenser 25pcs	1,5m 10-pack	2,5m 10-pack	5,5m Box Small Roll	22,5m Box Medium Roll	45,5m Box Large Roll
Tan	Extra Light	n/a	n/a	n/a	01-100501	n/a	01-104501
Yellow	Light	01-100102	01-101502	01-100202	01-100502	01-102202	01-104502
Red	Medium	01-100103	01-101503	01-100203	01-100503	01-102203	01-104503
Green	Heavy	01-100104	01-101504	01-100204	01-100504	01-102204	01-104504
Blue	Extra Heavy	01-100105	01-101505	01-100205	01-100505	01-102205	01-104505
Black	Special Heavy	01-100106	01-101506	01-100206	01-100506	01-102206	01-104506
Silver	Super Heavy	n/a	01-101507	01-100207	01-100507	01-102207	n/a
Gold	Ultra Heavy	n/a	n/a	n/a	01-100508	01-102208	n/a
Set of 5 resistances (yellow to black)		n/a	01-101511 (5-pack)	01-100211 (5-pack)	n/a		n/a

Resistance - MoVeS Latex-Free Band vs. Market Leader

Research and development is high on our list for MoVeS products. As an example, these continuous efforts have brought us the Snap-Stop feature for the regular MoVeS Band.

Now, we have improved our Latex-Free Bands. The MoVeS Latex-Free Band has been optimised by using an innovative mix of raw material. Following customer demand, the grip has been improved. Resistance levels are similar to the 'old' latex-free band.



MoVeS Band:
Best combination of resistance / linearity

Market Leader:
Non-linear behaviour

Other alternative:
Very low resistance

We are confident that MoVeS Latex-Free Bands are one of the best latex-free resistive exercise bands on the market. The elasticity and linearity are simply amazing. Finally people who don't support latex can exercise with the same high quality specifications. MoVeS Latex-Free bands start from high quality thermoplastic elastomer (TPE), so it's non-toxic and completely PVC and phthalate free.

The linearity from how the strength per elongation evolves is one of the most important specifications as it offers the possibility to apply a very extended range of exercises without risking an injury due to strange unexpected product behaviour. Latex-free alternatives will reach their 'elastic stop' more quickly; this is the point where the band feels like a regular string.

We notice a growing demand and trend for latex free products. For instance, more and more hospitals and centers aim to be completely latex free to avoid possible allergic reactions for their patients or staff. After extensive research and development, we were able to launch our MoVeS Latex-Free Bands. A product which we are extremely proud of.



MoVeS Latex-Free Band



Thanks to our technological leadership, we achieve a Latex-Free exercise band with characteristics that were never before so close to its latex counterpart. With an exceptional 'resistance vs. elongation' linearity, we are confident our MoVeS Latex-Free Band is the best product to fulfil your needs.

The 5 different resistance levels are determined by the thickness of the band. These bands are available in 3 lengths and come with a full colour manual. The band has a width of +/- 100mm.

- 1,5m 10 pre-cut bands of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header.
- 5,5m One small roll of our latex free band.
- 45,5m One large roll of our latex free band.

TPE
phthalate-free

COLOUR	TOTAL KG / ELONGATION %					
	50%	100%	150%	200%	250%	300%
Yellow	1,3	1,7	1,9	2,3	2,8	3,2
Red	1,4	1,9	2,2	2,6	3,0	3,5
Green	2,0	2,5	3,0	3,4	4,0	4,5
Blue	2,3	3,0	3,4	3,9	4,7	5,5
Black	2,4	3,4	3,9	4,3	5,1	6,1



Box with one Roll
5,5m



COLOUR	RESISTANCE	MOVES LATEX-FREE BAND		
		1,5m 10-pack	5,5m Box Small Roll	45,5m Box Large Roll
Yellow	Light	01-121502	01-120502	01-124502
Red	Medium	01-121503	01-120503	01-124503
Green	Heavy	01-121504	01-120504	01-124504
Blue	Extra Heavy	01-121505	01-120505	01-124505
Black	Special Heavy	01-121506	01-120506	01-124506
Set of 5 resistances (yellow to black)		01-121511 (5-pack)	n/a	n/a



Polybag (10-pack)
1,5m



Box with one Roll
45,5m



Latex-Free



MoVeS Tube



The MoVeS Tube is a high quality tube created by an intensive dipping process. The dipping process guarantees a multi-layered tube, which ensures a strong and safe training product. Extruded tubes are very susceptible to damage and can easily snap, this is not the case for our multi-layered dipped tubes. In general, dipped tubes are even more durable than bands, especially when used on a coarse surface.

MoVeS Tubes are available as a ready to use product - 120cm length with handles straight out of a polybag - and in rolls - 7,5m or 30m packed in a nicely designed box.

The table below shows the exceptional strength to elongation linearity from our tubes.

COLOUR	TOTAL KG / ELONGATION %					
	50%	100%	150%	200%	250%	300%
Yellow	1,0	1,5	1,9	2,3	2,9	3,6
Red	1,4	2,0	2,5	2,8	3,4	4,0
Green	2,4	3,4	4,0	4,9	5,5	6,5
Blue	2,9	3,9	4,7	5,6	5,9	7,4
Black	3,4	4,4	5,4	6,2	7,1	8,1
Silver	4,0	5,4	6,5	7,7	8,1	9,5

Combine a roll with our Safety Handles for superior grip (p22).



MoVeS Tubes with very useful for and muscle training. These high quality and very resistant tubes are created to train all muscles in your body and improve body flexibility. The lighter resistances are perfect for physiotherapy or a light work-out, while the highest resistances will even convince the more trained athlete.

Handles are physiotherapy. These high resistant



DIPPED
not extruded

Box with Rolls
7,5m / 30m



Box of 30m contains maximal 3 separate pieces.
Minimal length of 1 piece is 4m.



Polybag
120cm with Handles



COLOUR	RESISTANCE	MOVES TUBE		
		120cm with Handles	7,5m roll	30m roll
Yellow	Light	01-200102	01-200702	01-203002
Red	Medium	01-200103	01-200703	01-203003
Green	Heavy	01-200104	01-200704	01-203004
Blue	Extra Heavy	01-200105	01-200705	01-203005
Black	Special Heavy	01-200106	01-200706	01-203006
Silver	Super Heavy	01-200107	01-200707	01-203007




MoVeS Wide Loop



We offer wide loops (7,5cm) next to our regular loops (2,5cm) with the same resistances. Some users may prefer the wider loops for spreading the forces more evenly on their body parts, while others prefer the thicker narrow loops which won't roll-up while being stretched.

As your supplier for bands, tubes and loops, these wide loops fit perfectly in our range of resistive products to satisfy all users. Wide loops are also multi-layered and available in 6 resistances: from Extra Light (Tan) to Special Heavy (Black).



Multi-Layered

MoVeS Wide Loops:

High quality, cost effective and easy to use!

Lay-flat-length: 30cm

Width: 7,5cm

One unit contains 10 individually packed pcs of the same colour/resistance



COLOUR	RESISTANCE	MOVES WIDE LOOP	UNIT
Tan	Extra Light	01-100071	10 pcs
Yellow	Light	01-100072	10 pcs
Red	Medium	01-100073	10 pcs
Green	Heavy	01-100074	10 pcs
Blue	Extra Heavy	01-100075	10 pcs
Black	Special Heavy	01-100076	10 pcs



PLASTIC
FREE



MoVeS Loop



MoVeS Loops are extremely portable and unbelievably versatile products. They are very effective instruments for stretching and strengthening all muscle groups. The MoVeS Loop Set really is a "Gym on-the-go".

Our loops are also thicker loops

multi-layered which results in a more durable product. These won't roll-up during exercising, avoiding the uncomfortable feeling from thinner alternatives.



Your
Gym on-the-go

They are available in 6 resistances: from Extra Light (Tan) to Special Heavy (Black).

The loop set combines our most popular colours in a compact mesh bag. This affordable set is very useful to bring along wherever you go.

MoVeS Loops:

High quality, cost effective and easy to use!

Lay-flat-length: 30cm

Width: 2,5cm

One polybag contains 10 loose pcs of the same colour/resistance

COLOUR	RESISTANCE	MOVES LOOP	UNIT
Tan	Extra Light	01-100001	10 pcs
Yellow	Light	01-100002	10 pcs
Red	Medium	01-100003	10 pcs
Green	Heavy	01-100004	10 pcs
Blue	Extra Heavy	01-100005	10 pcs
Black	Special Heavy	01-100006	10 pcs
Loop Set of 4 resistances (yellow to blue) in mesh bag		01-100011	Set / 4 pcs

COLOUR	TOTAL KG / ELONGATION %					
	50%	100%	150%	200%	250%	300%
Tan	2.2	3.2	3.8	4.6	5.4	6.4
Yellow	2.8	4.0	4.8	5.8	6.8	8.0
Red	3.4	4.8	5.8	6.8	8.0	9.2
Green	4.4	6.4	7.8	9.2	10.8	12.4
Blue	5.0	7.2	9.0	10.8	12.6	14.4
Black	6.1	8.9	11.0	13.0	15.0	17.3



Multi-Layered




The production of MoVeS latex products start deep in the humid forests of Malaysia, where the best and strongest trees provide the natural latex rubber that will be processed from a milky fluid to your preferred resistive exercise band. When a tree is 5 - 6 years old, harvesting can begin: incisions are being made orthogonally in the bark of the tree, just deep enough to tap the latex vessels without harming the tree's growth. The fluid is collected in small buckets, a process known as rubber tapping.

Did you know that rubber trees combat greenhouse gases and thereby slowing down global warming? Additionally, latex harvesting helps prevent deforestation.

natural
100%
biodegradable

"This is why we are confident when using the words 'premium quality' when referring to MoVeS."

Only the best quality latex is being selected for MoVeS Bands. Production starts under the constant supervision of a team of scientists, ensuring an immaculate product quality. Latex sheets are produced, powdered, washed, dried and cut to the exact width. The amount of powder used is below the 2mg/dm² limit, which makes the MoVeS Bands officially powder-free. Each batch of latex sheeting that is produced, is subject to the meticulous inspection of our scientists by a series of laboratory experiments: colour-coding, dryness, tensile strength, tearing resistance, snapping resistance, ... During the rolling process, a 100% visual inspection is being conducted on both sides of the band. Indeed, before leaving the factory, every MoVeS product has been completely analysed, checked and examined to ensure premium product quality.



MoVeS FIT Band



MoVeS FIT, for athletes and home use.

MoVeS FIT Bands are low powder (in fact, less than 2mg/dm², so we're officially powder-free) latex resistive exercise band with the same premium quality as the original MoVeS Bands. The difference is purely the alternative, more youthful colour sequence.

The band comes in 6 different resistance levels which are determined by the thickness of the band, is available in 4 different lengths and comes with a full colour manual.

COLOUR	TOTAL KG / ELONGATION %					Compares to STANDARD
	50%	100%	150%	200%	250%	
Peach	1,0	1,3	1,5	1,9	2,3	Yellow
Orange	1,5	2,0	2,5	3,2	3,7	Red
Lime Green	1,6	2,5	3,0	3,5	4,1	Green
Blueberry	2,2	3,0	3,6	4,4	5,2	Blue
Plum	2,5	3,4	4,0	4,9	6,1	Black
Grey	4,1	5,8	7,2	7,9	9,5	Silver

Our high quality 1,5m and 2,5m MoVeS FIT Bands are packed in retail packs. 10 pre-cut pieces of the same resistance or 5 pieces of different resistance that come in separate polybags with header for easy dispensing. Full-colour instruction manual is included. These 1,5m and 2,5m bands are easy to use and cutting is not required. The available colours are peach, orange, lime green, blueberry, plum and grey.



		MOVES FIT BAND				
COLOUR	RESISTANCE	1,5m 10-pack	2,5m 10-pack	5,5m in Box Small Roll	22,5m in Box Medium roll	45,5m in Box Large Roll
Peach	Light	01-111502	01-110202	01-110502	n/a	01-114502
Orange	Medium	01-111503	01-110203	01-110503	n/a	01-114503
Lime Green	Heavy	01-111504	01-110204	01-110504	n/a	01-114504
Blueberry	Extra Heavy	01-111505	01-110205	01-110505	n/a	01-114505
Plum	Special Heavy	01-111506	01-110206	01-110506	n/a	01-114506
Grey	Super Heavy	01-111507	01-110207	01-110507	01-112207	n/a
Set of 5 resistances (peach to plum) (5-pack)		01-111511 (5-pack)	01-110211 (5-pack)	n/a	n/a	n/a



MoVeS FIT Tube



The MoVeS FIT Tube is a high quality tube created by an intensive dipping process. The dipping process guarantees a multi-layered tube, which ensures a strong and safe training product. Extruded tubes are very susceptible to damage and can easily snap, this is not the case for our multi-layered dipped tubes. In general, dipped tubes are even more durable than bands, especially when used on a coarse surface.

MoVeS FIT Tubes are available as a ready to use product in 30m rolls, packed in a great dispenser box. Safety Handles offer an unlimited amount of exercises and superior grip.

MoVeS FIT Tubes are available in 6 different colours. The table below shows the resistance and to which standard colour it relates.

COLOUR	TOTAL KG / ELONGATION %					Compares to STANDARD
	50%	100%	150%	200%	250%	
Peach	1,0	1,5	1,9	2,3	2,9	Yellow
Orange	1,4	2,0	2,5	2,8	3,4	Red
Lime Green	2,4	3,4	4,0	4,9	5,5	Green
Blueberry	2,9	3,9	4,7	5,6	5,9	Blue
Plum	3,4	4,4	5,4	6,2	7,1	Black
Grey	4,0	5,4	6,5	7,7	8,1	Silver

*DIPPED
not extruded*





Box of 30m contains maximal 3 separate pieces.
Minimal length of 1 piece is 4m.



COLOUR	RESISTANCE	MOVES FIT TUBE 30m roll
Peach	Light	01-213002
Orange	Medium	01-213003
Lime Green	Heavy	01-213004
Blueberry	Extra Heavy	01-213005
Plum	Special Heavy	01-213006
Grey	Super Heavy	01-213007



Combine a roll with our
Safety Handles for superior grip (p22).



MoVeS FIT Loop



*Your
Gym on-the-go*

Our high quality and cost effective loops are also available in FIT colours. MoVeS FIT Loops are easy to use and very versatile. You can stretch and strengthen all muscle groups.

Our loops are multi-layered, this results in a more durable product. Our thicker loops won't roll-up during exercising, avoiding the uncomfortable feeling from thinner alternatives.

Our loop set combines our most popular colours in a compact mesh bag. This affordable set is very useful to bring along wherever you go.

MoVeS FIT Loops:
High quality, cost effective and easy to use!
Lay-flat-length: 30cm
Width: 2,5cm
One Polybag contains 10 loose pcs

Multi-Layered

COLOUR	RESISTANCE	MOVES FIT LOOP	UNIT
Peach	Light	01-110002	10 pcs
Orange	Medium	01-110003	10 pcs
Lime Green	Heavy	01-110004	10 pcs
Blueberry	Extra Heavy	01-110005	10 pcs
Loop Set of 4 resistances (peach to blueberry) in mesh bag		01-110011	Set / 4 pcs



MoVeS FIT Superloop



COLOUR	SIZE IN MM	TOTAL KG / ELONGATION %				
		100%	200%	300%	400%	500%
Peach	1040 x 15 x 3,0	7,3	10,0	14,4	19,1	29,2
Orange	1040 x 22 x 4,0	13,0	20,7	27,1	35,0	66,3
Lime Green	1040 x 30 x 4,4	22,5	33,0	44,0	56,2	97,0
Blueberry	1040 x 40 x 4,4	25,6	40,0	53,0	76,0	132,0
Plum	1040 x 53 x 4,4	38,0	56,0	79,0	105,0	167,0

With the MoVeS FIT Superloop you can easily stretch and strengthen different muscle groups. It can be used for rehabilitation, fitness and strength training. The lay-flat-length from the MoVeS FIT Superloop is 104cm long (circumference is 208cm) and has an incredible elasticity.

The MoVeS FIT Superloops are multi-layered, this results in a more durable product. They are available in 5 different resistances: Peach, Orange, Lime Green, Blueberry and Plum.

COLOUR	RESISTANCE	MOVES FIT SUPERLOOP	UNIT
Peach	Light	01-110402	1 pc
Orange	Medium	01-110403	1 pc
Lime Green	Heavy	01-110404	1 pc
Blueberry	Extra Heavy	01-110405	1 pc
Plum	Special Heavy	01-110406	1 pc





MoVeS Bar



The MoVeS Bar is the ideal product to treat and resolve golfer's and tennis elbow.
This effective product is developed to focus on hand, arm and shoulder problems, this product allows flexion and extension exercises, adduction and abduction exercises and even soft tissue mobilisation, far more than the twisting and bending you would imagine.
The MoVeS Bar is available in 6 colours / resistances: extra light, light, medium, heavy, extra heavy and special heavy.

Size: 31cm x 4,5cm
Packaging: polybag with full-colour header
Material: dried natural rubber

1 size Ø; 6 Colours / Resistances

COLOUR	RESISTANCE	MOVES BAR
Tan	Extra Light	01-600101
Yellow	Light	01-600102
Red	Medium	01-600103
Green	Heavy	01-600104
Blue	Extra Heavy	01-600105
Black	Special Heavy	01-600106

Ideal for Golfer's & Tennis Elbow treatment



MoVeS Ring Tube



The MoVeS Cuff-Ring Tube is a versatile exercise tool to tone and strengthen all the muscles of the lower body. With velcro closure that fits most sizes.



For intensive training of isolated muscles. Efficient for all muscles of upper- and lower body. The centre of the MoVeS 8-Ring Tube offers lots of possibilities for short and long pull exercises.

Ideal for abdomen, bottom, hip and leg exercises. With the MoVeS O-Ring Tube you get an intensive training of each muscular group.



COLOUR	RESISTANCE	O-RING TUBE LOOP	FIGURE 8-RING TUBE LOOP	CUFF-RING TUBE LOOP*
Yellow	Light	01-200202	01-200302	01-200402
Red	Medium	01-200203	01-200303	01-200403
Green	Heavy	01-200204	01-200304	01-200404
Blue	Extra Heavy	01-200205	01-200305	01-200405
Black	Special Heavy	01-200206	01-200306	01-200406

DIPPED! not extruded



* This product cannot be sold in Germany.



Multi-Grip Exerciser

The Multi-Grip™ Exerciser is designed to provide unlimited exercise options to users who want to do various high impact performance exercises. Each exerciser contains both large and small loops that are easy to grip. Use the small loops as an anchor or for finger exercises. Use the larger loops around the hands and feet to perform high impact upper and lower body exercises. Easily change resistance by moving 1 loop up or down the exerciser, or using a different resistance progressively exercise with 8 different color-coded resistance levels.

The Multi-Grip™ Exerciser is made of an unscented cloth that can be washed and dried. It does not contain latex or a powder coating.



COLOUR	RESISTANCE	MULTI-GRIP EXERCISER 30 Yard Roll (27,4m)
Tan	XX-Light	01-703001
Yellow	X-Light	01-703002
Red	Light	01-703003
Green	Medium	01-703004
Blue	Heavy	01-703005
Black	X-Heavy	01-703006
Silver	XX- Heavy	01-703007
Gold	XXX-Heavy	01-703008



COLOUR	TOTAL KG / ELONGATION %	
	100%	200%
Tan	0,9	2,3
Yellow	1,4	3,4
Red	1,8	4,5
Green	2,3	5,7
Blue	2,7	6,8
Black	3,2	7,9
Silver	4,5	11,3
Gold	5,9	14,7



MoVeS Ankleciser




The MoVeS Ankleciser is a great ankle exerciser, complying to all the rigorous quality standards MVS In Motion demands of its MoVeS products. It's as easy in use as it's effective. The feet slip comfortably into the latex free elastic band which consists of two stitched loops. Available in 4 different resistances, this progressive exercise system enables the user to choose the desired resistance level to perform the exercises. Available in a set of 4 resistances, or in a 5-pack of the same resistance. Packed in a polybag with full-colour header and exercise chart.



COLOUR	RESISTANCE	MOVES ANKLECISER	UNIT
Yellow	Light	01-500102	5 pcs
Red	Medium	01-500103	5 pcs
Green	Heavy	01-500104	5 pcs
Blue	Extra Heavy	01-500105	5 pcs
Set of 4 resistances		01-500111	4 pcs



MoVeS Shoulder Rope Pulley - Metal Bracket




This MoVeS Shoulder Rope Pulley comes with a pulley, wooden handles, nylon cord, stopper system, and metal bracket.

Reference	Description	Unit
01-400102	MoVeS Shoulder Rope Pulley Metal Bracket Wooden Handles Stoppers	1 pc



MoVeS Shoulder Tube Pulley



Just like the MoVeS Shoulder Rope Pulley, the MoVeS Shoulder Tube Pulley helps increase range of shoulder motion. It can be used for prevention and rehabilitation of shoulder injuries.

The MoVeS Shoulder Tube Pulley has a light, medium or heavy tube. Each pack contains the selected resistance of tube, soft handles, pulley and door anchor.

COLOUR	RESISTANCE	SHOULDER TUBE PULLEY	UNIT
Yellow	Light	01-400202	1 pc
Red	Medium	01-400203	1 pc
Green	Heavy	01-400204	1 pc



MoVeS Knee Rope Pulley



The MoVeS Knee Rope Pulley helps increase range of knee motion. This easy-to-use product can be mounted in a doorway to help you recover from or prevent knee injuries.

While sitting on a chair, put your foot into the ankle strap. A rope starts from the ankle strap, goes through a pulley and ends at a wooden handle designed for a 2-handed grip. Pull your foot carefully upwards by moving the handle towards your abdomen which will generate a rotational movement at your knee.

The rope length can easily be adjusted at the handle while a print on the ankle strap helps you to secure your foot in the right position. The MoVeS Knee Rope Pulley comes in a full colour cardboard packaging with a pulley, wooden 2-handed handle, ankle strap, nylon cord and door anchor.

Reference	Description	Unit
01-400301	MoVeS Knee Rope Pulley Door Anchor 2-handed Wooden Handle Ankle Strap	1 pc




MoVeS Shoulder Rope Pulley - Door Anchor




The MoVeS Shoulder Rope Pulley is an economical, easy-to-use device that helps increase range of shoulder motion. It can be used for prevention and rehabilitation of shoulder injuries. This easy to use shoulder rope pulley can be mounted in a doorway.

NOTCH HANDLES m cb ex CE




This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and notch handles. The unique notch handles allow the user to easily adjust the length of the rope through the handle while the revolutionary stopper system is a great aid for consistent and progressive exercising.

ECO-PACK po ex



The ECO PACK is packed in a polybag, with notch handles but without stoppers on the nylon cord.

WOODEN HANDLES bl ex CE



This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and wooden handles.

HAND SUPPORT bl ex CE



This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and detachable hand support. This hand support allows patients without gripping power to use the shoulder pulley.

Reference	Description	Unit
01-400101	MoVeS Shoulder Rope Pulley Door Anchor Notch Handles Stoppers	1 pc
01-400104	MoVeS Shoulder Rope Pulley Door Anchor ECO PACK (notch handles, polybag, no stoppers)	1 pc
01-400105	MoVeS Shoulder Rope Pulley Door Anchor Wooden Handles Stoppers	1 pc
01-400103	MoVeS Shoulder Rope Pulley Door Anchor Hand Support Stoppers	1 pc

MoVeS Accessories



Reference	Description	Unit
01-300001	MoVeS Wall Mount (Plug & Play)	1 pc
01-300101	MoVeS Safety Handles	pair
01-300102	MoVeS Soft Handles	pair
01-300103	MoVeS Door Anchor	1 pc
01-300104	MoVeS Universal Assist	1 pc
01-300106	MoVeS Band & Tubing scissors NEW	1 pc



MoVeS Flex-Ion



The spring loaded buttons of the MoVeS Flex-Ion isolate each digit, so the strong fingers cannot compensate for the weaker ones. This results in a better, faster and easier to control rehab process. MoVeS Flex-Ion is available in 6 colour coded units with gradual resistance. This gradual resistance helps the therapist measure progress. The ergonomic base adds comfort and stability. Comes with a full colour and multilingual instruction manual.



COLOUR	RESISTANCE (1 digit all digits)	LEVEL	REFERENCE	UNIT
Tan	0,35kg 1,1kg	0	02-000101	1 pc
Yellow	0,7kg 2,3kg	1	02-000102	1 pc
Red	1,4kg 4,5kg	2	02-000103	1 pc
Green	2,3kg 7,3kg	3	02-000104	1 pc
Blue	3,2kg 10,4kg	4	02-000105	1 pc
Black	4,1kg 14,1kg	5	02-000106	1 pc
MoVeS Flex-Ion Set 6 Resistances 1 Display			02-000111	set
MoVeS Flex-Ion Empty Display Metal			02-000121	1 pc





Varigrip Therapy



Resistance adjustable per finger

The VariGrip Therapy has the ability to adjust tension in each individual finger, the resistance from each single finger knob can easily be adjusted. Resistance can be increased incrementally on the patient's need, which minimises overcompensation and strain. The VariGrip allows the patient to progress steadily between resistance levels.

5 different resistance levels are available. Each Varigrip's highest resistance ends where the next Varigrip resistance starts. Any other hand exerciser increases from Light (1,4kg) to Medium (2,3kg). That's almost a 70% increase gap, which can cause an overload on the muscles, resulting in injuries. The seamless progress from one Varigrip to the other, supported by the padded finger and palm area, guarantees a safe recovery and better feel for those with sensitivity issues.

Comparison Chart: VariGrip vs Digi-Flex® / Flex-Ion (kg per digit)

VariGrip	X-Light		Light		Medium		Heavy		X-Heavy	
Digi-Flex® / Flex-Ion	0,7	1,4	1,4	2,3	2,3	3,2	3,2	4,1	4,1	n/a

COLOUR	RESISTANCE (range per digit)	LEVEL	REFERENCE	UNIT
X-Light	0,7kg - 1,4kg	1	02-100102	1 pc
Light	1,4kg - 2,3kg	2	02-100103	1 pc
Medium	2,3kg - 3,2kg	3	02-100104	1 pc
Heavy	3,2kg - 4,1kg	4	02-100105	1 pc
X-Heavy	4,1kg - 5,4kg	5	02-100106	1 pc
VariGrip Set 5 Resistances 1 Display			02-100107	set
VariGrip Empty Display Plastic			02-100108	1 pc



Varigrip Sport



02-100304



02-100303

Resistance adjustable per finger


VariGrip Sport has the patented ability to adjust the tension level in each individual finger, allowing the user to create a custom programme to build and improve endurance, strength and dexterity more effectively for individual fingers, as well as the entire hand, wrist and forearm. VariGrip Sport is the only exerciser of its kind that includes a callus builder, ideal for climbers, gymnasts, gamers, musicians and athletes of all levels.

Reference	Description	Resistance (kg)	Unit
02-100303	VariGrip Sport	2,3 - 4,1	1 pc
02-100304	VariGrip Sport Pro	4,1 - 5,4	1 pc



Constant Force X-Tend

Unlike elastic-based exercisers that increase in resistance the further the unit is extended, Constant Force X-Tend uses retractable-based resistance, which keeps a constant level of resistance from start to finish. Ergonomic in design, the Constant Force X-Tend works the Intrinsic and Extrinsic Extensors of the hand and can isolate MP/PIP/DIP joint motion. It also allows each finger to exercise independently with smooth, consistent resistance, thus improving range of motion, endurance, strength and dexterity more effectively.



- Exercise each finger or thumb independently with smooth, consistent resistance to the intrinsic and extrinsic extensor muscles throughout their full range of motion.
- Improve range of motion, strength, dexterity, endurance, circulation for each finger as well as the hand, wrist and forearm.
- Ergonomic design fits comfortably in hand, and includes ability to isolate pip & dip joints.



COLOUR	RESISTANCE (per digit)	REFERENCE	UNIT
Light	0,11 kg	02-100403	1 pc
Medium	0,18 kg	02-100404	1 pc
Heavy	0,27 kg	02-100405	1 pc

not elastic = constant force



Varigrip Uno

Soft touch thumb button

The only thumb/individual finger exerciser with adjustable resistance for both flexion and extension. An ergonomic design with stabilisation points offers extra support.

Gradually adjust resistance and better track progression based on the patients needs. Benefits are improved strength, dexterity, endurance and range of motion.




Both flexion and extension

Light Spring	0,7kg - 2,3kg
Heavy Spring	2,3kg - 3,6kg



Adjustable resistance allows a smooth transition between resistance levels



Reference	Description	Unit
02-100204	VariGrip Uno	1 pc





Digi-Flex

Digi-Flex® develops isolated finger strength, flexibility and coordination. Suitable for rehabilitation after hand injury. The Digi-Flex® comes with an exercise leaflet.

COLOUR	RESISTANCE (1 digit - all digits)	REFERENCE	UNIT
Yellow	0,7kg - 2,3kg	02-010102	1 pc
Red	1,4kg - 4,5kg	02-010103	1 pc
Green	2,3kg - 7,3kg	02-010104	1 pc
Blue	3,2kg - 10,4kg	02-010105	1 pc
Black	4,1kg - 14,1kg	02-010106	1 pc
Digi-Flex Set 5 Resistances 1 Display		02-010111	set
Digi-Flex Empty Display Metal		02-010121	1 pc



Digi-Extend

The Digi-Extend® is a unique therapy device that can isolate the extension of all five fingers and each of the individual finger joints. The product strengthens the intrinsic and extrinsic muscle groups to maximise extensor tendon gliding, joint motion, coordination and muscle function in the hand and forearm.

The Digi-Extend® focuses on the extensor mechanism of one or all fingers as it targets the area of the hand that needs rehabilitation. Patients can easily change the resistance levels and vary their hand positions for specific exercises to help maximise the exercise goals and prevent injury. The Digi-Extend® can also be used for finger adduction, abduction, blocking and flexion exercises. It is compact in design, high-tech, easy and comfortable to use as well as effective in helping patients reach their therapy goals. Four colour coded band resistance levels allow each muscle or joint to work at a comfortable level and build as increased strength is achieved.

Reference	Description	Unit
02-010201	Digi-Extend Bands (Tan to Green) included	1 pc
02-010211	Digi-Extend Set 4 units 1 Display	set
02-010202	Digi-Extend Band Set 25 pcs 5x Tan, 5x Yellow, 5x Red, 5x Green, 5x Blue	set
02-010221	Digi-Extend Empty Display Metal	1 pc

Bands included:
- 5x Tan
- 4x Yellow
- 3x Red
- 2x Green



Prohands Series



VIA™ is an entry level hand exerciser designed with young athletes in mind. It has a round, soft look, wide finger pistons for smoother function & improved finger performance.

The VIA™ is available in two tensions.

RESISTANCE	VIA	UNIT
Light	02-060102	1 pc
Medium	02-060103	1 pc



PRO™ provides serious athletes with the most advanced hand exerciser in the world. The PRO™ features a larger profile, superior steel springs that resist wear, an integrated ergonomic palm bar for comfort.

The PRO™ is available in three tensions.

RESISTANCE	RESISTANCE (per digit)	PRO	UNIT
Medium	3,5kg	02-060303	1 pc
Heavy	4,5kg	02-060304	1 pc
Extra Heavy	5,5kg	02-060305	1 pc



Handmaster Plus



The Handmaster Plus™ strengthens the 9 muscles that close the hand as well as the 9 muscles that open the hand. The result is maximum strength, balance and blood flow to the hand, wrist and elbow.



3-pack set
= 3 balls for the price of 2!



Soft | Purple Ball - Yellow Cord | light exercise
Medium | Red Ball - Black Cord | medium exercise
Firm | Orange Ball - Grey Cord | strength, speed training



Reference	Description	Unit
02-050102	Handmaster Plus Soft Purple Ball - Yellow Cord	1 pc
02-050103	Handmaster Plus Medium Red Ball - Black Cord	1 pc
02-050104	Handmaster Plus Firm Orange Ball - Grey Cord	1 pc
02-050111	Handmaster Plus Set of 3 (one of each level)	Set





Power-Web Junior



The Power-Web® Junior is 18cm in diameter and weighs 225g. Just half the size of the Power-Web® Senior, it has all the same features as the larger product and delivers many of the same benefits when used in a hand and wrist exercise programme.



COLOUR	RESISTANCE	POWER-WEB JUNIOR
Yellow	Light	02-020302
Red	Medium	02-020303
Green	Heavy	02-020304



Power-Web Senior



The Power-Web® Senior is 36cm in diameter and weighs 450g. Its unique patented design allows a variety of exercises to be performed comfortably and effectively. Constructed of high quality rubber with special agents added for durability and strength, the Power-Web® will last for years. Now also available as latex-free!



COLOUR	RESISTANCE	POWER-WEB SENIOR	POWER-WEB LATEX-FREE SENIOR
Tan	Ultra Light	02-020201	02-022201
Yellow	Light	02-020202	02-022202
Red	Medium	02-020203	02-022203
Green	Heavy	02-020204	02-022204
Blue	Super Heavy	02-020205	02-022205
Black	Ultra Heavy	02-020206	02-022206



Power-Web Combo



The #1 hand exercise web in the world!

Each Power-Web comes with a 39 exercise poster



COLOUR / RESISTANCE		POWER-WEB COMBO
Tan/Ultra Light	Red/Medium	02-020101
Yellow/Light	Green/Heavy	02-020102
Red/Medium	Blue/Super Heavy	02-020103
Green/Heavy	Black/Ultra Heavy	02-020104



The Power-Web® Combo features a unique advantage of two resistance levels in one, offering greater flexibility and value. It is identical in size and weight to the Power-Web® Senior with the same high quality construction and materials.



Power-Web Flex-Grip



The creators of the best-selling Power-Web®, are proud to present the original Power-Web® Flex-Grip®. The Power-Web® Flex-Grip® is a unique exerciser that fits in your palm and provides the ability to work dexterity, extension and strengthening in your hands, fingers, and thumbs. Available in 3 different resistances, the Power-Web® Flex-Grip® has finger and thumb holes and has two sides: one side featuring a ball to squeeze and the other featuring buttons and a thumb pad to exercise fingers separately and increase dexterity and strength in your fingers and thumb. Latex-free and easy to use!



COLOUR	RESISTANCE	POWER-WEB FLEX-GRIP
Yellow	Light	02-021102
Red	Medium	02-021103
Green	Heavy	02-021104





MoVeS Squeeze Balls & Eggs



The MoVeS Squeeze Ball is developed to provide a progressive sequence of 5 colours. It strengthens hands, fingers, forearms and will work the muscles separately or in tandem. The MoVeS Squeeze Ball can be used for rehab as well as for sports training or stress relief.

The MoVeS Squeeze Egg is ergonomically shaped to provide highly effective exercise for the rehabilitation and strengthening of weakened or injured fingers, hands and wrists. Because it is small and easy to use, you can exercise anywhere, anytime.

Available colours/resistances: yellow, red, green, blue and black.



The MoVeS Squeeze Balls & Eggs can help to:


- Rehabilitate post-operative/post-traumatic hands and arms.
- Restore function in hands, arms and shoulders following stroke.
- Reduce pain and stiffness in arthritic fingers, hands and wrists.
- Relax tense, sore muscles and reduce stress.
- Strengthen aging hands and arms.

COLOUR	RESISTANCE	LEVEL	SQUEEZE BALL	SQUEEZE EGG
Yellow	Extra Soft	1	02-030102	02-030202
Red	Soft	2	02-030103	02-030203
Green	Medium	3	02-030104	02-030204
Blue	Firm	4	02-030105	02-030205
Black	Extra Firm	5	02-030106	02-030206



Why Anti-Microbial Agents?

Anti-microbial agents kill bacterium without additives and resist growth of odour-causing microbes. MoVeS Comfort Putty contains natural anti-microbial without addition of toxic chemicals. This additive ensures the MoVeS Comfort Putty material retains the clean, non-greasy feel from day one. Cheaper alternatives on the market do not contain these agents. As a small amount of bacteria and other micro-organisms will always be present on the patient's hands and because putty should be stored in room temperature conditions, anti-microbial agents are a must to avoid the putty becoming the ideal breeding ground for bacteria and fungi. Anti-microbial agents extend the product's lifespan and improves its cleanliness.



MoVeS Comfort Putty



The MoVeS Comfort Putty does everything a great putty should do: it strengthens weakened grasp through exercise and is used for rehabilitative exercise or stress relief. It is odourless, does not stain and will not leave any residue on the patient's hands. MoVeS Comfort Putty is non-toxic, greaseless, hypo-allergenic and gluten and latex-free. MoVeS Comfort Putty is available in 6 colours/resistances and comes pre-packed.

With anti-microbial agents



COLOUR	RESISTANCE	LEVEL	COMFORT PUTTY	
			85g	454g
White	Special Soft	-1	02-090101	02-090401
Tan	Extra Soft	0	02-090102	02-090402
Yellow	Soft	1	02-090103	02-090403
Red	Medium	2	02-090104	02-090404
Green	Firm	3	02-090105	02-090405
Blue	Extra Firm	4	02-090106	02-090406



Theraflex Putty



Theraflex™ Putty strengthens weakened grasp through exercise and is used for rehabilitative exercise or stress relief. It is odourless, does not stain and will not leave any residue on the patient's hands. Theraflex™ Putty is non-toxic, greaseless, hypo-allergenic and latex-free. Theraflex™ Putty is available in 5 colours/resistances.

With anti-microbial agents



COLOUR	RESISTANCE	LEVEL	THERAFLEX PUTTY					
			57g	85g	113g	170g	454g	2,3kg
Tan	Extra Soft	0	02-040001	02-040101	02-040201	02-040301	02-040401	02-040501
Yellow	Soft	1	02-040002	02-040102	02-040202	02-040302	02-040402	02-040502
Red	Medium	2	02-040003	02-040103	02-040203	02-040303	02-040403	02-040503
Green	Firm	3	02-040004	02-040104	02-040204	02-040304	02-040404	02-040504
Blue	Extra Firm	4	02-040005	02-040105	02-040205	02-040305	02-040405	02-040505



MoVeS Wrist Exerciser



Adjustable resistance



The MoVeS Wrist Exerciser is an adjustable wrist and forearm trainer which allows the user to exercise with progressive resistance. Its ergonomically designed forearm support keeps the wrist in the right place. Together with the soft handgrip, adjustable to 4 positions, it ensures maximum comfort while rehabilitating your wrist and forearm.

Reference	Description	Unit
02-070201	MoVeS Wrist Exerciser	1 pc




Adjustable Hand Grip



The Mambo Max Adjustable Hand Grip comes in two different versions. Whether it is the 5 to 20 kg light version (orange) or the 10 to 40 kg heavy version (grey), you can always adjust the required tension simply by turning the control knob.

Adjustable resistance



Reference	Description	Unit
02-070102	Mambo Max Adjustable Hand Grip - Light	1 pc
02-070104	Mambo Max Adjustable Hand Grip - Heavy	1 pc



Foam Hand Grip



The Mambo Max Foam Hand Grip assures a firm training grip. The foam handles provide a firm grip for the entire workout session. Take it with you to the office or just use it at home. Wherever you want to use it, the Mambo Max Foam Hand Grip assures a tight grasp.

Gripping power = 20kg

Reference	Description	Unit
02-080101	Mambo Max Foam Hand Grip	1 pair



Wheeleo



Reduced mobility presents 3 challenges:

- Keep the person moving;
- Reduce the risk of falling;
- Loss of autonomy.

Wheeleo, the rolling cane, is designed to offer a solution to these 3 challenges and offers additional benefits:

- Ease of handling, no special training required;
- **Walking** goes a lot **faster** than with a quadripod cane;
- **Safer walking** than with a quadripod cane, as there is constant ground contact;
- **One-handed use** while still maintaining constant ground contact, the other hand can be used to move an object, to open doors, ...
- **Minimal space required** and perfect for indoor use where a rollator is too big.
- **Reinforced stability**: the wheels offer a security system. A lateral push will result in a lateral movement, not a tilting/fall.



Faster and safer walking



Wheeleo is a unique quadripod rolling cane used to improve mobility. It is designed with the aim of offering help without stigmatising the disability of its user.

It will make walking trips smoother and safer which will allow the user to maintain or **improve autonomy**.

The Wheeleo is used with one hand. Thanks to its 4 wheels, you do not have to lift it, which improves stability and **reduces the risk of falling**.

The Wheeleo is super easy to use. Ultra-intuitive, its user will always position it in the right place to guarantee its stability and balance.

An important feature of the Wheeleo is that it **does not require brakes** to prevent its user from falling. Thanks to its rolling design, the Wheeleo will always be placed in the right position for its user.

When losing balance, the user will automatically push the Wheeleo in the right position that is comfortable for him, without lifting the Wheeleo.

Adding brakes would be the same as replacing the Wheeleo with a regular cane because the brakes would block the wheels. As a consequence the Wheeleo would tilt and its user would possibly fall.




Adjustable height

Rolling Cane Wheeleo

± 77,5 - 101,5cm
38 x 42cm
120kg

Reference	Description	Unit
03-060101	Wheeleo NEW	1 pc



MoVeS Oxycycle 1 - Active Magnetic



For ACTIVE training only

The OxyCycle 1 Pedal Exerciser is used for upper and lower body exercises. The **adjustable resistance** makes it a great product to progressively strengthen your arms and legs. The OxyCycle 1 Pedal Exerciser has a light-weight compact design, adjustable resistance and a multi-function LCD display readout of scan, time, distance, calories and RPM. A non-skid mat and interchangeable handgrip pedals are included.



- 50x46x38 cm
- 9 kg

Adjustable
RESISTANCE

The OxyCycle 1 comes with strap pedals and handgrip pedals.

Reference	Description	Unit
03-010101	MoVeS OxyCycle 1 Pedal Exerciser	1 pc



MoVeS Oxycycle Pedals



All pedals are interchangeable and fit OxyCycle 1, 2 & 3.

Strap Pedal

Handgrip Pedal



Full-Foot Pedal



Half-Foot Pedal

Reference	Description	Unit
03-020101	MoVeS Strap Pedals	pair
03-020102	MoVeS Half-Foot Pedals	pair
03-020103	MoVeS Full-Foot Pedals	pair
03-020201	MoVeS Handgrip Pedals	pair



MoVeS Oxycycle 2 - Passive



For PASSIVE training only

The MoVeS OxyCycle 2 Pedal Exerciser is a great tool for personalised, low-stress workouts at home. The **adjustable speed level** makes it a great product to progressively strengthen your arms and legs. This pedal exerciser provides low impact exercise to increase strength and flexibility, is versatile for both upper and lower body and it improves the blood circulation. It offers up to **60 Watt** of motorized assistance. It has enough power to maintain a constant rotational speed, whereas cheaper models only carry a 30 Watt motor. A safety stop will avoid injury in case of cramp or muscle spasm. The MoVeS OxyCycle 2 has a light-weight compact design, adjustable speed and a multi-function LCD display readout of scan, time, count, RPM-level and estimated calories burned. A non-skid mat and interchangeable handgrip pedals are included.



- 50x46x38 cm
- 8,5 kg
- 50 Hz
- 25 - 65 rpm

NEW
Comes with UK or EU plug!

Adjustable
ASSISTANCE

The OxyCycle 2 comes with half-foot pedals and handgrip pedals.

Reference	Description	Unit
03-010102	MoVeS OxyCycle 2 Pedal Exerciser EU-plug	1 pc
03-010102 UK	MoVeS OxyCycle 2 Pedal Exerciser UK-plug	1 pc



MoVeS Oxycycle 3 - Active & Passive



For ACTIVE and PASSIVE training

The MoVeS OxyCycle 3 Pedal Exerciser is a great tool for personalised, low-stress workouts at home. The **adjustable speed level and resistance** makes it a great product to progressively strengthen your arms and legs. This pedal exerciser provides low impact exercise to increase strength and flexibility, is versatile for both upper and lower body and it improves the blood circulation. It has a reliable **80 Watt** motor to offer enough power to maintain a constant rotational speed, whereas cheaper models only carry a 30 Watt motor. A safety stop will avoid injury in case of cramp or muscle spasm. The MoVeS OxyCycle 3 has a light-weight compact design, adjustable speed & resistance, a 15-minute timer switch and a multi-function LCD display readout of scan, time, count, total count and calories burned. A non-skid mat and interchangeable handgrip pedals are included.



- 50x46x38 cm
- 8,5 kg
- 50 Hz
- 25 - 65 rpm

NEW
Comes with UK or EU plug!

Adjustable
RESISTANCE &
ASSISTANCE

The OxyCycle 3 comes with full-foot pedals and handgrip pedals.

Reference	Description	Unit
03-010103	MoVeS OxyCycle 3 Pedal Exerciser EU-plug	1 pc
03-010103 UK	MoVeS OxyCycle 3 Pedal Exerciser UK-plug	1 pc





Curamotion

EXER 1
96 x 45 x 45 (cm) / 11,5 kg
Wall mount included



EXER 2
146 x 45 x 54 (cm) / 13,5 kg
Wall mount included



Improvement of blood circulation, muscle strength, joint range of motion and shoulder exercise. The CuraMotion Exer Series has an adjustable resistance and a multi-functional LCD-display readout scan which indicates rotation count, rotations per minute, time exercised and calories consumed. All are easy adjustable with a quick release knob.



EXER 3
37 x 38 x 70 (cm) / 9 kg
With foot stands



EXER 4
96 x 45 x 45 (cm) / 13 kg
Wall mount included

Reference	Description	Unit
03-030101	CuraMotion Exer 1	1 pc
03-030102	CuraMotion Exer 2	1 pc
03-030103	CuraMotion Exer 3	1 pc
03-030104	CuraMotion Exer 4	1 pc



MoVeS Pedal Exerciser



MoVeS Pedal Exercisers are ideal for restoring impaired muscle strength, improving blood circulation and give a better coordination to the legs or arms. It can be used on a table top as an arm exerciser or on the floor to give the legs a workout. A convenient resistance control knob allows a broad range of pedal resistance for a progressive exercise programme.

MoVeS Foldable Pedal Exerciser € €

- Foldable for compact shipping and storage, easy setup
- Pedal straps adjustable to two positions
- LCD Display (time, rotations per minute, total rotations and calories)
- In colour box



Foldable & LCD display

- 49x36x27 cm
- 2,3 kg



03-010203

MoVeS Standard Pedal Exerciser | Premium Chrome € €

- Pre-assembled
- Pedal straps adjustable to two positions
- Quality chrome finish
- In colour box or on specific demand in bulk with separate colour box



03-010201

MoVeS Standard Pedal Exerciser K/D € €


- Knock-down for compact shipping and storage (only 6 screws to assemble or disassemble)
- Pedal straps adjustable to two positions
- Budget solution
- In colour box



Budget Pack

03-010202

Reference	Description	Unit
03-010201	MoVeS Standard Pedal Exerciser Premium Chrome	1 pc
03-010202	MoVeS Standard Pedal Exerciser K/D White	1 pc
03-010203	MoVeS Foldable Pedal Exerciser + LCD Display	1 pc





Fit Stretch



The Mambo Max Fit Stretch is a fantastic, simple and portable tool to stretch your lower leg muscles. It is perfect for warm-ups and cool-downs. The use of a Mambo Max Fit Stretch helps increase your ankle, shin and calf flexibility. It is made of durable, moulded plastic and includes a non-skid foot pad.

Reference	Description	Unit
03-050101	Mambo Max Fit Stretch Single Leg	1 pc



03-010301: Saehan Cycle - Home

Saehan Cycle



The Saehan Cycles are professional pedal exercisers aimed at rehabilitation. Pedalling improves blood circulation and will increase strength and flexibility after an injury or for elderly. Both the Home and Clinical Saehan Cycles include attachments to fit the pedal exerciser to a chair or even a wheelchair. The pedal exerciser's length can be perfectly adjusted to the user when his feet are fitted in the pedals with foot straps. A convenient turning knob makes it possible to adjust the resistance (0 - 13 kg) which increases or decreases the effort, resulting in a progressive exercise programme.

The Clinical version from the Saehan Cycle has an even higher level of adjustability. The distance between pedal and crank can be adjusted and the resistance knob indicates the level selected. Together with the convenient grips to attach the Saehan Cycle to a (wheel)chair and different extra accessories, this product will be an essential item in any professional's practice room.

Because the Saehan Cycle - Clinical is a product aimed at professionals, we offer various accessories to increase its flexibility and user's safety. Be aware these do only apply to the clinical cycle, not to the home cycle.

Reference	Description	Unit
03-010311	Clinical Cycle Accessory Long front leg	1 pc
03-010312	Clinical Cycle Accessory Extension arm (right = left)	1 pc
03-010313	Clinical Cycle Accessory Hand strap - right	1 pc
03-010314	Clinical Cycle Accessory Hand strap - left	1 pc
03-010315	Clinical Cycle Accessory Long (straight) leg brace set	1 pc
03-010316	Clinical Cycle Accessory Bent leg brace	1 pc

Reference	Description	Unit
03-010301	Saehan Cycle - Home	1 pc
03-010302	Saehan Cycle - Clinical	1 pc



MoVeS Wooden Incline Board



The incline board increases extension in the lower and upper back while stretching the calves. It can be used:

- to improve the general posture,
- to stretch the calf muscles of athletes, football players and sportspersons in general,
- to stretch the upper-leg muscles (Ischias),
- to increase extension in the lower and upper back.

MoVeS Wooden Incline Board - Ladder



The MoVeS Wooden Incline Board - Ladder is made out of natural birch wood and is very resistant. It is adjustable in 4 different angles (10°, 15°, 20° and 25°). After using the MoVeS Wooden Incline board, it can easily be folded to achieve its compact form when not in use

MoVeS Wooden Incline Board - Fix



58cm x 30cm x 16cm
150kg

The new MoVeS Wooden Incline Board - Fix is perfect for who's looking for an easy and budget-friendly Incline Board. It's made of polymerised wood for a strong feel and high user weight. The inclination (17,5°) is perfect for physiotherapy exercises.

Reference	Description	Unit
03-040001	MoVeS Wooden Incline Board - Ladder	1 pc
03-040003	MoVeS Wooden Incline Board - Fix	1 pc



Plastic Incline Board

Apart from stretching the calf and upper-leg muscles, the incline board is a fantastic device to improve general posture and extension in the lower back.

Made of a composite plastic for a maximum load of 250kg. Adjusts to 4 levels and folds flat (8cm) to store easily. Angles are 5°, 15°, 25° and 35°. The last slot on the board makes an excellent hand-hold for stretching. Non-slip top measures 35x35 cm.

250kg

Reference	Description	Unit
03-040102	Incline Board Standard	1 pc



Jumping Fitness Trampoline



The improved Mambo Max Jumping Fitness Trampoline is equipped with bungee cords and a handlebar. Bungee cords protect and support joints and ligaments by providing a deeper and gentler bounce. This allows the muscles to work more, and for a longer period of time, burning calories efficiently without placing unnecessary strain on the body.

Similar to our complete range of trampolines, the Mambo Max Jumping Trampoline helps to flush the lymphatic system, where most toxins are stored. Rebounding also provides an ideal balance of gravity and weightlessness to stimulate the organs and work major and minor muscles in the body. It also serves other core functions, from digestion to mental health.

Finally, bungee trampolines tend to be quieter, sturdier and more reliable than spring trampolines, making upkeep and maintenance a breeze.

Ø 106cm / ↕ 32cm
 ↕ 74-109cm (handrail height from mat)
 200kg

Ready to use



Metal hooks

Bungee cords

Foldable Legs



Reference	Description	Unit
03-111102	Mambo Max Jumping Fitness Trampoline	1 pc

B-Pod Trampoline



Ø trampoline 96cm / Ø C-Hoop: 58cm
 ↕ mat: 26cm
 ↕ C-hoop: 50-80cm
 (C-hoop height from mat)
 120kg

The Health Bounce Pod™ is the only device on the market able to provide a safe and suitable environment for elderly and injured people to exercise. It is robust, balanced, adjustable, and comfortable. It has been designed by professional trampoline athletes in cooperation with specialists in the field of extreme trauma and accidents.

The adjustable C-hoop, connected to the steel frame by 4 symmetrically placed springs makes it friendly to use by people with sight loss, who can exercise freely and up to their own capability. No other form of exercise can be so easily available with true benefits, with the extra bonus of being fun!

Reference	Description	Unit
03-110202	B-Pod Trampoline	1 pc

Triggerpoint Rollers



Massage rollers are used to put pressure on the muscles' triggerpoints and to loosen the tension in the muscles. Trigger points are often the cause of different pains.

SPIKY FOAM ROLLER



↕ 33cm
 Ø 13cm



The Mambo Max Spiky Foam Roller is designed for deeper massage & acupressure therapy. The spikes make it possible to get deeper into the muscles to relieve tension and soreness and to improve mobility. It also has a bi-directional zone for comfort when rolling along the spine. The Spiky Foam Roller is mainly used for stretching, self-massage, body support and stability exercises.

HOLLOW FOAM ROLLER


The Mambo Max Hollow Foam Roller can easily help people who are experiencing long term pains. Its three-dimensional pattern with different shapes replicates the massage from a therapist's hands. The larger surfaces allow a superficial massage while the smaller shapes are designed for a deeper massage.



↔ 33cm
 Ø 14cm



ERGONOMIC FOAM ROLLER




↔ 33cm
 Ø 15cm



The Mambo Max Ergonomic Foam Roller is a curved hollow foam roller, which enables more exercises. Except rolling back and forth, it's possible to make sideways movements. More muscles will be addressed through which the body will be toned even more. The soft foam in combination with the nipples is perfect for trigger point massages. The central groove will reduce pressure on the spine.

3-IN-1 FOAM ROLLER

The Mambo Max 3-in-1 foam roller combines 3 different foam rollers in one package. Thanks to its smooth surface, the outer part is perfect for a superficial massage. The second part has grooves which can be used for a deeper massage and to exercise higher pressure on the trigger points. Finally, with the inner massage stick it's possible to exercise high pressure on sore muscles to improve blood circulation and the muscles' recovery time. This last item is available as a separate product too: The Mambo Max Massage Stick.



↔ 46cm
 Ø 15cm

↔ 46cm
 Ø 10cm

↔ 46cm
 Ø 4cm

MASSAGE STICK



Reference	Description	Unit
04-050101	Mambo Max Hollow Foam Roller	1 pc
04-050102	Mambo Max Ergonomic Foam Roller	1 pc
04-050103	Mambo Max Spiky Foam Roller NEW	1 pc
04-050301	Mambo Max 3-in-1 Foam Roller	1 set
04-050201	Mambo Max Massage Stick	1 pc



Yoga



This yoga set consists of:
 - 1 pvc yoga mat (173 x 61 x 0,4 cm)
 - 2 EVA foam yoga blocks (7,5 x 15 x 23 cm)
 - 1 yoga strap with buckle (180 x 3,8 cm)

This set is packed in a full colour box.
 Items are also sold separately.



Reference	Description	Unit
04-010201	Mambo Max Yoga Mat	1 pc
04-010202	Mambo Max Yoga Block	5-pack
04-010203	Mambo Max Yoga Strap	1 pc
04-010211	Mambo Max Yoga Set 1 Mat 2 Blocks 1 Strap	Set



Stepping Stones



New Colour Combination

The non-slip 16 cm Mambo Max Stepping Stone is an amazing tool to spice up your workout. They can be used by both children and adults for indoor and outdoor workouts to improve dynamic balance, body awareness and coordination. It can also serve as a foot massage tool while seated.



The Stepping Stones now come in bright and vivid colours: yellow, orange, red, purple, blue and green.
 They are packed in cardboard boxes per 6 pieces and each piece has a different colour.



Reference	Description	Unit
04-030141	Mambo Max Stepping Stone	6-pack



Pilates



SOFT-OVER-BALL



The Mambo Max Soft-Over-Ball is an inflatable ball used for resistance work in Pilates. It is a wonderful tool for adding focus and energy to your workouts. In addition, when slightly deflated, it can be used as a soft prop under the head, neck, feet, back... Now available in 3 different colours and sizes, giving a range of 9 different Soft-Over-Balls.



Reference	Description	Unit
04-010105	Mambo Max Pilates Soft-Over-Ball - 17-19 cm Green	1 pc
04-010113	Mambo Max Pilates Soft-Over-Ball - 17-19 cm Blue	1 pc
04-010107	Mambo Max Pilates Soft-Over-Ball - 17-19 cm Black	1 pc
04-010115	Mambo Max Pilates Soft-Over-Ball - 21-23 cm Green	1 pc
04-010123	Mambo Max Pilates Soft-Over-Ball - 21-23 cm Blue	1 pc
04-010117	Mambo Max Pilates Soft-Over-Ball - 21-23 cm Black	1 pc
04-010125	Mambo Max Pilates Soft-Over-Ball - 25-27 cm Green	1 pc
04-010103	Mambo Max Pilates Soft-Over-Ball - 25-27 cm Blue	1 pc
04-010127	Mambo Max Pilates Soft-Over-Ball - 25-27 cm Black	1 pc



FOAM ROLLER



The Mambo Max Foam Roller has been designed for Pilates training and spine exercises (back aches) and is made of EVA. Ideal for balance and stabilisation exercises. Foam rollers are used for stretch and strength exercises, as well as for massage and relaxation.

Ø 15 cm



Reference	Description	Unit
04-010102	Mambo Max Pilates Foam Roller - EVA 90cm	1 pc
04-010104	Mambo Max Pilates Foam Roller - EVA 45cm	1 pc



RING



Foam covered fiberglass resistive exercise ring with foam padded grips. For exercises focusing on strength, balance, flexibility and posture. The Mambo Max Pilates Ring adds light resistance to pilates exercises and helps target specific muscle groups. Often used to help tone the arms and thighs, and increase the core strength challenge of an exercise.

Ø 38 cm

Reference	Description	Unit
04-010101	Mambo Max Pilates Ring	1 pc





Sitting Cushions




Mambo Max Sitting Cushions are air-filled cushions for active and dynamic sitting, recommended by therapists. These easy-to-carry air cushions activate intrinsic trunk muscles to support the spine and are suitable for everybody. The surface from these sitting cushions can easily be cleaned with soap and water.



Soft and comfortable surface

Comfort Cushion 33
Ø 33cm
↓ 7cm
250kg



one smooth | one senso side

Standard Cushion 37
Ø 37,5cm
↓ 7cm
200kg



Packed in colour box with needle pump or in polybag only



Needle pump

Faster Blaster pump



one smooth | one senso side

Standard Cushion 33
Ø 33cm
↓ 7cm
200kg



one smooth | one senso side

Sitting Wedge
□ 34cm
↓ 7cm
250kg



one smooth | one senso side


SensAir Cushion
Ø 32cm
↓ 7cm
200kg

Reference	Description	Unit
04-020301	Mambo Max Comfort Cushion 33 Blue + Needle Pump	1 pc
04-020302	Mambo Max Comfort Cushion 33 Red + Needle Pump	1 pc
04-020303	Mambo Max Comfort Cushion 33 Black + Needle Pump	1 pc
04-021301	Mambo Max Comfort Cushion 33 Blue in Polybag	1 pc
04-021302	Mambo Max Comfort Cushion 33 Red in Polybag	1 pc
04-021303	Mambo Max Comfort Cushion 33 Black in Polybag	1 pc
04-020101	Mambo Max Standard Cushion 33 + Needle Pump	1 pc
04-020102	Mambo Max SensAir Cushion + Needle Pump	1 pc
04-020103	Mambo Max Standard Cushion 37 + Faster Blaster Pump	1 pc
04-020201	Mambo Max Sitting Wedge + Faster Blaster Pump	1 pc



Massage Tools

The **Knobble® II** massage tool is made of durable polymer with a band of non-latex. The Knobble® II gives your hand more power when you need it. Use either the tip or the base for applying pressure on yourself or others. The Knobble® II doesn't replace your healing hands, it just helps them do more and last longer.



The **Jacknobber® II** massage tool can be used in many ways to apply deep pressure to myofascial trigger points and areas of muscle tension. The smaller knobs provide more intense pressure while the larger knobs provide a more diffuse pressure. Use one knob at a time to press deeply into taut bands with embedded trigger points; use two knobs at a time on either side of the spine.

The **Index Knobber® II** massage tool is a simple device designed to let you apply deep pressure to muscles within easy reach of your hand. It is the ideal instrument for clinic or home use allowing for firm, sustained pressure on targeted muscle. It's unique design allows it to be used in several hand positions and by the right or left hands comfortably and interchangeably.



The **Palmassager™** is one of the most efficient massage therapy aids available anywhere. Hold it firmly in the palm of the hand and use one, two or three of its round knobs to massage sore muscles. It is used for back, shoulder, arms and legs.

Reference	Description	Unit
04-031101	Jacknobber II	1 pc
04-031102	Index Knobber II	1 pc
04-031103	Palmassager	1 pc
04-031104	Backnobber II	1 pc
04-031105	Knobble II NEW	1 pc

The **Backnobber® II** massage tool is allows the user to apply as much pain-releasing pressure as long as required to achieve remarkable effects. A locking mechanism allows it to come apart for convenient travel and storage.



Massage Ball, Roll & Ring

The Mambo Max Massage Balls, Roll & Ring are used for relaxing muscles and reducing tension by relieving trigger points. The Mambo Max Massage Ball helps to strengthen hands and feet and to improve coordination and flexibility. It can be used for massage as well as to improve dexterity skills.



Reference	Description	Unit
04-030103	Mambo Max Massage Ball 6 cm	1 pc
04-030105	Mambo Max Massage Ball 7 cm	1 pc
04-030107	Mambo Max Massage Ball 8 cm	1 pc
04-030109	Mambo Max Massage Ball 9 cm	1 pc
04-030111	Mambo Max Massage Ball 10 cm	1 pc
04-030121	Mambo Max Massage Roll 16 cm	1 pc
04-030131	Mambo Max Massage Ring	1 pc






Gym Mat



Comfort Gym Mat
■ ■ ■ 140 x 60 x 1 (cm)




Comfort

1cm

Latex-free


Xtra Comfort Gym Mat
■ ■ ■ 180 x 60 x 1,5 (cm)



Xtra Comfort

1,5cm

Xtra Comfort Gym Mat
■ ■ ■ 180 x 100 x 1,5 (cm)



Xtra Comfort

1,5cm

The Mambo Max Gym Mats provide a soft touch and comfortable exercise environment. They are perfect for gym & fitness. Non-slip, well-cushioned surface with a unique closed cell structure. The mats are equipped with plastic rings 40cm apart so they can be suspended while not in use. Each mat comes with a handy storage/carry strap.

The Comfort Gym Mat is a new addition to our line-up and comes in 3 popular colours: black, blue and magenta. Our well-known 1,5cm-thick exercise mats are now called 'Xtra Comfort Gym Mat' to highlight the extra comfort they offer compared to thinner mats. They continue to be available in widths of 60cm and 100cm.



With rings
and carry strap.



May be used in water

The Mambo Max Gym Mat
Hanger holds up to 10 mats.
Distance between bars is 40cm.

Reference	Description	Unit
05-030001	Mambo Max Gym Mat Wall Hanger for 10 mats	1 pc
05-030303	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) Blue Rings Strap	1 pc
05-030304	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) Black Rings Strap	1 pc
05-030305	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) Magenta Rings Strap	1 pc
05-030101	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (cm) Red Rings Strap	1 pc
05-030102	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (cm) Green Rings Strap	1 pc
05-030103	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (cm) Blue Rings Strap	1 pc
05-030201	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Red Rings Strap	1 pc
05-030202	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Green Rings Strap	1 pc
05-030203	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Blue Rings Strap	1 pc



Fabric Elastic Band



The Mambo Max Fabric Elastic Band consists of 10 consecutive loops. This great toning and work-out tool is suitable for early training up to professional athletes. Three different resistances are available, while the 10 consecutive loops allow you to spice things up. Physical training can be intensified by using a lower amount of consecutive loops resulting in a higher resistance.



The Mambo Max Fabric Elastic Band is durable and extremely versatile. It can be used for stretching, toning or even as an alternative for dumbbells as a weight training accessory. The progressive resistance offers you a broad range of exercises.



COLOUR	RESISTANCE	REFERENCE	UNIT
Light	7 kg	05-120107	1 pc
Medium	10 kg	05-120110	1 pc
Heavy	15 kg	05-120115	1 pc



Dynadome



Half sphere, inflatable exercise dome. This versatile product enables the user to perform exercises on both sides of the product, allowing it to be used for a full body work-out as well as for balance exercises. Basically it's Mambo Max' economical answer to the BOSU® system, but with extra tubing and handles!



Ø 59cm
± 21-25cm
175kg

Tubes with
Handles included




Reference	Description	Unit
05-040401	Mambo Max Dynadome Tubing with handles Faster Blaster Pump	1 pc



Aerobic Step



Step aerobics is a form of aerobics that is practised worldwide by millions of people. The reason for its popularity is simple, it offers a great cardio workout without the need for extensive space or equipment. While the aim of aerobics is to increase the body's consumption of oxygen, step aerobics builds on that goal by offering a more intense workout. The Mambo Max Aerobic Step provides this extra boost by incorporating the action of stepping onto the platform to intensify the workout. It is low impact, and doesn't stress the joints as much as running or jogging.



96cm x 33cm
 15-20-25cm
 100kg



Reference	Description	Unit
05-070101	Mambo Max Aerobic Step	1 pc



Ball Accessories



Electric Compressor Pump



Faster Blaster Pump



Needle Pump




Foot Pump



Hand Pump

Ball Stacker
 with velcro
 closure



Ball Stabiliser



The Ball Stabiliser facilitates balance. Its notches enable combined AB Gym Ball / MoVeS exercises.

Reference	Description	Unit
05-020101	Mambo Max Ball Stabiliser	1 pc
05-020102	Mambo Max Ball Stacker	3 pcs
05-020201	Mambo Max Faster Blaster Pump 2x 120mm	1 pc
05-020202	Mambo Max Hand Pump 2x 2000mm	1 pc
05-020203	Mambo Max Needle Pump 8" - 100mm	1 pc
05-020204	Mambo Max Foot Pump	1 pc
05-020205	Electric Compressor Pump NEW	1 pc



AB Gym Ball



ANTI-BURST

350kg
 Yes, 350 kg!!!

The Mambo Max AB (Anti-Burst) Gym Ball increases flexibility, improves strength, fitness and posture. Exercising with this ball also enhances your balance and coordination. All AB Gym Balls are Anti-Burst: if the ball gets punctured by a sharp object the ball won't burst, thus reducing risk of injury.

This extremely durable Mambo Max AB Gym Ball comes in a Mambo Max cardboard retail pack, including AB Gym Ball, Faster Blaster Pump, Exercise Chart, Plug and Plug Remover.



Faster Blaster Pump

The Mambo Max AB Gym Ball is also available in the budget pack: a polybag containing only the ball and plug.

Budget Pack

Ecofriendly pvc. Heavy metal content in compliance with EN 71 Standards. Phthalates content meets EU-directive 2005/84/EC.

Ball Size	Colour	Body Height	Mambo Max AB Gym Ball Box + Faster Blaster Pump	Mambo Max AB Gym Ball Polybag	Unit
45cm	Yellow	< 125 cm	05-010102	n/a	1 pc
55cm	Red	125-150 cm	05-010103	05-010203	1 pc
65cm	Green	150-175 cm	05-010104	05-010204	1 pc
75cm	Blue	175-185 cm	05-010105	05-010205	1 pc
85cm	Black	185-200 cm	05-010106	n/a	1 pc
95cm	Silver	> 200 cm	05-010107	n/a	1 pc



AB Peanut Ball



ANTI-BURST



Faster Blaster Pump

The Mambo Max AB (Anti-Burst) Peanut ball helps users who struggle with the instability of a regular gym ball. The peanut shape offers extra sideways stability, which gives more confidence to start working out with Gym Balls.



150kg

The Mambo Max AB Peanut Ball offers a more predictable therapy or fitness training as it can only roll in one direction. This certainly is a huge benefit for children, elderly or people who suffered an injury. In addition, Peanut Balls are very popular amongst pregnant women and are very helpful devices during labour.

Similar to our regular Mambo Max AB Gym Ball, the Peanut Ball is Anti-Burst: if the ball gets punctured by a sharp object the ball won't burst, thus reducing risk of injury.

The Mambo Max AB Peanut Ball comes in a Mambo Max cardboard retail pack, including Faster Blaster Pump and Plug.

Peanut Ball Size	Colour	Mambo Max AB Peanut Ball + Faster Blaster Pump	Unit
50 x 100cm	Red	05-011103	1 pc



Balance Board



The Mambo Max Balance Board is an effective tool for balance training and fitness exercises. Stand on the rounded platform and start balancing. It is also a great tool to strengthen your ankles and thigh muscles. For a more severe balance training, add the included screw-on top to your Mambo Max Balance Board.

Ø 42cm
 ↓ 5-8cm
 100kg



Reference	Description	Unit
05-040201	Mambo Max Balance Board	1 pc



Wooden Balance Boards

The Mambo Max Wooden Balance Board is your perfect partner if you want to improve your balance on an ecologic and high-end product. The Wooden Balance board is available in 3 versions: Uni-Directional, Standard & Sport. Their different form and balance point make it possible to have 3 products with a very different character and level of difficulty.

Children and elderly could experience a balance board with pivot point as too challenging. The Uni-directional Balance Board offers more confidence. This enables the user to improve body stability without the fear of falling.

A Uni-directional Balance Board offers a more predictable training as it can only pivot in one direction. When turning the Balance board upside-down, the arches can be used as push-up grips. Once the user is ready for the next level of training, the Balance Board Rehab or Sport are ideal next steps in the process.



The Balance Board Standard is lower than the Balance Board Sport. In addition with its smoother pivot point, it reacts completely different than the Sport. It's a great middle step between the Balance Board Uni-Directional and Balance Board Sport.

The most challenging Balance Board is our Sport version, which, as the name specifies, makes it a real sport. Finding your balance and avoiding the edges touching the ground is hard and requires you to use all your muscles in your body. It's a great and fun way to exercise a full-body workout.

Reference	Description	Unit
05-040001	Mambo Max Wooden Balance Board - Standard	1 pc
05-040002	Mambo Max Wooden Balance Board - Sport	1 pc
05-040003	Mambo Max Wooden Balance Board - Uni-directional	1 pc

Uni: 40cm x 40cm / ↓ 10cm
 150kg



Standard: 40cm x 40cm / ↓ 6cm
 Max inclination: 12,5°
 150kg



Sport: Ø 40cm / ↓ 9cm
 Max inclination: 20°
 150kg



Balance Pad



This foam balance pad provides the comfortable exercise environment needed to perform reinforcing balance exercises. Its soft and thus unstable surface challenges the user to work the foot-, ankle-, and leg muscles.

Reduced in price

Available in three versions:

- 1 rectangular right angle shaped pad, 47 x 39 x 6 (cm), blue
- 1 rectangular chamfered shaped pad, 47 x 39 x 6 (cm), anthracite
- 2 oval shaped pads, 37 x 22 x 6 (cm) - sold per pair



Reference	Description	Unit
05-040301	Mambo Max Balance Pad Rectangular Blue	1 pc
05-040303	Mambo Max Balance Pad Rectangular Anthracite	1 pc
05-040302	Mambo Max Balance Pad Oval Blue	pair



Balance Trainer



The Mambo Max Balance Trainer is designed for balance training. Use it to improve your sense of equilibrium, coordination, reaction and motor skills. It is also ideal for muscle exercise. The Mambo Max Balance Trainer can be adjusted by inflating and deflating the air through the valve. It can be used for lying, kneeling, sitting and stepping movement and it can be incorporated with other exercise equipment to spice up your rehabilitation.

Can be used as Sitting Cushion

Ø 45cm
 ↓ 16cm
 100kg



Ø 60cm
 ↓ 14cm
 100kg



Reference	Description	Unit
05-040102	Mambo Max Balance Trainer 45	1 pc
05-040103	Mambo Max Balance Trainer 60	1 pc



Balance Board Freeman



The Mambo Max Balance Board Freeman is an effective tool for balance training. Stand on the rounded platform and start balancing. This is a great way to strengthen your ankles and thigh muscles. The handgrips on the bottom are used for upper body balancing and workouts.



Ø 40cm
 ↓ 8cm
 100kg

Reference	Description	Unit
05-040202	Mambo Max Balance Board Freeman	1 pc

Body Twister



The Mambo Max Body Twister is your ideal partner to slim your thighs, waist and hips. While twisting your body, the disc below your feet turns and glutes, abdominals and legs are being trained. Continue to twist your body from side to side to tone all muscles in your body. The Mambo Max Body Twister is also an effective product to enhance balance and coordination.

Ø 25cm
 ↑↓ 3cm
 100kg



Reference	Description	Unit
05-100101	Mambo Max Body Twister	1 pc





Ekilibre

BEGINNER



The Mambo Max Balance Plate “Ekilibre” is designed for balance training. Combining balance pad with balance board, it offers the best of both worlds. It improves your sense of equilibrium, coordination, reaction and motor skills. It can be used for lying, kneeling, sitting and stepping movement and it can be incorporated with other exercise equipment.




 60cm x 39cm
 ↓ 9cm
 100kg



INTERMEDIATE



Ø 39cm
 ↑↓ 9cm
 100kg



Reference	Description	Unit
05-040501	Mambo Max Ekilibre Beginner *	1 pc
05-040502	Mambo Max Ekilibre Intermediate *	1 pc

* Can not be sold in Canada.



FOR A QUICK RECOVERY

NAQI®
SPORTS LOTIONS

RECOVERY GEL

OPTIMAL RECOVERY



EXERCISE CAREFREE

- Reduces **muscle fatigue, stiffness** and **cramps**
- Improves the **circulation** and **reduces waste materials**
- Lightens **heavy legs**
- **Cooling** and **analgesic** effect



AFTER SPORTS

NAQI nv Stadsbeemd 1037, 3545 Halen, Belgium www.naqi.com Tel.+32(0)13.460.190

Wooden Work-Out Station



With the Mambo Max Wooden Work-Out Station an even higher amount of exercises is suddenly available. The cleverly spread rings allow you to attach a resistive exercise band or tubing for various vertical and diagonal exercises.

The Work-Out Station is designed to fit an Aerobic Step or a set of Oval Balance Pads. This mix of various products results in a complete and challenging workout. Combining balance with resistive exercise triggers many different muscles in your body. The Mambo Max Wooden Work-Out Station truly offers a complete and three-dimensional workout.



Combine with Balance Pad, Aerobic Step, Stepping Stones... for a huge variety of exercises

Outer Dimension: 103 x 40 (cm)
Inner Dimension: 96 x 33 (cm)



Reference	Description	Unit
05-110001	Mambo Max Wooden Work-Out Station + Tube with Handles + Carabiner	Set
05-110099	Mambo Max Pair of Tubes with Handles + Carabiner	Pair



FAIRS



MEDICA
Düsseldorf - Germany
16 - 19 November 2020

MEDICA has been the largest medical trade fair in the world and the event on the agenda of medical experts for over 40 years. In 2019, the fair welcomed 5,100 exhibitors, nearly 130,000 visitors from 144 countries and had a satisfaction rating for closed deals of 96%.



FIBO
Cologne - Germany
8 - 11 April 2021

FIBO Global Fitness is an expo for sports, nutrition, and healthy lifestyle. It brings global top decision makers of the industry together looking for the latest innovations and trends for their gyms, wellness facilities, doctors and physiotherapy practises, etc...



Therapie Hamburg
Hamburg - Germany
22 - 23 October 2021

Therapie Hamburg is the new platform in northern Germany for physiotherapists, occupational therapists and sports therapists as well as speech therapists, masseurs and lifeguards. Innovations, products and industry trends are shown here.



ARAB HEALTH
Dubai - UAE
21 - 24 June 2021

Arab Health is the largest gathering of healthcare and trade professionals in the MENA region. With over 55,000 local and international visitors, 4,250 exhibitors and participants from over 150 countries, Arab Health is truly where the healthcare world comes to do business.

Swim Noodle




Noodles can be used as resistance equipment for muscle strengthening or as floating device. They are made of ultra-flexible foam that can bend into any shape. Ideal for aquagym, for fitness exercises and a perfect way to strengthen muscles in the water, especially the abdominals. Also great for relaxation and fun in the pool.

Ø 7cm
↔ 150cm



Reference	Description	Unit
05-050101	Mambo Max Swim Noodle	1 pc



Coated Dumbbell

The Mambo Max Dumbbells are cast iron dumbbells, coated with vinyl or neoprene. The vinyl coating is easy washable or can be treated with a water based disinfectant. The neoprene coating has an added resin to give them a scrubby feel and superior grip.

Available in 7 weights with matching colour code.
The Dumbbell Rack holds 6 Dumbbells.



Description	Weight	Vinyl	Neoprene	Unit
Mambo Max Dumbbell	0,5kg	06-010102	06-010202	pair
Mambo Max Dumbbell	1,0kg	06-010103	06-010203	pair
Mambo Max Dumbbell	1,5kg	06-010104	06-010204	pair
Mambo Max Dumbbell	2,0kg	06-010105	06-010205	pair
Mambo Max Dumbbell	3,0kg	06-010107	06-010207	pair
Mambo Max Dumbbell	4,0kg	06-010108	06-010208	pair
Mambo Max Dumbbell	5,0kg	06-010109	06-010209	pair
Mambo Max Dumbbell rack	n/a	06-000101		1 pc



SoftMeds


Mambo Max SoftMeds are colour coded hand weighted balls that are small and flexible enough to grasp with one hand, yet big enough to hold with both hands. The soft, pliable vinyl cover makes it easy to hold and is a safe and effective way to add weight training to your workout.

Mambo Max SoftMeds are fully inflated and ready to use, straight out-of-the box. When needed, air pressure can be adjusted with a needle pump.



Reference	Weight	Description	Unit
06-040102	0,5kg	Mambo Max SoftMeds Tan	1 pc
06-040103	1,0kg	Mambo Max SoftMeds Yellow	1 pc
06-040104	1,5kg	Mambo Max SoftMeds Red	1 pc
06-040105	2,0kg	Mambo Max SoftMeds Green	1 pc
06-040106	2,5kg	Mambo Max SoftMeds Blue SOLD OUT	1 pc
06-040107	3,0kg	Mambo Max SoftMeds Black LAST UNITS	1 pc





Wrist & Ankle Weight



0,5 - 2kg

To quickly build stamina and muscle, simply add Mambo Max Wrist and Ankle Weights to your workout. They are filled with iron sand and will effectively shape and tone your body and can be worn anywhere. Wear them around the office, out shopping or at the gym. Eventually, your body will adjust to running, walking and living with this extra weight.

3 - 5kg



Reference	Description	Weight	Unit
06-020102	Mambo Max Wrist & Ankle Weights	0,5kg	pair
06-020103	Mambo Max Wrist & Ankle Weights	1kg	pair
06-020104	Mambo Max Wrist & Ankle Weights	1,5kg	pair
06-020105	Mambo Max Wrist & Ankle Weights	2kg	pair
06-020107	Mambo Max Wrist & Ankle Weights	3kg	pair
06-020108	Mambo Max Wrist & Ankle Weights	4kg	pair
06-020109	Mambo Max Wrist & Ankle Weights	5kg	pair



Thumb Lock Wrist Weight



The Mambo Max Thumb Lock Wrist Weight enables you to fix the weight to your thumb, ensuring a fixed position during your entire workout. Being it boxing, running or Zumba, the weights will remain in place.

Reference	Description	Weight	Unit
06-020202	Mambo Max Thumb Lock Wrist Weights	0,5kg	pair
06-020203	Mambo Max Thumb Lock Wrist Weights	1kg	pair



Medicine Ball



A medicine ball (also known as an exercise ball, a med ball, or a fitness ball) is often used for strength training as it plays an important role in the field of sports medicine. Medicine ball training is one of the oldest forms of strength and conditioning training – the first reference to wrestlers training with sand filled bladders appears in Persia nearly 3000 years ago. In ancient Greece the physician Hippocrates had his patients throw them back and forth for injury prevention and rehabilitation.

Reference	Description	Ø	Unit
06-030103	Mambo Max Medicine Ball - 1kg	20cm	1 pc
06-030104	Mambo Max Medicine Ball - 1,5kg	20cm	1 pc
06-030105	Mambo Max Medicine Ball - 2kg	20cm	1 pc
06-030107	Mambo Max Medicine Ball - 3kg	22cm	1 pc
06-030108	Mambo Max Medicine Ball - 4kg	22cm	1 pc
06-030109	Mambo Max Medicine Ball - 5kg	27cm	1 pc
06-000301	Mambo Max Medicine Ball rack (6 balls)		1 pc



28cm x 57cm x 183cm
(l x w x h)



Kettlebell

The Mambo Max Kettlebells are an important part of any fitness or CrossFit training. Made from high quality iron surrounded with a colourful vinyl coating, they will not damage any surfaces when being stored. Starting from 2 kg up to 24 kg, we also offer you a complete set, including a premium quality Kettlebell Rack, which can hold up to 150kg per shelf (total of 300kg).



120cm x 50cm x 70cm (l x w x h)

Available in: 2 - 4 - 6 - 8 - 10 - 12 - 16 - 20 - 24kg

Reference	Weight	Description	Unit
06-050102	2kg	Mambo Max Kettlebell Sky Blue	1 pc
06-050104	4kg	Mambo Max Kettlebell Pink	1 pc
06-050106	6kg	Mambo Max Kettlebell Yellow	1 pc
06-050108	8kg	Mambo Max Kettlebell Red	1 pc
06-050110	10kg	Mambo Max Kettlebell Green	1 pc
06-050112	12kg	Mambo Max Kettlebell Blue	1 pc
06-050116	16kg	Mambo Max Kettlebell Black	1 pc
06-050120	20kg	Mambo Max Kettlebell Silver	1 pc
06-050124	24kg	Mambo Max Kettlebell Orange	1 pc
06-051001		Mambo Max Kettlebell Rack	1 pc
06-051101		Mambo Max Kettlebell Rack with 1 Kettlebell of each weight	set



Para Pro - Paraffin Bath



The Para Pro is a professional paraffin wax heating bath which is made to be used on 220/230 Volt power sockets and which can be used 24/7. In combination with our premium paraffin pearls it only needs between 2 and 3 hours to completely melt 2,5kg of paraffin. The temperature is adjustable from 40-60°C and is easy to control due to the digital thermometer.



When ordering the Para Pro, you receive 2,5kg of white paraffin pearls or 6 blocks (2,72kg) of Tea Tree Green Paraffin free of charge.

Para Pro Specifications:

- Outer dimensions: 39cm x 25cm x 18cm (l x w x h)
- Inner dimensions: 32,5cm x 18cm x 13cm (l x w x h)
- Weight: 3,1kg
- 220 V: No transformer needed, which will lengthen the durability.
- Low power consumption once the paraffin has melted
- 24-hours a day fully operational
- **Fast melting between 2 - 3 hours!**
- Adjustable temperature (**40 - 60°C**)
- Digital thermometer

Paraffin Blocks:

- Green fully refined paraffin wax
- Melting point: 46°C
- Tea Tree scented



NEW
Comes with UK or EU plug!

Paraffin Pearls:

- White fully refined paraffin wax
- Congealing point: 53°C
- Melting point: +/- 54°C
- Meets the purity requirements of BGA1, Ph.EUR.2 and FDA3

1. Bundesgesundheitsamt recommendation No. XXV B II;
2. European Pharmacopoeia, current edition;
3. Food and Drug Administration, § 172.886 Petroleum Wax



Reference	Description	Unit
07-100101	MoVeS Para Pro Bath + 2,5kg Paraffin Pearls - White EU-plug	1 pc
07-100101 UK	MoVeS Para Pro Bath + 2,5kg Paraffin Pearls - White UK-plug	1 pc
07-100206	MoVeS Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg) - Tea Tree Green EU-plug	1 pc
07-100206 UK	MoVeS Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg) - Tea Tree Green UK-plug	1 pc
07-110102	MoVeS Paraffin Pearls - 2,5kg - White	1 pc
07-110125	MoVeS Paraffin Pearls - 25kg - White	1 pc
07-110206	MoVeS Paraffin Blocks 46° Set of 6 x 453g (2,72 kg) - Tea Tree Green	1 pc
07-109801	MoVeS Insulating Gloves for Para Pro Bath Blue	1 pair
07-109802	MoVeS Disposable Gloves for Para Pro Bath	10 pairs



PearlPack



PearlPacks are reusable Hot/Cold Packs filled with small pearls. Hot and Cold Therapy can be applied to relieve pain and inflammations associated with muscle strain and sprain. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The small pearls absorb heat or cold, which they can retain for 20 minutes. During these 20 minutes, these packs will relieve a part of the pain.

Not only do they look good, the products are very convenient in use. PearlPacks won't feel wet nor drip after freezing or heating them, avoiding a mess. All materials used are non-toxic and odourless.

Besides relieving pain from an injury or muscle strain, hot/cold packs are often used to reduce swelling after an operation. PearlPacks truly are multifunctional products.



07-040301: PearlPack Knee/Elbow

This offers a solution for Hot/Cold therapy on the knee or elbow joints. The large opening ensures that your knee or elbow remains free and can continue to move freely, while the Pearl technology sends direct relief to the pain. The flexible velcro strap keeps the pack firmly in place.



07-041101: PearlPack Eye Mask

This mask adapts to the patient's face and has a fabric strap. Hot treatment can be applied when eyes feel dry or if there's an infection of the eyelid. Cold treatment will help when suffering from headaches, swollen eyes or sinusitis.

07-040201: PearlPack Back Wrap

This large pack divided in 3 segments is equipped with velcro strap. It fits perfectly around the patient's waist for an optimal pain relief on the back. The pearls can be pushed from one segment to the other if one area of the back needs less or more thermal treatment.



07-040101: PearlPack Sport Pack

This portable classic pack can be applied on any part of the body.



Reference	Description	Dimension	Unit
07-040101	MoVeS PearlPack - Sport Pack	19cm x 12,5cm	1 pc
07-040201	MoVeS PearlPack - Back Wrap	43cm x 17,5cm	1 pc
07-041101	MoVeS PearlPack - Eye Mask	22,5cm x 7,5cm	1 pc
07-040301	MoVeS PearlPack - Knee/Elbow NEW	37cm x 25,6cm	1 pc





Roll-On Pack



MoVeS Roll-On packs combine hot/cold treatment with practicality. The tubular shape and flexible material allow an easy roll-on application, even after being heated or frozen. The various sizes will fit perfectly over finger, wrist, ankle or knee. This design ensures that the whole area will benefit from the hot or cold treatment.

From the same factory as the original Torex® products. Reusable, non-toxic and latex free. Available in 5 sizes.



Reference	Description	Inside Ø	Circumference	Uses:	Unit
07-050101	MoVeS Roll-On Pack Finger	2,5cm	One size fits all	Finger	1 pc
07-050102	MoVeS Roll-On Pack Small	10,0cm	Fits up to 25cm	Hand, Wrist, Arm, Elbow	1 pc
07-050103	MoVeS Roll-On Pack Medium	12,5cm	Fits 25cm - 38cm	Hand, Wrist, Arm, Elbow, Ankle	1 pc
07-050104	MoVeS Roll-On Pack Large	15,0cm	Fits 38cm - 53cm	Elbow, Ankle, Foot, Calf, Knee	1 pc
07-050105	MoVeS Roll-On Pack Extra Large	17,5cm	Fits 53cm - 58cm	Knee, Thigh	1 pc



Torex Pro Ankle Boot



Torex® Pro Ankle Boot
A high performance ankle hot/cold and compression brace designed by PRO Athletic Trainers.



Reference	Description	Unit
07-030201	TorexHot/Cold Pro Ankle Boot	1 pc



Hot/Cold Pack



The MoVeS Hot/Cold Packs are usable for hot as well as cold therapy. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The versatile, reusable packs deliver hot and cold therapy effectively and economically. The use of the MoVeS Hot/Cold Pack is unlimited and its material is non-toxic. Even below -20°C it remains soft. The MoVeS Hot/Cold Pack can be bent in any angle, allowing maximum coverage of joint and muscle. All Hot/Cold Packs can be heated in the microwave and can be frozen in the freezer. Temperature holds for approximately 30 minutes.

MoVeS Hot/Cold Pack | Classic - Gel

The classic (and economical) transparent Hot/Cold Pack. The special gel is non-toxic and biodegradable.

MoVeS Hot/Cold Pack | Standard - Gel

Our standard Hot/Cold Pack is a nylon covered gelpack. The special gel is non-toxic and biodegradable.

MoVeS Hot/Cold Pack | Soft Touch - Clay

The Soft Touch is improved with a more user friendly surface. The frost-free surface allows usage in direct contact with the skin. Filled with clay, a 100% natural product.

Sold per 10 pieces



Sold per 10 pieces



Sold per 10 pieces



100% polyester protection sleeve for hot/cold packs. Machine washable. Available in 3 sizes.



07-040201: XXL Hot/Cold Pack

This is ideal for pain relief on shoulders and back and has a nylon cover.
Sold per single piece!



07-010311: Tri-Sectional Hot/Cold Pack



Dimension		MoVeS Hot/Cold Packs			Protection Sleeve	Unit
		Classic (gel)	Standard (gel)	Soft Touch (clay)	100% polyester	
15cm x 25cm	Small	07-010102	07-010202	07-010302	07-020102	10 pcs
20cm x 30cm	Medium	07-010103	07-010203	07-010303	07-020103	10 pcs
25cm x 35cm	Large	07-010104	07-010204	07-010304	07-020104	10 pcs
33cm x 47cm	XXL	n/a	07-010206	n/a	n/a	1 pc
40cm x 20cm	Tri-Sectional	n/a	07-010211	07-010311	n/a	10 pcs

Infrared Thermometer

EAR & FOREHEAD thermometer

The infrared thermometer is used for measuring body temperature by two-in-one functions: Ear and Forehead temperature. It is suitable for medical unit and home use.

- Fast and easy temperature measurement by advanced infrared technology
- Range between 34.0°C~43.0°C
Accuracy; <35.0°C and >42.0°C = ±0.3°C
35.0°C~42.0°C = ±0.2°C
- Tri-colour backlit;
 - <37.4°C = Green
 - 37.4°C~38.0°C = Yellow
 - >38.0°C = Red
- Memory mode
- Suitable for both infants and adults
- Easy operating and cleaning
- Automatic power off after 60 seconds
- 2 AAA alkaline batteries



SHOTGUN thermometer

Measures your body temperature without direct contact with the skin. Can measure between 5-15 cm from the skin.

- Detailed 34 memories with time/date
- Response time: 1 second
- Effective distance: 5-15cm
- User guidance by infrared beam
- Fever alarm



Reference	Description	Unit
08-070301	Infrared Thermometer Ear & Forehead auto-sensing 2 x AAA	1 pc
08-070303	Infrared Thermometer Shotgun Forehead auto-sensing 2 x AAA NEW	1 pc



3-colour
backlit



Blood Pressure Monitor



The blood pressure monitor is used for measuring blood pressure and pulse rate. It is suitable for people older than 12 years.



Features:

- Large LCD display: 72mm x 43,5mm
- 2 users, 99 memories per user
- Irregular heart beating sign
- Indicating blood pressure level in 4 colours
- Average of last 3 records
- Low power
- Over pressure protection
- Automatic switch-off



Reference	Description	Unit
08-070501	Blood Pressure Monitor + Cuff NEW	1 pc

Fingertip Pulse Oximeter

A Pulse Oximeter is a very important and common device to check patient blood-oxygen saturation (SpO₂) level and pulse rate. As a small, compact, simple, reliable and durable physiological monitoring device, Fingertip Pulse Oximeter greatly enhances patient care.



ADULT oximeter

Features:

- Dual colour OLED displays SpO₂, PR, pulse bar and waveform
- SpO₂; measurement range between 70%~99%
Accuracy; ±2% on the stage of 80%~99%
- PR; measurement range between 30BPM~240BPM
Accuracy; ±1BPM or ±1% (the larger one)
- Low power indicator
- 6 display modes with adjustable brightness
- Automatically power off
- AAA alkaline batteries



JUNIOR oximeter

Rechargeable oximeter with OLED-screen, specifically designed for kids.

- Rechargeable via USB-cable
- Dual colour OLED displays SpO₂, PR, pulse bar and waveform
- 4-direction & 6-mode display for convenient readings
- Set-up alarm range of SpO₂ and pulse rate
- Menu-function setting (beep sounds, etc.)
- Lithium battery; automatically power off
- Gravity function, automatic rotating reading option



Reference	Description	Unit
08-070201	Fingertip Pulse Oximeter Adult OLED 2 x AAA	1 pc
08-070202	Fingertip Pulse Oximeter Junior OLED Rechargeable NEW	1 pc



Hydraulic 3-piece hand evaluation kit



Portable carrying case contains all instruments necessary to evaluate the strength and range-of-motion of the hand and fingers.

Set includes:

- hydraulic hand dynamometer
- hydraulic pinch gauge
- stainless steel finger goniometer 15cm - 180° per 5°

Quality Hard Case Box

Reference	Description	Unit
08-010113	Saehan Hydraulic Hand Evaluation Kit 3 pieces	Set



Hydraulic 7-piece hand evaluation kit



Portable carrying case contains all instruments necessary to evaluate the strength and range-of-motion of the hand and fingers.

Set includes:

- hydraulic hand dynamometer
- hydraulic pinch gauge
- stainless steel finger goniometer 15cm - 180° per 5°
- 2-point discriminator with 3rd point
- Wartenberg pinwheel
- finger gauge
- functional finger motion gauge

Quality Hard Case Box

Reference	Description	Unit
08-010117	Saehan Hydraulic Hand Evaluation Kit 7 pieces	Set



Hydraulic hand dynamometer



The Hydraulic Hand Dynamometer gives accurate grip strength readings without the patient being able to "feel" the handle move. The adjustable handle has 5 positions to accommodate to any hand size. The indicator remains at the patient's maximum reading until it is reset. Has a scale in kilograms. Comes in a hard carrying case.

Range: 0 - 90kg (200lb)

Quality Hard Case Box



Reference	Description	Unit
08-010101	Saehan Hydraulic Hand Dynamometer	1 pc



Digital hand dynamometer with Data Transfer Software- DHD-1



Designed for reliable grip-force evaluations from 0 to 90kg (200lb). It has an easy-to-read LCD display. This hand dynamometer features an electronic zero calibration system and up to 40 test data can be stored in memory. It is powered by a rechargeable Li-Ion polymer battery and data will not be lost when changing the battery. The Data Transfer Software can transfer the test results into .txt-files, which can be easily accessed using Notepad or MS-Word for editing and printing.

This Digital Hand Dynamometer hard carrying case includes a Digital Hand Dynamometer, the DTS-software and a USB cable.

Range: 0 - 90kg (200lb)

Quality Hard Case Box



Reference	Description	Unit
08-010201	Saehan Digital Hand Dynamometer - DHD-1	1 pc



Digital hand dynamometer with G-STAR™ Software - DHD-3



The DHD-3 Digital Hand Dynamometer (SH1003) provides the same features as the DHD-1 when used as a stand-alone device, but adds automated Rapid Exchange Grip Testing. Using Auto Mode, up to 10 strength tests for each hand as well as the time allowed to switch hands can be configured to automate the testing procedure.

The Grip Strength Testing and Research (G-STAR) is an easy to use software tool which automates grip strength testing and endurance testing in both clinical and research environments. All patient information and their test results are stored in an Access database.

Range: 0 - 90kg (200lb) / Minimum requirements: Windows XP + USB port

Reference	Description	Unit
08-010202	Saehan Digital Hand Dynamometer - DHD-3 (incl. G-STAR software)	1 pc



Quality Hard Case Box



Pinch Gauge

SAEHAN®

08-020102: Mechanical Pinch Gauge

Used for hand-strength testing by measuring tip, key and palmar pinch force. Made of durable painted stainless steel. Mechanical properties will not degrade over time. It can serve as a diagnostic tool or for monitoring the progress of rehabilitation.



08-020102

08-020101: Hydraulic Pinch Gauge

Uses the improved hydraulic system of the hand dynamometer to assure convenience, product reliability, measurement accuracy and repeatability. This results in a more accurate test result for all pinch tests.




08-020101

Range:
0 - 20kg (45lb)

Reference	Description	Unit
08-020101	Saehan Hydraulic Pinch Gauge	1 pc
08-020102	Saehan Mechanical Pinch Gauge	1 pc

CE

Push-Pull Dynamometer



Range: 0 - 45kg (100lb)

Lift platform and chain to perform lifting evaluation not included.



Quality Hard Case Box



Reference	Description	Unit
08-010501	Baseline Hydraulic Push-Pull Dynamometer	1 pc
08-010502	Dual Grip for Push-Pull Dynamometer	1 pc

CE

Squeeze Dynamometer

SAEHAN®



Quality Soft Pouch Case

The Squeeze (bulb) Dynamometer is the inexpensive device to reliably measure hand and finger strength. The dynamometer has an indicator dial that remains at the maximum reading until reset. Unit is calibrated in pounds and kilograms.


Range: 0 - 70kg (150lb)

Reference	Description	Unit
08-010301	Saehan Squeeze Dynamometer 70kg	1 pc

CE

Smedley Hand Dynamometer

SAEHAN®



Ergonomic Handle

The Smedley spring-type Dynamometer has been successfully used for many years to measure grip strength. The ergonomic handle is adjustable to fit all hand sizes. Unit is calibrated in kilograms and pounds. Stored in a soft carrying case.

Range: 0 - 100kg (220lb)

Quality Soft Pouch Case



Reference	Description	Unit
08-010401	Saehan Smedley Hand Dynamometer - Ergonomic Handle	1 pc

CE

Back-Leg-Chest Dynamometer

SAEHAN®

Measure strength of back, leg and chest. Base provides secure footing. Chain length is adjustable to accommodate for height differences or to vary the point of force application. Shows pounds and kilograms. Pointer remains at maximum until reset. Comes fully assembled with 1,5m chain. Range: 0 - 300kg (660lb)



Reference	Description	Unit
08-010601	Saehan Back-Leg-Chest Dynamometer	1 pc

CE

Fat Caliper

Professional light-weight aluminium caliper. Measurement arms with floating tips to assure accurate and reproducible results. 65mm scale is large and easy to read. Corrective scale for sex and age included in multilingual manual.

- Professional, very accurate unit with spring loaded levers and low-pressure bearings for constant pressure.
- Scale permits reading up to 65mm
- Aluminium casted to insure years of reliable use.
- Includes carrying case
- CE marked and FDA approved



08-040101
CE
SAEHAN

Quality Soft Pouch Case



08-040104
CE

The Plastic Skinfold Caliper from Baseline is the inexpensive alternative for skinfold analysis. Spring loaded arm assures reproducible “pinch”. Read directly from the 50mm scale. Ideal instrument for screening large population. The Plastic Skinfold Caliper comes with a booklet and tables.

Enjoy the ability to measure body fat easily, by yourself, in the privacy of your own home, with the reliability and accuracy that you expect from world class training tools. Recommended in Body-for-LIFE and endorsed by the World Natural Bodybuilding Federation, Accu-Measure® Personal Body Fat Testers have gold standard accuracy to within 1.1% of underwater weighing results. Each Accu-Measure® Fitness 3000 model comes with:

- An original Accu-Measure® Personal Body Fat Tester
- Body Fat Measurement Instructions
- Men’s and Women’s Body Fat Measurement Charts
- Body Fat Tracking Charts with the Right Way to Weigh
- Questions and Answers on Body Fat and Calipers and Body Fat Tips



08-040103



08-040102

The FatTrack® FITNESS TK Body Transformation Toolkit enables you to measure (and track) Body Fat %, Lean Body Mass, your workouts, macronutrient information and more by yourself - Easily, accurately and in the privacy of your own home!

- Each FatTrack FITNESS TK Includes:
- Updated FatTrack® Fitness Digital Body Fat Caliper
 - User Manual with instructions
 - Body Tracker Software
 - MyoTape Body Tape Measure
 - FREE copy of “Improve Your Measuring Skill” with practical pointers for getting the best measurement readings and valuable information on ideal weight and muscle mass.

Reference	Description	Unit
08-040101	Saehan Skinfold Fat Caliper Aluminium - Green	1 pc
08-040104	Baseline Skinfold Fat Caliper Plastic - Black	1 pc
08-040103	Accu-Measure Skinfold Fat Caliper Plastic - White	1 pc
08-040102	FatTrack Skinfold Fat Caliper Digital	1 pc

Inclinometer

Baseline’s Bubble and Digital Inclinometers are used to measure range-of-motion (ROM).



Reference	Description	Unit
08-060101	Baseline Bubble Inclinometer	1 pc
08-060201	Baseline Digital Inclinometer	1 pc



Quality Hard Case Box



CE

Scolimeter



Scolimeters are inclinometers designed to measure trunk asymmetry or Angle of Truncal Rotation (ATR) used in scoliosis screening programmes.

Reference	Description	Unit
08-060102	Pedi-Scolimeter	1 pc
08-060104	Baseline Plastic Scolimeter	1 pc

CE

Measurement Tape



Measurement Tape with Gulick Attachment



Finger Circumference Gauge



Measurement Tape with Hands-free Attachment


Reference	Description	Unit
08-080101	Baseline Measurement Tape with Gulick Attachment - 150cm	1 pc
08-080102	Baseline Measurement Tape with Gulick Attachment - 180cm	1 pc
08-080103	Baseline Measurement Tape with Hands-free Attachment - 150cm	1 pc
08-080201	Baseline Finger Circumference Gauge - 15cm	1 pc

CE

Baseline & Lafayette Goniometer



30cm Plastic Goniometer
08-030401




30cm Plastic Goniometer - HiRes™
08-030402



Absolute + Axis® Attachment
08-030411



9cm Plastic Finger Goniometer
08-030601



Lafayette Extendable Goniometer
08-030701

Reference	Description	Unit
08-030401	Baseline Plastic Absolute+Axis Goniometer - 30cm 360°	1 pc
08-030402	Baseline Plastic Absolute+Axis Goniometer - HiRes™ - 30cm 360°	1 pc
08-030411	Baseline Plastic Absolute+Axis - Attachment Only for 30cm Goniometer	1 pc
08-030601	Baseline Plastic Finger Goniometer - 9cm Flexion to Hyper-Extension	1 pc
08-030701	Lafayette Gollehon Extendable Goniometer with Magnifying Lens - 24cm - 66cm 180°	1 pc



Goniometer Kit



The Plastic set includes:

- Plastic Goniometer 15 cm | 0° to 360° per 1°
- Plastic Goniometer 20 cm | 0° to 360° per 1°
- Plastic Goniometer 30 cm | 0° to 360° per 1°
- Rulong Goniometer 18 cm | 0° to 360° per 2°
- Pocket Goniometer 15 cm | 180° x2
- Plastic Flexion-Hyperextension goniometer

Quality Soft
Pouch Case

The Stainless Steel set includes:

- Robinson pocket goniometer 15cm - 180°
- Stainless steel goniometer 18cm - 180°
- Stainless steel goniometer 35cm - 180° x2
- Stainless steel goniometer 35cm - 360°
- Stainless steel finger goniometer 15cm
- X-ray goniometer 20cm - 180°

Reference	Description	Unit
08-030211	Saehan Goniometer Set Stainless Steel 6 pieces	Set
08-030121	Saehan Goniometer Set Plastic 6 pieces	Set



Plastic Goniometer



Clear plastic permits observation of joint's axis of motion, and its range of motion. The 360° head has three scales calibrated to be used with ISOM (International Standard Of Measurements, STFR) system.



15cm Goniometer
08-030101



20cm Goniometer
08-030102



30cm Goniometer
08-030103



Rulong Goniometer
08-030111




Pocket Goniometer
08-030112

Reference	Description	Unit
08-030101	Saehan Plastic Goniometer 15cm 0° to 360° per 1°	1 pc
08-030102	Saehan Plastic Goniometer 20cm 0° to 360° per 1°	1 pc
08-030103	Saehan Plastic Goniometer 30cm 0° to 360° per 1°	1 pc
08-030111	Saehan Rulong Goniometer 18cm 0° to 360° per 2°	1 pc
08-030112	Saehan Pocket Goniometer 15cm 180° x2 per 5°	1 pc



Stainless Steel Goniometer




15cm / 180° scale / 5° increments
08-030201



20cm / 1° increments
2 opposing 180° scales
08-030202



35cm / 2 opposing 180° scales / 1° increments
Arm tension is controlled by thumb knob
08-030203



Extra long arms for measuring all joints
35cm / 360° scale / 1° increments
08-030204

Reference	Description	Unit
08-030201	Saehan Stainless Steel Robinson Pocket Goniometer 15cm 180°	1 pc
08-030202	Saehan Stainless Steel Goniometer 20cm 180° x 2	1 pc
08-030203	Saehan Stainless Steel Goniometer 35cm 180° x 2	1 pc
08-030204	Saehan Stainless Steel Goniometer 35cm 360°	1 pc



Finger Goniometer



Finger and Flexion/Hyper-Extension Goniometers measure range-of-motion of finger joints. Head has two opposing 180° scales marked in 5° increments. A linear scale in centimetres is on arm.

Reference	Description	Unit
08-030301	Saehan Finger Goniometer 9cm Stainless Steel 180° per 5°	1 pc
08-030302	Saehan Finger Goniometer 15cm Stainless Steel 180° per 5°	1 pc



Body Fat Scale



Baseline body fat scale is made with tempered glass. This scale measures weight, body fat and water % and memorises data up to 8 persons. Measurements range from 4% - 45% body fat; 37.8% - 66% body water and 330lb/150kg weight capacity. The Required 1 x 3V battery is included.

Resolution:
- 0.1% bodyfat
- 0.1% water

Dimensions: 35cm x 35cm x 6cm
Product weight: 2kg


Reference	Description	Unit
08-070101	Baseline BMI Body Fat Scale	1 pc



Sit & Reach Unit



When using the Sit & Reach unit you sit with your feet against the "sit and reach box" and flex forward moving the slide bar. Flexion is read directly from the calibrated bar.



Reference	Description	Unit
08-060103	Sit & Reach Unit	1 pc



Monofilament



Monofilament evaluators are used to measure cutaneous sensory perception threshold.

Reference	Description	Unit
08-050401	Baseline Tactile Monofilament - LEAP Program - Disposable - 5.07 - 10 g	1 pc
08-050421	Baseline Tactile Monofilament - 2-piece Set - Normal - 2.83, 3.61 g	Set
08-050425	Baseline Tactile Monofilament - 2-piece Set - Protective Sensation - 4.56, 5.07 g	Set



Tuning Fork

These aluminium Tuning Forks are developed to measure sensitivity to vibration which can be an indicator of nerve damage. It is also used for foot exam by diabetic people to check for signs of neuropathy.

Reference	Description	Unit
08-050301	Tuning Fork 128cps	1 pc
08-050302	Tuning Fork 256cps	1 pc




Neurological Reflex Instruments




08-050101: Buck Hammer

Two heads for cutaneous sensation evaluation. Comes with a needle and a brush.



08-050102: Babinski Hammer

Pointed bottom can be used as gross aesthesiometer.



08-050103: Taylor Hammer

Use to elicit responses during examination
Latex-free!



08-050104: Dejerine Hammer

This hammer is used for testing deep tendon reflexes in a diagnostic evaluation. It features a double-sided head composed of two rounded rubber surfaces for increased patient comfort.



08-050201: Wartenberg Pinwheel

Elicit cutaneous sensory and pain perception responses.



08-050105: Troemner Hammer

This hammer is shaped like a two-headed mallet. The larger mallet is used to elicit tendon stretch reflexes, and the smaller mallet is used to elicit percussion myotonia.



Reference	Description	Unit
08-050101	MoVeS Buck Hammer	1 pc
08-050102	MoVeS Babinski Hammer	1 pc
08-050103	MoVeS Taylor Hammer	1 pc
08-050104	MoVeS Dejerine Hammer NEW	1 pc
08-050105	MoVeS Troemner Hammer NEW	1 pc
08-050201	MoVeS Wartenberg Pinwheel	1 pc



Discriminator



08-050203
Touch-Test 2-point Discriminator


Reference	Description	Unit
08-050203	Touch-Test 2-point Discriminator	1 pc
08-050204	Baseline 2-point Discriminator	1 pc
08-050205	Baseline 3-point Discriminator	1 pc
08-050206	Baseline 2-point Discrim-A-Gon, 2 Disk Set	1 pc



08-050204
Baseline 2-point Discriminator



08-050205
Baseline 3-point Discriminator



08-050206
Baseline 2-point Discrim-A-Gon, 2 Disk Set



TENS device



- The Neuropex TENS is accepted to be far more accurate than a standard TENS machine.
- Highly user friendly with 3 programs and modes:
 - B (burst)
 - C (continuous)
 - M (Pulse Rate Modulation)
 - Regulated pulse width control for adjustable intensity and stimulation strength
 - Clear plastic window allows controls to be viewed during operation
 - Easy to use and very compact
 - Adjustable intensity and stimulation strength
 - Protective case for easy electrode storage and transportation
 - Complete set of electrodes, in-built timer, complimentary spare battery



- Neuropex Long Lasting Self-Adhesive Electrodes
- Resealable pouch for safe storage
 - 4 Electrodes per pack
 - 50mm x 50mm
 - Ideal for use with all Neuropex TENS and Interferential devices



Reference	Description	Unit
08-070401	2-Channel Professional TENS unit	1 pc
08-070411	Replacement TENS Electrode - 4 piece set	set



PRIVATE LABEL



We have nice products, but you would like your own brand or your customer's name on the band, loop, superloop...? We'll do it! Together with our manufacturers from all the products in our range, we can create your desired product with the look you want. We can replace our logo with yours, change the colour, even the packaging can be tweaked to your needs. Private Labels will allow you to design the product the way you want to match it to your company's or customer's identity.



Private Label products have the same quality as all our products, but are putting your brand in the spotlights. This is perfect if you need a homogeneous range of products or want to increase your brand awareness. Or maybe you're looking for a one time opportunity: a promotional product during your company's special occasion, a gift for visitors or an event you're supporting?

Please contact us for more information about pricing and volume.

Product description	Minimum Order Quantity
Resistive Exercise Band	5.000m per colour
Resistive Exercise Loop	500 pcs per colour
Soft-Over-Ball	3.000 pcs per colour
Hollow Foam Roller	300 pcs
AB Gym Ball	500 pcs
NBR Comfort Gym Mat	300 pcs per colour
NBR Xtra Comfort Gym Mat	
Pearl Packs	500 pcs
Plastic Goniometer	2.000 pcs
Various other products	Contact us for more info



MVS IN MOTION CALIBRATION SERVICE

Measurement Instruments are calibrated at manufacturing with a full one year guarantee, before sending them to users all over the world. Ideally, these instruments should be recalibrated every year.

In order to make this happen in an easy and inexpensive way, Hydraulic Hand Dynamometers (Saehan, Jamar, Baseline...) and Hydraulic Pinch Gauges are recalibrated by our experts at the MVS In Motion Centre in Tisselt - Belgium. This operation consists of re-establishing the exact relationship between the force applied to a dynamometer and the measurement it shows. Small repairs can be done as well.



Simply send your request to our customer service department for a return authorisation. After each recalibration, you will receive an updated calibration certificate, certifying accurate and repeatable test results.

If you have many calibrations to do, you might consider buying your own calibration unit.

Reference	Description	Unit
08-019901	Calibration Unit	1 pc
08-019900	Calibration & Repair	-



MVS IN MOTION REPAIR CENTRE



Since machines are designed and assembled by humans, they are not flawless. That is why MVS In Motion installed an inhouse repair centre, to swiftly and correctly solve any problems that may occur with our products.

No matter what the cause, be it faulty fabrication or misuse, our technician will shine his light on the problem.

Send an email to our customer service department to inform us on the type of problem and the serial number/production date of the device. After our authorisation, please send over the patient, accompanied by your name and client number and a brief description of the symptoms to info@mvs-in-motion.com. We'll do our best to have him cured in no time.

3-in-1 Foam Roller	41	Flex-Ion	23	Roll-On Packs	60
Accessories	22	Foam Hand Grip	32	Saehan Cycle	38
Adjustable Hand Grip	32	Foot Pump	48	Safety Handles	22
Aerobic Step	48	Goniometers	70-71	Scolimeter	69
Ankleciser	19	Gym Ball	49	Shoulder Rope Pulleys	20,21
Babinski Hammer	73	Gym Mat	46	Shoulder Tube Pulley	20
Back Wrap	59	Hand Pump	48	Sit & Reach Unit	72
Back-Leg-Chest Dynamometer	67	Handmaster Plus	27	Sitting Cushions	44
Balance Boards	50,51	Hollow Foam Roller	41	Smedley Hand Dynamometer	67
Balance Pads	51	Hot/Cold Packs	61	Soft Handles	22
Balance Trainers	51	Hydraulic Hand Dynamometer	65	SoftMeds	55
Ball Accessories	48	Hydraulic Hand Evaluation Kits	64	Soft-Over-Ball	43
Ball Stabiliser	48	Incline Boards	39	Spiky Foam Roller	41
Ball Stacker	48	Inclinometer	69	Sport Pack	59
Band	6-9,15	Kettlebells	57	Squeeze Balls & Eggs	30
Band & Tubing Scissors	22	Knee Rope Pulley	20	Squeeze Dynamometer	67
Band Latex-Free	11	Knee/Elbow Wrap	59	Stepping Stones	42
Bar	18	Loop	13,17	Superloop	17
Blood Pressure Monitor	63	Massage Ball, Roll & Ring	45	Swim Noodle	54
Body Fat Scale	72	Massage Stick	41	Taylor Hammer	73
Body Twister	52	Massage Tools	45	TENS Device	74
Buck Hammer	73	Measurement Tape	69	Theraflex Putty	31
Comfort Putty	31	Medicine Ball	57	Thermometers	62
Constant Force X-Tend	25	Monofilament	72	Thumb Lock Wrist Weight	56
Cuff-Ring Tube Loop	18	Multi-Grip Exerciser	19	Torex Pro Ankle Boot	60
Curamotion	36	Needle Pump	48	Trampolines	40
Dejerine Hammer	73	Neurological Instruments	73	Triggerpoint Rollers	41
Digi-Extend	26	O-Ring Tube Loop	18	Troemner Hammer	73
Digi-Flex	26	Oximeters	63	Tube	12,16
Digital Hand Dynamometers	65	Oxycycles	34-35	Universal Assist	22
Discriminators	73	Paraffin Bath	58	Varigrip Sport	24
Door Anchor	22	Peanut Ball	49	Varigrip Therapy	24
Dumbbells	55	PearlPacks	59	Varigrip Uno	25
Dynadome	47	Pedal Exercisers	37	Wall Mount	22
Ekilibres	52	Pilates Foam Roller	43	Wartenberg Pinwheel	73
Electric Compressor Pump	48	Pilates Ring	43	Wheeleo	33
Ergonomic Foam Roller	41	Pinch Gauge	66	Wide Loop	13
Eye Mask	59	Power-Webs	28-29	Work-Out Station	53
Fabric Elastic Band	47	Prohands Series	26	Wrist & Ankle Weight	56
Faster Blaster Pump	48	Pumps	48	Wrist Exerciser	32
Fat Caliper	68	Push-Pull Dynamometer	66	Yoga Block	42
Figure 8-Ring Tube Loop	18	Puttys	31	Yoga Mat	42
Fit Stretch	38	Ring Tubes	18	Yoga Strap	42

TERMS & CONDITIONS ///

1. Terms: Prepayment, unless otherwise specified. Buyer should bear all bank charges.
2. Minimum purchase value: An administration fee of €10,00 will be added on all orders below the minimum purchase value of €250,00.
3. Freight Terms: Prices are ex-works Tisselt, Belgium. All shipping/delivery costs are at charge of the buyer, unless otherwise specified. Visible damage to goods delivered, not duly noted on the consignment note of the carrier, cannot be compensated in any form.
4. Prices: All prices for merchandise purchased are based on seller's price list in effect on the time of the order. Prices can be subject to change without notice.
5. Claims: Unless otherwise provided for herein, claims of any kind in goods sold by seller shall be made by buyer to seller, in writing, within a period of seven days after such goods are delivered to buyer.
6. Credit and Exchange: Written authorisation must be obtained before returning merchandise for credit or exchange! Prior to return, buyer shall furnish seller an itemized inventory of the products to be returned, and any other information seller may request, such as product lot number, manufacturing date, expiration date, ... In no event shall seller accept products for return more than thirty (30) days after delivery by seller to buyer and products which are damaged, used or not in their original and unopened packaging. Seller will charge a handling fee of ten percent (10%) of the returned goods.
7. Finance Charge: A finance charge of 1% per month (12 % per annum) will be added to all amounts past due. A penalty of 15% on value, with a minimum of €40,00, for late payment will be added. Terms of payment are from date seller's invoice and not contingent upon delivery.
8. Governing Law: All disputes arising out of all sales by seller to buyer shall be governed by the laws of Belgium and place of jurisdiction is Antwerp, section Mechelen.
9. Warranties: Seller warrants that it has good title to the goods purchased by buyer and that its goods will conform to buyer's specifications, if any, except as otherwise provided herein. All other warranties, whether express or implied, and including, but not limited to, the warranty of merchantability and fitness for a particular purpose are excluded. In all events, buyer's sole and exclusive remedy shall be limited to the amount of the purchase price of any product on which a claim is made as its exclusive remedy for breach of the warranty expressed herein and for negligence and otherwise. Seller reserves the right to repair or replace warranted products in its discretion. Under no circumstances shall seller be liable to buyer or any other person for consequential or incidental damages.
10. Property: Seller has the exclusive property on the delivered products until buyer has paid the full amount due. This also implies that if buyer does not pay in time, seller has the right to, without prenotice, claim back the delivered goods.




Your dependable partner **FERROX**
for dependable products & service.



The choice of physiotherapists for 50 years,
with know-how and expertise at your service.

www.ferrox.it
info@ferrox.it • +39 (0)4 38 77 70 99

Our
Own
Brands



Master
Distributor
For



WWW.MOVES-YOU.COM
WWW.MOVES-BAND.COM
WWW.MAMBOMAX.COM

Visit our websites for:
- Product information
- Exercise charts and videos
- Shop locator

MVS In Motion - Westdijk 150, 2830 Tisselt (Willebroek) - Belgium