

CATALOGUE 2021 2022

We are proud to present to you our new catalogue 2021-2022!

Since 1996, our company has been known as a worldwide manufacturer and master distributor of a broad range of physiotherapy, rehabilitation and sports & fitness products. Our products reach the end-users through a global network of dealers and distributors.

Our own branded MoVeS and Mambo Max products, as well as other well-known brands in our market, are all presented in this catalogue.

On page 7, you can find our new products: our MoVeS Kinesiology Tape, Mambo Max Balance Pad Pro, and our updated MoVeS PearlPack Sport

All our medical devices, marked with the MD symbol, comply with all the requirements of the EU's Medical Devices Regulation (MDR), which came into force May 26, 2021.

For all your inquiries for Private Label products, please contact us. We will be happy to assist you working out your private label projects.

Unfortunately, Covid-19 is still around, which has made 2021 a challenging year as we are currently experiencing the economic impact of this pandemic. Shortages of raw materials & containers as well as available vessel spaces cause high product price and transport price increases, as well as many and long delays on deliveries. We have adapted our systems to this present-day reality to minimise stock shortages.

This year is also special in the company's history; I herewith take the opportunity to introduce you to the new owners. It is with great confidence that I hand over the company to Evelyn and Tom, who will continue in the same spirit and with enormous enthusiasm. I wish them lots of success!

Finally, I would of course like to thank YOU for your confidence in the company and its products! Stay safe, stay healthy!

Sincerely, Sandra Vrancaert Former CEO/Owner



Sandra Vrancaert Former CEO/Owner sandra.vrancaert@mvs-in-motion.com



Evelyn Deceuninck CEO/Owner



Tom Vanhecke CEO/Owner tom.vanhecke@mvs-in-motion.com



Tom Lambrechts Business Development/Sales/QA



Sales tom.benoy@mvs-in-motion.com



Monique Breuer **Customer Service**



Patrick Verlinden Administration



Natascha Jooris **Customer Service**



Julie Persyn Marketing



Logistics logistics@mvs-in-motion.com



Patrick Steenackers Warehousing



Cem Iserio Warehousing



















marketing@mvs-in-motion.com

TABLE OF CONTENT ///

OUR OWN BRANDS ONLINE	6	3. UPPER & LOWER BODY	3
LEGEND NEW PRODUCTS	6 7	Fit Stretch	3
THE ENVIRONMENT	7	Saehan Cycle	3
		Wheeleo	3
4 DECICENT	4.4	MoVeS Oxycycle 1 - Active Magnetic	3
1. RESISTIVE	11	MoVeS Oxycycle 2 - Passive	3
MoVeS Band	11	MoVeS Oxycycle 3 - Active & Passive	3
MoVeS Latex-Free Band	13	MoVeS Oxycycle Pedals	3
MoVeS Tube	14	Curamotion	3
MoVeS Wide Loop	15	MoVeS Pedal Exerciser	3
MoVeS Loop	15	MoVeS Foldable Pedal Exerciser	3
MoVeS F!T Band	17	MoVeS Standard Pedal Exerciser Premium Chrome	3
MoVeS F!T Tube	18	MoVeS Standard Pedal Exerciser K/D	3
MoVeS F!T Loop	19	MoVeS Wooden Incline Board	4
MoVeS F!T Superloop	19	MoVeS Wooden Incline Board - Ladder	4
MoVeS Bar	20	MoVeS Wooden Incline Board - Fix	4
MoVeS Ring Tube	20	Plastic Incline Board	4
Multi-Grip Exerciser	21		
MoVeS Ankleciser	21		
MoVeS Shoulder Rope Pulley - Metal Bracket	22	4. WELLNESS	4
MoVeS Shoulder Tube Pulley	22	Wheelshor Torr	
MoVeS Knee Rope Pulley	22	Kinesiology Tape	4
MoVeS Shoulder Rope Pulley - Door Anchor	23	Yoga	4
Notch Handles	23	Stepping Stones	4
Value-Pack	23	Triggerpoint Rollers	4
	_	Spiky Foam Roller	4
Wooden Handles	23	Hollow Foam Roller	4
Hand Support	23	Ergonomic Foam Roller	4
MoVeS Accessories	24	3-in-1 Foam roller	4
		Massage Stick	4
2. HAND	25	Sitting Cushions	4
Z. HAND	23	Pilates	4
MoVeS Flex-Ion	25	Soft-Over-Ball	4
Digi-Flex	26	Foam Roller	4
Digi-Extend	26	Ring	4
Prohands Series	26	Massage Tools	4
Handmaster Plus	27	Massage Ball, Roll & Ring	4
Constant Force X-Tend	27		
Varigrip Therapy	28		
Varigrip Sport	28	5. FITNESS	4
Varigrip Uno	28		
Power-Web Combo	29	Gym Mat	4
Power-Web Flex-Grip	29	Ball Accessories	4
Power-Web Senior	30	Aerobic Step	4
Power-Web Junior	30	Fabric Elastic Band	4
MoVeS Squeeze Balls & Eggs	31	Dynadome	4
MoVeS Wrist Exerciser	31	Balance Board Freeman	5
Adjustable Hand Grip	32	Balance Board	5
Foam Hand Grip	32	Wooden Balance Boards	5
MoVeS Comfort Putty	33	AB Gym Ball	5
Theraflex Putty	33	AB Peanut Ball	5
,		Body Twister	5
		Ekilibre	5
		Wooden Work-Out Station	5
		Balance Pad Pro	5
		Balance Pad	5
		Balance Trainer	5
		Swim Noodle	5

FAIRS	54	Hydraulic hand dynamometer	6
		Digital hand dynamometer with Data Transfer Software-	
C CTRENCTU	55	DHD-1	6
6. STRENGTH	33	Digital hand dynamometer with G-STAR™ Software -	
Coated Dumbbell	55	DHD-3	6
SoftMeds	55	Smedley Hand Dynamometer	6
Wrist & Ankle Weight	56	Push-Pull Dynamometer	6
Thumb Lock Wrist Weight	56	Body Fat Scale	6
Medicine Ball	57	Pinch Gauge	6
Kettlebell	57	Squeeze Dynamometer	6
		Back-Leg-Chest Dynamometer	6
		Fat Caliper	6
7. THERMAL	58	Measurement Tape	6
		Tuning Fork Inclinometer	69
Para Pro - Paraffin Bath	58	Scoliometer	69
PearlPack	59	Goniometers	69 70
PearlPack Knee/Elbow	59		-
PearlPack Eye Mask	59	Baseline	7
PearlPack Back Wrap	59	Lafayette	7
PearlPack Sports Pack + strap	59	Goniometer Kit	7
Roll-On Pack	60	Plastic Goniometer	7
Torex Pro Ankle Boot	60	Stainless Steel Goniometer	7
Hot/Cold Pack	61	Finger Goniometer	7
		Monofilament	7
		Discriminator	7
8. EVALUATION	62	Neurological Reflex Instruments	7
lufusus d The sure street	63	Sit & Reach Unit	7
Infrared Thermometer	62	TENS device	7
Ear & Forehead Thermometer	62		
Shotgun Thermometer	62	DDIVATE LADEI	-
Blood Pressure Monitor	63	PRIVATE LABEL AFTER-SALES	7! 7(
Fingertip Pulse Oximeter	63	INDEX	7
Adult Oximeter	63	TERMS & CONDITIONS	78
Junior Oximeter	63		
Hydraulic 3-niece hand evaluation kit	64		



Hydraulic 7-piece hand evaluation kit

OUR OWN BRANDS ///



MoVeS products are aimed towards physiotherapy and rehabilitation, predominantly resistive exercise. Resistive exercise products have various levels of resistance, which also makes them perfect for fitness. MoVeS products are generally MD certified.

MoVeS

Our Mambo Max products are the go-to-products for wellness, fitness and strength training at an affordable price tag. Combine various products for a complete work-out, without leaving your home.

ONLINE

www.moves-you.com www.moves-band.com

The virtual home of our MoVeS and MoVeS F!T product line. This B2C platform connects our distributors directly to the end users, with product information, a shop locator and numerous free downloadable exercises. The website has been upgraded and is mobile friendly.

www.mambomax.com

Completely revitalised! This is the virtual home of our Mambo Max product line. This B2C platform connects our products directly to the end users, with product information and numerous free downloadable exercises.

www.mvs-in-motion.com

MVS In Motion's B2B site, with all available products and digital catalogue. Ask for your login to download all high resolution pictures, product information, logos, videos and many more marketing tools.

Social Media in 🖸 f 💿

Here you will find our latests and useful information. Like or follow our pages for a higher brand awareness to create more demand. Our pages are continuously monitored to answer your questions as soon as possible.

We are active on LinkedIn and YouTube. Soob, you will also be able to find us on Facebook and Instagram


LEGEND

In order to keep all the information in this catalogue as clear and easy accessible as possible, we use a system of icons that are shown next to the product where applicable.



Contains multilingual instructions.









Plastic free



Packed in a shrinkfoil



Product Diameter



Comes with a pump.

Packaging contains


Downloadable exercises

available on our websites.



Product Height


Product Width




printed exercises.



Maximum User Weight







European Conformity (Health, safety & environment)

NEW PRODUCTS







THE ENVIRONMENT



We, at MVS In Motion, continue our strive to reduce our carbon footprint on the planet. Because after all, there is no planet B!

By making exercises available on our websites, not only do we reduce our ecological footprint by cutting down on the number of unused prints, it also allows us to keep costs low and offer you better value-for-money.





We have also been making continuous efforts over the past years to replace our plastic packaging with the more eco-friendly and recyclable cardboard packaging. And we will continue these efforts in the future!





LEADER IN TECHNOLOGY



One-Stop-Shop for all your bands and tubes

MVS In Motion is your one-stop-shop for all resistive exercise products. With MoVeS we offer both original colour and F!T colour sequence. In addition, we offer a wide range of Latex-Free resistive exercise bands in original colours. Our tubes and loops are also available in the original and F!T colour sequence.





The latex bands, tubes and (super)loops are made of high quality latex from huge rubber tree plantations in Malaysia.

All products have an extraordinary elasticity and linearity. The Snap-Stop® feature on our latex bands, dipped production process for our tubes and the multi-layered structure of our (super)loops makes them your ideal partner for your physiotherapy trainings or intensive work-outs. The different resistances allow you to use MoVeS for rehabilitation as well as for fitness.

MoVeS Latex-Free is made from TPE, to the same high quality standards as the regular MoVeS bands. Never before were the characteristics from Latex-Free bands so close to those from their latex counterpart.

MoVeS Band: 8 resistances (Tan - Gold) & 5 lengths (1,5m - 2,5m - 5,5m - 22,5m - 45,5m) MoVeS F!T Band: 6 resistances (Peach - Grey) & 4 lengths (1,5m - 2,5m - 5,5m - 45,5m) MoVeS Latex-Free Band: 5 resistances (Yellow - Black) & 3 lengths (1,5m - 5,5m - 45,5m)



Original Colours


Extra Light
Light
Medium
Heavy
Extra Heavy
Special Heavy
Super Heavy
Ultra Heavy





F!T Colours

Peach	Light
Orange	Medium
Lime Green	Heavy
Blueberry	Extra Heavy
Plum	Special Heavy
Grey	Super Heavy



- Packaging with two tone layout brings more variety
- MoVeS brings premium quality to an affordable price range.





cb po m ex ex

MoVeS Band MD < €



MoVeS Bands are premium quality, low powder (in fact, less than 2mg/dm², so we're officially powder-free) latex resistive exercise bands. The unique Snap-Stop® fabrication process makes this product extremely tear resistant. Using the MoVeS Bands for resistive exercises provides positive as well as negative forces on the muscles and improves the range of motion, the strength and the cooperation of the muscle groups. The 8 different resistance levels are determined by the thickness of the band. MoVes bands are available in 5 lengths and come with a full colour manual. The band has a width of +/- 127mm.

Dispenser 25 bands of 1,5m in a nicely designed OTC dispenser box. Each band comes in a separate polybag.


- 1,5m 10 pre-cut bands of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header and barcode.
- 2,5m 10 pre-cut pieces of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header and barcode.
- 5,5m One small roll of our premium band.
- 22,5m One medium roll of our premium band.

Dispenser 1,5m

45,5m One large roll of our premium band.



TOTAL KG / ELONGATION %








PLASTIC

				MOVES	BAND		
COLOUR	RESISTANCE	1,5m Dispenser 25pcs	1,5m 10-pack	2,5m 10-pack	5,5m Box Small Roll	22,5m Box Medium Roll	45,5m Box Large Roll
Tan	Extra Light	n/a	n/a	n/a	01-100501	n/a	01-104501
Yellow	Light	01-100102	01-101502	01-100202	01-100502	01-102202	01-104502
Red	Medium	01-100103	01-101503	01-100203	01-100503	01-102203	01-104503
Green	Heavy	01-100104	01-101504	01-100204	01-100504	01-102204	01-104504
Blue	Extra Heavy	01-100105	01-101505	01-100205	01-100505	01-102205	01-104505
Black	Special Heavy	01-100106	01-101506	01-100206	01-100506	01-102206	01-104506
Silver	Super Heavy	n/a	01-101507	01-100207	01-100507	01-102207	n/a
Gold	Ultra Heavy	n/a	n/a	n/a	01-100508	01-102208	n/a
	5 resistances w to black)	n/a	01-101511 (5-pack)	01-100211 (5-pack)	n/a	n/a	30-St




150%

200%

9,0 8,0 8,0 7,0 5,0 1,0 50% 100% 150% 200% 250% 300% Percentage Elongation

Special Heavy



Other alternative:

Very high resistance

MoVeS Bands can be stretched up to 600% elongation! MoVeS Bands are therefore your perfect partner for a rehabilitation programme or a full body work-out. Even after an effective training with amazing results, MoVeS Bands allow you to use the same band again. When you prestretch your band, you will be able to achieve higher forces without having to worry about the 'elastic stop', the point where the band feels like a regular string. Most competitor's bands are not to be stretched more than 3 times their resting length, while a MoVeS Band achieves easily the double, maintaining its amazing elasticity.

These specifications apply to all MoVeS resistances. The high-quality material enables you to have a much wider array of resistances and possibilities compared to the competitors. At the start of your exercise, the band already gives some resistance, which allows you to achieve much better results than when the band is at its resting length.

Note

- 50% elongation = 1,5 x lay-flat length
- 100% elongation = 2 x lay-flat length
- 150% elongation = 2,5 x lay-flat length





Research and development is high on our list for MoVeS products. As an example, these continuous efforts have brought us the Snap-Stop feature for the regular MoVeS Band.


Resistance - MoVeS Latex-Free Band vs. Market Leader -


Now, we have improved our Latex-Free Bands. The MoVeS Latex-Free Band has been optimised by using an innovative mix of raw material. Following customer demand, the grip has been improved.

Resistance levels are similar to the 'old' latex-free band.











Best combination of resistance / linearity

Market Leader: Non-linear behaviour Very low resistance

We are confident that MoVeS Latex-Free Bands are one of the best latex-free resistive exercise bands on the market. The elasticity and linearity are simply amazing. Finally people who don't support latex can exercise with the same high quality specifications. MoVeS Latex-Free bands start from high quality thermoplastic elastomer (TPE), so it's non-toxic and completely PVC and phthalate free.

The linearity from how the strength per elongation evolves is one of the most important specifications as it offers the possibility to apply a very extended range of exercises without risking an injury due to strange unexpected product behaviour. Latex-free alternatives will reach their 'elastic stop' more quickly; this is the point where the band feels like a regular string.

We notice a growing demand and trend for latex free products. For instance, more and more nospitals and centers aim to be completely latex free to avoid possible allergic reactions for their patients or staff. After extensive research and development, we were able to launch our MoVeS Latex-Free Bands. A product which we are extremely proud of.

Thanks to our technological leadership, we achieve a Latex-Free exercise band with characteristics that were never before so close to its latex counterpart. With an exceptional 'resistance vs. elongation' linearity, we are confident our MoVeS Latex-Free Band is the best product to fulfil your needs.

The 5 different resistance levels are determined by the thickness of the band. These bands are available in 3 lengths and come with a full colour manual. The band has a width of +/- 100mm.

1,5m 10 pre-cut bands of the same resistance or 5 pieces of different resistance (yellow, red, green, blue and black). Each band comes in a separate polybag with header.

MOVES LATEX-FREE BAND

5,5m Box

Small Roll

01-120502

01-120503

01-120504

01-120505

01-120506

n/a

45,5m Box

Large Roll

01-124502

01-124503

01-124504

01-124505

01-124506

n/a

1,5m

10-pack

01-121502

01-121503

01-121504

01-121505

01-121506

01-121511

(5-pack)

5,5m One small roll of our latex free band.

45,5m One large roll of our latex free band.

COLOUD		TOTAL KG / ELONGATION %						
COLOUR	50%	100%	150%	200%	250%	300%		
Yellow	1,3	1,7	1,9	2,3	2,8	3,2		
Red	1,4	1,9	2,2	2,6	3,0	3,5		
Green	2,0	2,5	3,0	3,4	4,0	4,5		
Blue	2,3	3,0	3,4	3,9	4,7	5,5		
Black	2,4	3,4	3,9	4,3	5,1	6,1		



RESISTANCE

Light

Medium Heavy

Extra Heavy

Special Heavy

Set of 5 resistances

(vellow to black)

COLOUR

Yellow

Red

Blue

Black





cb ex

MoVeS Tube MD < €

The MoVeS Tube is a high quality tube created by an intensive dipping process. The dipping process guarantees a multi-layered tube, which ensures a strong and safe training product. Extruded tubes are very susceptible to damage and can easily snap, this is not the case for our multi-layered dipped tubes. In general, dipped tubes are even more durable than bands, especially when

MoVeS Tubes are available as a ready to use product - 120cm length with handles straight out of a polybag - and in rolls - 7,5m or 30m packed in a nicely designed box.

The table below shows the exceptional strength to elongation linearity from our tubes.

COLOUR -	TOTAL KG / ELONGATION %					
COLOUR	50%	100%	150%	200%	250%	300%
Yellow	1,0	1,5	1,9	2,3	2,9	3,6
Red	1,4	2,0	2,5	2,8	3,4	4,0
Green	2,4	3,4	4,0	4,9	5,5	6,5
Blue	2,9	3,9	4,7	5,6	5,9	7,4
Black	3,4	4,4	5,4	6,2	7,1	8,1
Silver	4,0	5,4	6,5	7,7	8,1	9,5

Combine a roll with our Safety Handles for superior grip (p22).

MoVeS Tubes with Handles are very useful for physiotherapy

and muscle training. These high quality and very resistant tubes are created to train all muscles in your body and improve body flexibility. The lighter resistances are perfect for physiotherapy or a light work-out, while the highest resistances will even convince the more trained athlete.



Box of 30m contains maximal 3 separate pieces. Minimal length of 1 piece is 4m



- 111.	1
	Polybag 120cm with Handles
	MoVes Value Value

COLOUR RECICTANCE		MOVES TUBE					
COLOUR	RESISTANCE	120cm with Handles	7,5m roll	30m roll			
Yellow	Light	01-200102	01-200702	01-203002			
Red	Medium	01-200103	01-200703	01-203003			
Green	Heavy	01-200104	01-200704	01-203004			
Blue	Extra Heavy	01-200105	01-200705	01-203005			
Black	Special Heavy	01-200106	01-200706	01-203006			
Silver	Super Heavy	01-200107	01-200707	01-203007			



MoVeS Wide Loop MD (€ | | | |



We offer wide loops (7,5cm) next to our regular loops (2,5cm) with the same resistances. Some users may prefer the wider loops for spreading the forces more evenly on their body parts, while others prefer the thicker narrow loops which won't roll-up while being stretched.

As your supplier for bands, tubes and loops, these wide loops fit perfectly in our range of resistive products to satisfy all users. Wide loops are also multi-layered and available in 6 resistances: from Extra Light (Tan) to Special Heavy (Black).

MoVeS Wide Loops:

High quality, cost effective and easy to use!

Lay-flat-length: 30cm

Width: 7,5cm

One unit contains 10 individually packed pcs of the same colour/resistance

COLOUR	RESISTANCE	MOVES WIDE LOOP	UNIT
Tan	Extra Light	01-100071	10 pcs
Yellow	Light	01-100072	10 pcs
Red	Medium	01-100073	10 pcs
Green	Heavy	01-100074	10 pcs
Blue	Extra Heavy	01-100075	10 pcs
Black	Special Heavy	01-100076	10 pcs







MoVeS



MoVeS Loop MD (€ | | | |

MoVeS Loops are extremely portable and unbelievably versatile products. They are very effective instruments for stretching and strengthening all muscle groups. The MoVeS Loop Set really is a "Gym on-the-go".

Our loops are also multi-layered which results in a more durable product. These thicker loops won't rollup during exercising, avoiding the uncomfortable feeling from thinner alternatives.

They are available in 6 resistances: from Extra Light (Tan) to Special Heavy (Black).

The loop set combines our most popular colours in a compact mesh bag. This affordable set is very useful to bring along wherever you

Multi-Layered MoVeS Loops:

High quality, cost effective and easy to use! Lay-flat-length: 30cm

Width: 2,5cm

One polybag contains 10 loose pcs of the same colour/resistance



_	TOTAL KG / ELONGATION %							
COLOUR	50%	100%	150%	200%	250%	300%		
Tan	2.2	3.2	3.8	4.6	5.4	6.4		
Yellow	2.8	4,0	4,8	5,8	6,8	8,0		
Red	3,4	4,8	5,8	6,8	8,0	9,2		
Green	4,4	6,4	7,8	9,2	10,8	12,4		
Blue	5,0	7,2	9,0	10,8	12,6	14,4		
Black	6,1	8,9	11,0	13,0	15,0	17,3		

14

Box with Rolls

7,5m / 30m

Made in Malaysia –









The production of MoVeS latex products start deep in the humid forests of Malaysia, where the best and strongest trees provide the natural latex rubber that will be processed from a milky fluid to your preferred resistive exercise band. When a tree is 5 - 6 years old, harvesting can begin: incisions are being made orthogonally in the bark of the tree, just deep enough to tap the latex vessels without harming the tree's growth. The fluid is collected in small buckets, a process known as rubber tapping.

Did you know that rubber trees combat 🌡 greenhouse gases and thereby slowing down global warming? Additionally, latex harvesting helps prevent deforestation.



"This is why we are confident when using the words 'premium quality' when referring to MoVeS."

Only the best quality latex is being selected for MoVeS Bands. Production starts under the constant supervision of a team of scientists, ensuring an immaculate product quality.

Latex sheets are produced, powdered, washed, dried and cut to the exact width. The amount of powder used is below the 2mg/dm2 limit, which makes the MoVeS Bands officially powder-free. Each batch of latex sheeting that is produced, is subject to the meticulous inspection of our scientists by a series of laboratory experiments: colour-coding, dryness, tensile strength, tearing resistance, snapping resistance, ...

During the rolling process, a 100% visual inspection is being conducted on both sides of the band. Indeed, before leaving the factory, every MoVeS product has been completely analysed, checked and examined to ensure premium product quality.







MoVeS F!T, for athletes and home use.

MoVeS F!T Bands are low powder (in fact, less than 2mg/dm², so we're officially powder-free) latex resistive exercise band with the same premium quality as the original MoVeS Bands. The difference is purely the alternative, more youthful colour sequence. The band comes in 6 different resistance levels which are determined by the thickness of the band, is available in 4 different lengths and comes with a full colour manual.

TOTAL KG / ELONGATION %						Compares to
COLOUR	50%	100%	150%	200%	250%	STANDARD
Peach	1,0	1,3	1,5	1,9	2,3	Yellow
Orange	1,5	2,0	2,5	3,2	3,7	Red
Lime Green	1,6	2,5	3,0	3,5	4,1	Green
Blueberry	2,2	3,0	3,6	4,4	5,2	Blue
Plum	2,5	3,4	4,0	4,9	6,1	Black
Grey	4,1	5,8	7,2	7,9	9,5	Silver

Our high quality 1,5m and 2,5m MoVeS F!T Bands are packed in retail packs. 10 pre-cut pieces of the same resistance or 5 pieces of different resistance that come in separate polybags with header for easy dispensing. Full-colour instruction manual is included. These 1,5m and 2,5m bands are easy to use and cutting is not required. The available colours are peach, orange, lime green, blueberry, plum and grey.













The MoVeS F!T Tube is a high quality tube created by an intensive dipping process. The dipping process guarantees a multi-layered tube, which ensures a strong and safe training product. Extruded tubes are very susceptible to damage and can easily snap, this is not the case for our multi-layered dipped tubes. In general, dipped tubes are even more durable than bands, especially when used on a coarse surface.

MoVeS F!T Tubes are available as a ready to use product in 30m rolls, packed in a great dispenser box. Safety Handles offer an unlimited amount of exercises and superior grip.

MoVeS F!T Tubes are available in 6 different colours. The table below shows the resistance and to which standard colour it relates.

COLOUR -	٦	TOTAL KG	Compares to			
COLOUR	50%	100%	150%	200%	250%	STANDARD
Peach	1,0	1,5	1,9	2,3	2,9	Yellow
Orange	1,4	2,0	2,5	2,8	3,4	Red
Lime Green	2,4	3,4	4,0	4,9	5,5	Green
Blueberry	2,9	3,9		5,6	5,9	Blue
Plum	3,4	4,4	5,4	6,2	7,1	Black
Grey	4,0	5,4	6,5	7,7	8,1	Silver



Box of 30m contains maximal 3 separate pieces.
Minimal length of 1 piece is 4m.



COLOUR	RESISTANCE	MOVES F!T TUBE
COLOUR	RESISTANCE	30m roll
Peach	Light	01-213002
Orange	Medium	01-213003
Lime Green	Heavy	01-213004
Blueberry	Extra Heavy	01-213005
Plum	Special Heavy	01-213006
Grev	Super Heavy	01-213007





MoVeS F!T Loop MD (MoVeS F!\(\sigma\)

Our high quality and cost effective loops are also available in F!T colours. MoVeS F!T Loops are easy to use and very versatile. You can stretch and strengthen all muscle groups.

Our loops are multi-layered, this results in a more durable product. Our thicker loops won't roll-up during exercising, avoiding the uncomfortable feeling from thinner alternatives.

Our loop set combines our most popular colours in a compact mesh bag. This affordable set is very useful to bring along wherever you go.

Width: 2,5cm

MoVeS F!T Loops: High quality, cost effective and easy to use! Lay-flat-length: 30cm

One Polybag contains 10 loose pcs

Multi-Layered

COLOUR	RESISTANCE	MOVES F!T LOOP	UNIT
Peach	Light	01-110002	10 pcs
Orange	Medium	01-110003	10 pcs
Lime Green	Heavy	01-110004	10 pcs
Blueberry	Extra Heavy	01-110005	10 pcs
	4 resistances perry) in mesh bag	01-110011	Set / 4 pcs



J. 185				TI F					
ex ex			MoVe	S F!T	Supe	rloop	MD (€		
COLOUD	CIZE IN MANA		TOTAL KO	G / ELONG	SATION %			MoVeS F!	
COLOUR	SIZE IN MM	100%	200%	300%	400%	500%			
Peach	1040 x 15 x 3,0	7,3	10,0	14,4	19,1	29,2			
Orange	1040 x 22 x 4,0	13,0	20,7	27,1	35,0	66,3			
Lime Green	1040 x 30 x 4,4	22,5	33,0	44,0	56,2	97,0			
Blueberry	1040 x 40 x 4,4	25,6	40,0	53,0	76,0	132,0		A PART OF THE PART	
Plum	1040 x 53 x 4 4	38.0	56.0	79.0	105.0	167.0			

With the MoVeS F!T Superloop you can easily stretch and strengthen different muscle groups. It can be used for rehabilitation, fitness and strength training. The lay-flat-length from the MoVeS F!T Superloop is 104cm long (circumference is 208cm) and has an incredible elasticity.

The MoVeS F!T Superloops are multi-layered, this results in a more durable product. They are available in 5 different resistances: Peach, Orange, Lime Green, Blueberry and Plum.

COLOUR	RESISTANCE	MOVES F!T SUPERLOOP	UNIT
Peach	Light	01-110402	1 pc
Orange	Medium	01-110403	1 pc
Lime Green	Heavy	01-110404	1 pc
Blueberry	Extra Heavy	01-110405	1 pc
Plum	Special Heavy	01-110406	1 pc









The MoVeS Bar is the ideal product to treat and resolve golfer's and tennis

This effective product is developed to focus on hand, arm and shoulder problems, this product allows flexion and extension exercises, adduction and abduction exercises and even soft tissue mobilisation, far more than the twisting and bending you would imagine.

The MoVeS Bar is available in 6 colours / resistances: extra light, light, medium,

Size: 31cm x 4,5cm

KG-FORCE* 2,5 - 3,5

3,5 - 4,5

4,7 - 5,7 6,7 - 7,7

9,8 - 10,8

12,0 - 14,0

Packaging: polybag with full-colour header

Material: dried natural rubber

1 size Ø; 6 Colours / Resistances

RESISTANCE

Extra Light

Light

Medium

Heavy

Extra Heavy Special Heavy

	Tenr
MOVES BAR	
01-600101	
01-600102	
01-600103	
01-600104	
01-600105	
01-600106	

*force applied when bending to a U-shape


COLOUR

Tan

Yellow

Green

Blue



heavy, extra heavy and special heavy.



For intensive training of isolated muscles.

Efficient for all muscles of upper- and

MoVeS 8-Ring Tube offers lots of

possibilities for short and long pull

lower body. The centre of the


MoVeS Ring Tube MD < €

exercises.



is a versatile exercise tool to tone and strengthen all the muscles of the lower body. With velcro closure that fits

The MoVeS Cuff-Ring Tube most sizes.



exercises. With the MoVeS O-Ring Tube you							
get an	intensive	training	of	each	muscular		
group.	١						
)						

	Annual Control			
COLOUR	RESISTANCE	O-RING	FIGURE 8-RING	CUFF-RING
		TUBE LOOP	TUBE LOOP	TUBE LOOP*
Yellow	Light	01-200202	01-200302	01-200402
Red	Medium	01-200203	01-200303	01-200403
Green	Heavy	01-200204	01-200304	01-200404
Blue	Extra Heavy	01-200205	01-200305	01-200405
Black	Special Heavy	01-200206	01-200306	01-200406



DIPPEded not extruded

* This product cannot be sold in Germany.



Multi-Grip Exerciser MD < €

The Multi-Grip™ Exerciser is designed to provide unlimited exercise options to users who want to do various high impact performance exercises. Each exerciser contains both large and small loops that are easy to grip. Use the small loops as an anchor or for finger exercises. Use the larger loops around the hands and feet to perform high impact upper and lower body exercises. Easily change resistance by moving 1 loop up or down the exerciser, or using a different resistance progressively exercise with 8 different colour-coded resistance levels.

The Multi-Grip™ Exerciser is made of an unscented cloth that can be washed and dried. It does not contain latex or a powder coating.



COLOUR	RESISTANCE	MULTI-GRIP EXERCISER
COLOUR	RESISTANCE	30 Yard Roll (27,4m)
Tan	XX-Light	01-703001
Yellow	X-Light	01-703002
Red	Light	01-703003
Green	Medium	01-703004
Blue	Heavy	01-703005
Black	X-Heavy	01-703006
Silver	XX- Heavy	01-703007
Gold	XXX-Heavy	01-703008

COLOUR	TOTAL KG / EL	ONGATION %
COLOUR -	100%	200%
Tan	0,9	2,3
Yellow	1,4	3,4
Red	1,8	4,5
Green	2,3	5,7
Blue	2,7	6,8
Black	3,2	7,9
Silver	4,5	11,3
Gold	5,9	14,7
	Yellow Red Green Blue Black Silver	Tan 0,9 Yellow 1,4 Red 1,8 Green 2,3 Blue 2,7 Black 3,2 Silver 4,5

MoVeS Ankleciser <€







	COLOUR RESISTANCE		MOVES	UNIT
	COLOUR RESISTANCE	ANKLECISER	UNII	
	Yellow	Light	01-500102	5 pcs
Red	Medium	01-500103	5 pcs	
	Green	Heavy	01-500104	5 pcs
	Blue	Extra Heavy	01-500105	5 pcs
	Set of 4	resistances	01-500111	4 pcs

The MoVeS Ankleciser is a great ankle exerciser, complying to all the rigorous quality standards MVS In Motion demands of its MoVeS products. It's as easy in use as it's effective. The feet slip comfortably into the latex free elastic band which consists of two stitched loops. Available in 4 different resistances, this progressive exercise system enables the user to choose the desired resistance level to perform the exercises.

Available in a set of 4 resistances, or in a 5-pack of the same resistance. Packed in a polybag with full-colour header and exercise chart.



03-UPPER & LOWER BODY

MoVeS Shoulder Rope Pulley - Metal Bracket MD (€ MoVeS





MoVeS








This MoVeS Shoulder Rope Pulley comes with a pulley, wooden handles, nylon cord, stopper system, and metal bracket.

Reference	Description	Unit
01-400102	MoVeS Shoulder Rope Pulley Metal Bracket Wooden Handles Stoppers	1 pc




MoVeS Shoulder Tube Pulley MD < €

Just like the MoVeS Shoulder Rope Pulley, the MoVeS Shoulder Tube Pulley helps increase range of shoulder motion. It can be used for prevention and rehabilitation of shoulder injuries.

The MoVeS Shoulder Tube Pulley has a light, medium or heavy

Each pack contains the selected resistance of tube, soft handles, pulley and door anchor.

CO	_OUR	RESISTANCE	SHOULDER TUBE PULLEY	UNIT
Ye	llow	Light	01-400202	1 pc
F	ed	Medium	01-400203	1 pc
Gr	een	Heavy	01-400204	1 pc





MoVeS Knee Rope Pulley MD (€





The MoVeS Knee Rope Pulley helps increase range of knee motion. This easy-to-use product can be mounted in a doorway to help you recover from or prevent knee injuries.


While sitting on a chair, put your foot into the ankle strap. A rope starts from the ankle strap, goes through a pulley and ends at a wooden handle designed for a 2-handed grip. Pull your foot carefully upwards by moving the handle towards your abdomen which will generate a rotational movement

The rope length can easily be adjusted at the handle while a print on the ankle strap helps you to secure your foot in the right position. The MoVeS Knee Rope Pulley comes in a full colour cardboard packaging with a pulley, wooden 2-handed handle, ankle strap, nylon cord and door anchor.

Reference	Description	Unit
01-400301	MoVeS Knee Rope Pulley Door Anchor 2-handed Wooden Handle Ankle Strap	1 pc

MoVeS Shoulder Rope Pulley - Door Anchor

The MoVeS Shoulder Rope Pulley is an economical, easy-to-use device that helps increase range of shoulder motion. It can be used for prevention and rehabilitation of shoulder injuries. This easy to use shoulder rope pulley can be mounted in a doorway.



This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and notch handles.

The unique notch handles allow the user to easily adjust the length of the rope through the handle while the revolutionary stopper system is a great aid for consistent and progressive exercising.



The VALUE PACK is packed in a polybag, with notch handles but without stoppers on the nylon cord.







This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and wooden handles.

Hand Support








This MoVeS Shoulder Rope Pulley comes with a pulley, door anchor, nylon cord, stopper system and detachable hand support.


This hand support allows patients without gripping power to use the shoulder pulley.

Reference	Description	Unit
01-400101	MoVeS Shoulder Rope Pulley Door Anchor Notch Handles Stoppers	1 pc
01-400104	MoVeS Shoulder Rope Pulley Door Anchor VALUE PACK (notch handles, polybag, no stoppers)	1 pc
01-400105	MoVeS Shoulder Rope Pulley Door Anchor Wooden Handles Stoppers	1 pc
01-400103	MoVeS Shoulder Rope Pulley Door Anchor Hand Support Stoppers	1 pc

MoVeS Accessories < €









Reference	Description	Unit
01-300001	MoVeS Wall Mount (Plug & Play)	1 pc
01-300101	MoVeS Safety Handles	pair
01-300102	MoVeS Soft Handles	pair
01-300103	MoVeS Door Anchor	1 pc
01-300104	MoVeS Universal Assist	1 pc
01-300106	MoVeS Band & Tubing scissors	1 pc







MoVeS Flex-Ion MD (€

M₀V_eS

The spring loaded buttons of the MoVeS Flex-Ion isolate each digit, so the strong fingers cannot compensate for the weaker ones. This results in a better, faster and easier to control rehab process. MoVeS Flex-Ion is available in 6 colour coded units with gradual resistance. This gradual resistance helps the therapist measure progress. The ergonomic base adds comfort and stability. Comes with a full colour and multilingual instruction manual.



COLOUR	RESISTANCE (1 digit all digits)	LEVEL	REFERENCE	UNIT
Tan	0,35kg 1,1kg	0	02-000101	1 pc
Yellow	0,7kg 2,3kg	1	02-000102	1 pc
Red	1,4kg 4,5kg	2	02-000103	1 pc
Green	2,3kg 7,3kg	3	02-000104	1 pc
Blue	3,2kg 10,4kg	4	02-000105	1 pc
Black	4,1kg 14,1kg	5	02-000106	1 pc
MoVeS Flex	k-Ion Set 6 Resistances 1 Displa	ny	02-000111	set
MoVeS Flex	k-Ion Empty Display Metal		02-000121	1 pc



Digi-Flex MD (€

Digi-Flex® develops isolated finger strength, flexibility and coordination. Suitable for rehabilitation after hand injury. The Digi-Flex® comes with an exercise leaflet.

COLOUR	RESISTANCE (1 digit all digits)	REFERENCE	UNIT
Yellow	0,7kg 2,3kg	02-010102	1 pc
Red	1,4kg 4,5kg	02-010103	1 pc
Green	2,3kg 7,3kg	02-010104	1 pc
Blue	3,2kg 10,4kg	02-010105	1 pc
Black	4,1kg 14,1kg	02-010106	1 pc
Digi-Flex Set 5	Resistances 1 Display	02-010111	set
Digi-Flex Empty	02-010121	1 pc	
	·	·	





Digi-Extend MD (€

The Digi-Extend* is a unique therapy device that can isolate the extension of all five fingers and each of the individual finger joints. The product strengthens the intrinsic and extrinsic muscle groups to maximise extensor tendon gliding, joint motion, coordination and muscle function in the hand and forearm.

The Digi-Extend* focuses on the extensor mechanism of one or all fingers as it targets the area of the hand that needs rehabilitation. Patients can easily change the resistance levels and vary their hand positions for specific exercises to help maximise the exercise goals and prevent injury. The Digi-Extend* can also be used for finger adduction, abduction, blocking and flexion exercises. It is compact in design, high-tech, easy and comfortable to use as well as effective in helping patients reach their therapy goals. Four colour coded band resistance levels allow each muscle or joint to work at a comfortable level and build as increased strength is achieved.

Reference	Description	Unit
02-010201	Digi-Extend Bands (Tan to Green) included	1 pc
02-010211	Digi-Extend Set 4 units 1 Display	set
02.040202	Digi-Extend Band Set 25 pcs	sot
02-010202	5x Tan, 5x Yellow, 5x Red, 5x Green, 5x Blue	set
02-010221	Digi-Extend Empty Display Metal	1 pc



Bai	nds	included:
-	5x	Tan
-	4x	Yellow
-	3x	Red
-	2x	Green



Prohands Series <€



VIA[™] is an entry level hand exerciser designed with young athletes in mind. It has a round, soft look, wide finger pistons for smoother function & improved finger performance.

The VIA[™] is available in two tensions.

RESISTANCE	RESISTANCE (per digit)	VIA	UNIT
Light	1,5 kg	02-060102	1 pc
Medium	2,5 kg	02-060103	1 pc



PRO™ provides serious athletes with the most advanced hand exerciser in the world. The PRO™ features a larger profile, superior steel springs that resist wear, an integrated ergonomic palm bar for comfort.

The $\text{PRO}^{\mbox{\tiny M}}$ is available in three tensions.

RESISTANCE	RESISTANCE	PRO	UNIT
RESISTANCE	(per digit)		
Medium	3,5kg	02-060303	1 pc
Heavy	4,5kg	02-060304	1 pc
Extra Heavy	5,5kg	02-060305	1 pc



is maximum strength, and balance to the hand, wrist and elbow.

Handmaster Plus

The Handmaster Plus[™] strengthens the 9 muscles that close the hand as well as the 9 muscles that open the hand. The result

 ϵ

3-pack set

3-pack set

3 balls for the price of 2!





Soft | Purple Ball - Yellow Cord | light exercise Medium | Red Ball - Black Cord | medium exercise Firm | Orange Ball - Grey Cord | strength, speed training







Reference	Description	Unit
02-050102	Handmaster Plus Soft Purple Ball - Yellow Cord	1 pc
02-050103	Handmaster Plus Medium Red Ball - Black Cord	1 pc
02-050104	Handmaster Plus Firm Orange Ball - Grey Cord	1 pc
02-050111	Handmaster Plus Set of 3 (one of each level)	Set



Constant Force X-Tend

CE

Unlike elastic-based exercisers that increase in resistance the further the unit is extended, Constant Force X-Tend uses retractable-based resistance, which keeps a constant level of resistance from start to finish.

Ergonomic in design, the Constant Force X-Tend works the Intrinsic and Extrinsic Extensors of the hand and can isolate MP/PIP/DIP joint motion.

It also allows each finger to exercise independently with smooth,

consistent resistance, thus improving range of motion, endurance, strength and dexterity more effectively.



- Exercise each finger or thumb independently with smooth, consistent resistance to the intrinsic and extrinsic extensor muscles throughout their full range of motion.
- Improve range of motion, strength, dexterity, endurance, circulation for each finger as well as the hand, wrist and forearm.
- Ergonomic design fits comfortably in hand, and includes ability to isolate pip & dip joints.

COLOUR RESISTANCE (per digi		REFERENCE	UNIT
Light	0,11 kg	02-100403	1 pc
Medium	0,18 kg	02-100404	1 pc
Heavy	0,27 kg	02-100405	1 pc



po m ex

Varigrip Therapy <€



The VariGrip Therapy has the ability to adjust tension in each individual finger, the resistance from each single finger knob can easily be adjusted. Resistance can be increased incrementally on the patient's need, which minimises overcompensation and strain. The VariGrip allows the patient to progress steadily between resistance levels.

5 different resistance levels are available.

Comparison Chart: VariGrip vs Digi-Flex® / Flex-Ion (kg per digit)										
VariGrip	X-L	ight	Lig	ht	Med	lium	He	avy	Х-Не	eavy
Digi-Flex® / Flex-Ion	0,7	1,4	1,4	2,3	2,3	3,2	3,2	4,1	4,1	n/a

Each Varigrip's highest resistance ends where the next VariGrip resistance starts. Any other hand exerciser increases from Light (1,4kg) to Medium (2,3kg). That's almost a 70% increase gap, which can cause an overload on the muscles, resulting in injuries. The seamless progress from one VariGrip to the other, supported by the padded finger and palm area, guarantees a safe recovery and better feel for those with sensitivity issues.


COLOUR	COLOUR RESISTANCE LEVEL (range per digit)		REFERENCE	UNIT
X-Light	0,7kg - 1,4kg	1	02-100102	1 pc
Light	1,4kg - 2,3kg	2	02-100103	1 pc
Medium	Medium 2,3kg - 3,2kg 3 Heavy 3,2kg - 4,1kg 4		02-100104	1 pc
Heavy			02-100105	1 pc
X-Heavy	4,1kg - 5,4kg	5	02-100106	1 pc
VariGrip Set	5 Resistances 1 Disp	02-100107	set	
VariGrip Emp	VariGrip Empty Display Plastic			





po m ex

Varigrip Sport <€



a callus builder, ideal for climbers, gymnasts, gamers, musicians and athletes of all levels

02-100304

Reference	Description	Resistance (kg)	Unit
02-100303	VariGrip Sport	2,3 - 4,1	1 pc
02-100304	VariGrip Sport Pro	4,1 - 5,4	1 pc

Resistance ad-

justable

per finger




02-100303

Soft touch

thumb button

Varigrip Uno <€





The only thumb/individual finger exerciser with adjustable resistance for both flexion and extension. An ergonomic design with stabilisation points offers extra support.

Gradually adjust resistance and better track progression based on the patients needs. Benefits are improved strength, dexterity, endurance and range of motion.

Light Spring 0,7kg - 2,3kg 2,3kg - 3,6kg **Heavy Spring**

The VariGrip Uno includes two easily interchangeable springs. The red spring allows a light resistance for early therapy or light training. Once a higher resistance is required the red spring can be swapped by the green spring, allowing a higher resistance.

Reference	Description	Unit
02-100204	VariGrip Uno	1 pc





Power-Web Combo MD (€





The #1 hand exercise web in the world!



Each Power-Web comes with a 39 exercise poster



COLOUR /	POWER-WEB COMBO	
Tan/Ultra Light	Red/Medium	02-020101
Yellow/Light	Green/Heavy	02-020102
Red/Medium	Blue/Super Heavy	02-020103
Green/Heavy	Black/Ultra Heavy	02-020104

The Power-Web® Combo features a unique advantage of two resistance levels in one, offering greater flexibility and value. It is identical in size and weight to the Power-Web® Senior with the same high quality construction and materials.



Power-Web Flex-Grip MD < €



COLOUR	RESISTANCE	POWER-WEB
COLOUR	RESISTANCE	FLEX-GRIP
Yellow	Light	02-021102
Red	Medium	02-021103
Green	Heavy	02-021104

The creators of the best-selling Power-Web[®], are proud to present the original Power-Web® Flex-Grip®. The Power-Web® Flex-Grip® is a unique exerciser that fits in your palm and provides the ability to work dexterity, extension and strengthening in your hands, fingers, and thumbs. Available in 3 different resistances, the Power-Web® Flex-Grip® has finger and thumb holes and has two sides: one side featuring a ball to squeeze and the other featuring buttons and a thumb pad to exercise fingers separately and increase dexterity and strength in your fingers and thumb.

Latex-free and easy to use!



\\\ cb m ex

Power-Web Senior MD (€




The Power-Web® Senior is 36cm in diameter and weighs 450g. Its unique patented design allows a variety of exercises to be performed comfortably and effectively. Constructed of high quality rubber with special agents added for durability and strength, the Power-Web® will last for years.

Now also available as latex-free!



COLOUR	RESISTANCE	POWER-WEB	POWER-WEB
COLOUR	RESISTANCE	SENIOR	LATEX-FREE SENIOR
Tan	Ultra Light	02-020201	02-022201
Yellow	Light	02-020202	02-022202
Red	Medium	02-020203	02-022203
Green	Heavy	02-020204	02-022204
Blue	Super Heavy	02-020205	02-022205
Black	Ultra Heavy	02-020206	02-022206

Power-Web Junior MD (€





COL	OUR	RESISTANCE	POWER-WEB
COL	OUR	RESISTANCE	JUNIOR
Yel	low	Light	02-020302
R	ed	Medium	02-020303
Gre	een	Heavy	02-020304



NoVeS Squeeze Balls & Eggs MD (€





5 colours. It strengthens hands, fingers, forearms and will work the muscles separately or in tandem. The MoVeS Squeeze Ball can be used for rehab as well as for sports training or stress relief. The MoVeS Squeeze Egg is ergonomically shaped to provide highly effective

exercise for the rehabilitation and strengthening of weakened or injured fingers, hands and wrists. Because it is small and easy to use, you can exercise anywhere, anytime.


Available colours/resistances: yellow, red, green, blue and black.



The MoVeS Squeeze Balls & Eggs can help to:

- Rehabilitate post-operative/post-traumatic hands and arms.
- Restore function in hands, arms and shoulders following stroke.
- Reduce pain and stiffness in arthritic fingers, hands and wrists.
- Relax tense, sore muscles and reduce stress.
- Strengthen aging hands and arms.

COLOUR	RESISTANCE	LEVEL	SQUEEZE BALL	SQUEEZE EGG
Yellow	Extra Soft	1	02-030102	02-030202
Red	Soft	2	02-030103	02-030203
Green	Medium	3	02-030104	02-030204
Blue	Firm	4	02-030105	02-030205
Black	Extra Firm	5	02-030106	02-030206










The MoVeS Wrist Exerciser is an adjustable wrist and forearm trainer which allows the user to exercise with progressive resistance. Its ergonomically designed forearm support keeps the wrist in the right place. Together with the soft handgrip, adjustable to 4 positions, it ensures maximum comfort while exercising your wrist and forearm.

Reference	Description	Unit
02-070201	MoVeS Wrist Exerciser	1 pc

Adjustable Hand Grip



The Mambo Max Adjustable Hand Grip comes in two different versions. Whether it is the 5 to 20 kg light version (orange) or the 10 to 40 kg heavy version (grey), you can always adjust the required





The Mambo Max Foam Hand Grip assures a firm training grip. The foam handles provide a firm grip for the entire workout session. Take it with you to the office or just use it at home. Wherever you want to use it, the Mambo Max Foam Hand Grip assures a tight

Gripping power = 20kg

grasp.

Reference	Description	Unit
02-080101	Mambo Max Foam Hand Grip	1 pair



Why Anti-Microbial Agents? _


Anti-microbial agents kill bacterium without additives and resist growth of odour-causing microbes. MoVeS Comfort






MoVeS Comfort Putty MD C€













Theraflex Putty MD (€



microbial agents

Theraflex[™] Putty strengthens weakened grasp through exercise and is used for rehabilitative exercise or stress relief. It is odourless, does not stain and will not leave any residue on the patient's hands. Theraflex™ Putty is nontoxic, greaseless, hypo-allergenic and latex-free. Theraflex[™] Putty is available in 5 colours/resistances.



COLOUR RESISTANCE		RESISTANCE LEVEL THERAFLEX PUTTY						
		LEVEL	57g	85g	113g	170g	454g	2,3kg
Tan	Extra Soft	0	02-040001	02-040101	02-040201	02-040301	02-040401	02-040501
Yellow	Soft	1	02-040002	02-040102	02-040202	02-040302	02-040402	02-040502
Red	Medium	2	02-040003	02-040103	02-040203	02-040303	02-040403	02-040503
Green	Firm	3	02-040004	02-040104	02-040204	02-040304	02-040404	02-040504
Blue	Extra Firm	4	02-040005	02-040105	02-040205	02-040305	02-040405	02-040505








The Mambo Max Fit Stretch is a fantastic, simple and portable tool to stretch your lower leg muscles. It is perfect for warm-ups and cooldowns. The use of a Mambo Max Fit Stretch helps increase your ankle, shin and calf flexibility. It is made of durable, moulded plastic and includes a non-skid foot pad.

Reference	Description	Unit
03-050101	Mambo Max Fit Stretch Single Leg	1 pc



Saehan Cycle (€ §SAEHAN®





The Clinical version from the Saehan Cycle has an even higher level of adjustability. The distance between pedal and crank

can be adjusted and the resistance knob indicates the level

sistance (0 - 13 kg) which increases or decreases the effort,

resulting in a progressive exercise programme.

03-010301: Saehan Cycle - Home Because the Saehan Cycle - Clinical

is a product aimed at professionals,

we offer various accessories to

increase its flexibility and user's safety. Be aware these do only apply to the clinical cycle, not to

the home cycle.

Reference	Description	Unit
03-010311	Clinical Cycle Accessory Long front leg	1 pc
03-010312	Clinical Cycle Accessory Extension arm (right = left)	1 pc
03-010313	Clinical Cycle Accessory Hand strap - right	1 pc
03-010314	Clinical Cycle Accessory Hand strap - left	1 pc
03-010315	Clinical Cycle Accessory Long (straight) leg brace set	1 pc
03-010316	Clinical Cycle Accessory Bent leg brace	1 pc

selected.

Reference	Description	Unit
03-010301	Saehan Cycle - Home	1 pc
03-010302	Saehan Cycle - Clinical	1 pc



Yesterday, you had the choice between a classic cane which is impractical because it is difficult to coordinate, does not help when it is in the air. No one likes having to use it because it highlights their difficulties... Or a rollator which requires two hands to use it, which stigmatizes as much as the cane, which is cumbersome and not easy to handle, let alone indoors (try walking through a door with a rollator!). Therefore, the walker is parked at the entrance of the house or apartment and leaves its user without help. Now there is the Wheeleo.

Reduced mobility presents 3 challenges:

- Keep the person moving;
- Reduce the risk of falling;
- Loss of autonomy.

Rolling Cane

Wheeleo, the rolling cane, is designed to offer a solution to these 3 challenges and offers additional benefits:

- Ease of handling, no special training required;
- Walking goes a lot faster than with a quadripod cane;
- Safer walking than with a quadripod cane, as there is constant ground contact;
- One-handed use while still maintaining constant ground contact, the other hand can be used to move an object, to open doors, ...
- Minimal space required and perfect for indoor use where a rollator is too big.
- Reinforced stability: the wheels offer a security system. A lateral push will result in a lateral movement, not a tilting/fall.



Wheeleo is a unique quadripod rolling cane used to improve mobility. It is designed with the aim of offering help without stigmatising the disability of its It will make walking trips smoother and safer which will allow the user to

The Wheeleo is used with one hand. Thanks to its 4 wheels, you do not have to lift it, which improves stability and reduces the risk of falling.


The Wheeleo is super easy to use. Ultra-intuitive, its user will always position it in the right place to guarantee its stability and balance.

An important feature of the Wheeleo is that it does not require brakes to prevent its user from falling. Thanks to its rolling design, the Wheeleo will always be placed in the right position for its user.

When losing balance, the user will automatically push the Wheeleo in the right position that is

comfortable for him, without lifting the Wheeleo. Adding brakes would be the same as replacing the Wheeleo with a regular cane because the brakes would block the wheels. As a consequence the Wheeleo would tilt and its user would possibly fall.

maintain or improve autonomy.





Adjustable

height




03-UPPER & LOWER BODY

MoVeS Oxycycle 1 - Active Magnetic MD (€









The OxyCycle 1 Pedal Exerciser is used for upper and lower body exercises. The adjustable resistance makes it a great product to progressively strengthen your arms and legs. The OxyCycle 1 Pedal Exerciser has a light-weight compact design, adjustable resistance and a multi-function LCD display readout of scan, time, distance, calories and RPM. A non-skid mat and interchangeable handgrip pedals are



The OxyCycle 1 comes with strap pedals and handgrip pedals.



Reference	Description	Unit
03-010101	MoVeS OxyCycle 1 Pedal Exerciser	1 pc








The OxyCycle 2 comes with

half-foot pedals and

handgrip pedals.

MoVeS Oxycycle 2 - Passive < €

For PASSIVE training only





For home use only

• 50x46x38 cm

- 8,5 kg
- 50 Hz

36

• 25 - 65 rpm



Comes with UK or EU plug!









MoVeS Oxycycle 3 - Active & Passive <€







For home use only

03-010103 MoVeS OxyCycle 3 Pedal Exerciser | EU-plug

03-010103 UK MoVeS OxyCycle 3 Pedal Exerciser | UK-plug



Reference Description



- 50x46x38 cm Comes with UK or EU plug!
- 8,5 kg
- 50 Hz
- 25 65 rpm



handgrip pedals are included. Adjustable RESISTANCE & ASSISTANCE



The OxyCycle 3 comes with full-foot pedals and handgrip pedals.

The MoVeS OxyCycle 3 Pedal Exerciser is a great tool for

personalised, low-stress workouts at home. The adjustable

speed level and resistance makes it a great product to

progressively strengthen your arms and legs. This pedal

exerciser provides low impact exercise to increase strength

and flexibility, is versatile for both upper and lower body.

It has a reliable 80 Watt motor to offer enough power to

maintain a constant rotational speed, whereas cheaper

models only carry a 30 Watt motor. A safety stop will

avoid injury in case of cramp or muscle spasm. The MoVeS OxyCycle 3 has a light-weight compact design, adjustable

speed & resistance, a 15-minute timer switch and a multi-

function LCD display readout of scan, time, count, total count and calories burned. A non-skid mat and interchangeable



For ACTIVE and PASSIVE training





MoVeS Oxycycle Pedals ∈ €

1 pc







Half-Foot Pedal

Reference	Description	Unit
03-020101	MoVeS Strap Pedals	pair
03-020102	MoVeS Half-Foot Pedals	pair
03-020103	MoVeS Full-Foot Pedals	pair
03-020201	MoVeS Handgrip Pedals	pair

03-UPPER & LOWER BODY

Curamotion C€

96 x 45 x 45 (cm) / 11,5 kg Wall mount included



EXER 2 146 x 45 x 54 (cm) / 13,5 kg Wall mount included




Improvement of muscle strength and joint range of motion. The CuraMotion Exer Series has an adjustable resistance and a multi-functional LCD-display readout scan which indicates rotation count, rotations per minute, time exercised and calories

All are easy adjustable with a quick release knob.





37 x 38 x 70 (cm) / 9 kg With foot stands







96 x 45 x 45 (cm) / 13 kg Wall mount included



MoVeS Pedal Exerciser MD (€

MoVeS Pedal Exercisers are ideal for restoring impaired muscle strength, improving blood circulation and give a better coordination to the legs or arms. It can be used on a table top as an arm exerciser or on the floor to give the legs a workout. A convenient resistance control knob allows a broad range of pedal resistance for a progressive exercise programme.

MoVeS Foldable Pedal Exerciser

- Foldable for compact shipping and storage, easy setup
- Pedal straps adjustable to two positions
- LCD Display (time, rotations per minute, total rotations and calories)
- In colour box





MoVeS Standard Pedal Exerciser | Premium Chrome

- Pre-assembled
- Pedal straps adjustable to two positions
- Quality chrome finish
- In colour box or on specific demand in bulk with separate colour box





MoVeS Standard Pedal Exerciser K/D

- Knock-down for compact shipping and storage (only 6 screws to assemble or disassemble)
- Pedal straps adjustable to two positions
- **Budget solution**



Reference	Description	Unit
03-010201	MoVeS Standard Pedal Exerciser Premium Chrome	1 pc
03-010202	MoVeS Standard Pedal Exerciser K/D White	1 pc
03-010203	MoVeS Foldable Pedal Exerciser + LCD Display	1 pc



MoVeS Wooden Incline Board MD < € The incline board increases extension in the lower and upper back while stretching the calves.



• to improve the general posture,

It can be used:

- to stretch the calf muscles of athletes, football players and sportspersons in general,
- to stretch the upper-leg muscles (Ischias),
- to increase extension in the lower and upper back.



MoVeS Wooden Incline Board - Fix



The MoVeS Wooden Incline Board - Fix is perfect for who's looking for an easy and budget-friendly Incline Board. It's made of polymerised wood for a strong feel and high user weight. The inclination (17,5°) is perfect for physiotherapy exercises.

150kg

Reference	Description	Unit
03-040001	MoVeS Wooden Incline Board - Ladder	1 pc
03-040003	MoVeS Wooden Incline Board - Fix	1 pc

Plastic Incline Board MD (€

Apart from stretching the calf and upper-leg muscles, the incline board is a fantastic device to improve general posture and extension in the lower back.

Made of a composite plastic for a maximum load of 250kg. Adjusts to 4 levels and folds flat (8cm) to store easily. Angles are 5°, 15°, 25° and 35°. The last slot on the board makes an excellent hand-hold for stretching. Non-slip top measures 35x35 cm.



Reference	Description	Unit
03-040102	Incline Board Standard	1 pc





Kinesiology Tape MD < €





MoVeS Kinesiology Tape is made from quality cotton fibres and is used for joint and muscle support. This skin-friendly Kinesiology Tape is hypoallergenic, latexfree, breathable and water resistant. The adhesive side of the Kinesiology Tape is made of acrylic which prevents sweat and moisture from penetrating the tape. This tape with a wave-like structure can be stretched up to 170% lengthways.

Available in 4 colours.





Reference	Description		Unit
04-060101	MoVeS Kinesiology Tape 5cm x 5m Tan	NEW	1 pc
04-060102	MoVeS Kinesiology Tape 5cm x 5m Black	NEW	1 pc
04-060103	MoVeS Kinesiology Tape 5cm x 5m Pink	NEW	1 pc
04-060104	MoVeS Kinesiology Tape 5cm x 5m Light blue	NEW	1 pc



NACI Massage Lotions

- Excellent Lubrication
- Economical usage
- Rinses off with water Skin-friendly
- Hypo allergen perfume
 - Prolonged hydration



Non greasy leaves skin smooth, soft and dry

www.nagi.com

NAQI° nv Stadsbeemd 1037, 3545 Halen, Belgium



This yoga set consists of:

- 1 pvc yoga mat (173 x 61 x 0,4 cm)
- 2 EVA foam yoga blocks (7,5 x 15 x 23 cm)
- 1 yoga strap with buckle (180 x 3,8 cm)

This set is packed in a full colour box. Items are also sold separately.



Mambo Max

Reference	Description	Unit
04-010201	Mambo Max Yoga Mat	1 pc
04-010202	Mambo Max Yoga Block	5-pack
04-010203	Mambo Max Yoga Strap	1 pc
04-010211	Mambo Max Yoga Set 1 Mat 2 Blocks 1 Strap	Set

Stepping Stones <€

The non-slip 16 cm Mambo Max Stepping Stone is an amazing tool to spice up your workout. They can be used by both children and adults for indoor and outdoor workouts to improve dynamic balance, body awareness and coordination. It can also serve as a foot massage tool while seated.



The Stepping Stones now come in bright and vivid colours: yellow, orange, red, purple, blue and

They are packed in cardboard boxes per 6 pieces and each piece has a different colour.



Reference	Description	Unit
04-030141	Mambo Max Stepping Stone	6-pack



™ ** Triggerpoint Rollers** C€

Massage rollers are used to put pressure on the muscles' triggerpoints and to loosen the tension in the muscles. Trigger points are often the cause of different pains.

Spiky Foam Roller



The Mambo Max Spiky Foam Roller is designed for deeper massage & acupressure therapy. The spikes make it possible to get deeper into the muscles to relieve tension and soreness and to improve mobility. It also has a bi-directional zone for comfort when rolling along the spine. The Spiky Foam Roller is mainly used for stretching, self-massage, body support and stability exercises.

Hollow Foam Roller

The Mambo Max Hollow Foam Roller can easily help people who are experiencing long term pains. Its three-dimensional pattern with different shapes replicates the massage from a therapist's hands. The larger surfaces allow a superficial massage while the smaller shapes are designed for a deeper massage.



Ergonomic Foam Roller



The Mambo Max Ergonomic Foam Roller is a curved hollow foam roller, which enables more exercises. Except rolling back and forth, it's possible to make sideways movements. More muscles will be addressed through which the body will be toned even more. The soft foam in combination with the nipples is perfect for trigger point massages. The central groove will reduce pressure on the spine.

3-in-1 Foam roller

The Mambo Max 3-in-1 foam roller combines 3 different foam rollers in one package. Thanks to its smooth surface, the outer part is perfect for a superficial massage. The second part has grooves which can be used for a deeper massage and to exercise higher pressure on the trigger points. Finally, with the inner massage stick it's possible to exercise high pressure on sore muscles to improve blood circulation and the muscles' recovery time. This last item is available as a separate product too: The Mambo Max Massage Stick.







02-HAND

Sitting Cushions <€



Mambo Max Sitting Cushions are air-filled cushions for active and dynamic sitting, recommended by therapists. These easy-to-carry air cushions activate intrinsic trunk muscles to support the spine and are suitable for everybody. The surface from these sitting cushions can easily be cleaned with soap and water.

Comfort Cushion 33





Packed in colour box with needle pump or in polybag only

Soft and comfortable surface


Needle pump

Standard Cushion 33

SensAir Cushion

Standard Cushion 37

Sitting Wedge





Faster Blaster pump



Ø 32cm

☐ 34cm Ф 7с 0 250kg



oth one senso side	

Reference	Description	Unit
04-020301	Mambo Max Comfort Cushion 33 Blue + Needle Pump	1 pc
04-020302	Mambo Max Comfort Cushion 33 Red + Needle Pump	1 pc
04-020303	Mambo Max Comfort Cushion 33 Black + Needle Pump	1 pc
04-021301	Mambo Max Comfort Cushion 33 Blue in Polybag	1 pc
04-021302	Mambo Max Comfort Cushion 33 Red in Polybag	1 pc
04-021303	Mambo Max Comfort Cushion 33 Black in Polybag	1 pc
04-020101	Mambo Max Standard Cushion 33 + Needle Pump	1 pc
04-020102	Mambo Max SensAir Cushion + Needle Pump	1 pc
04-020103	Mambo Max Standard Cushion 37 + Faster Blaster Pump	1 pc
04-020201	Mambo Max Sitting Wedge + Faster Blaster Pump	1 pc
04-020201	Mambo Max Sitting Wedge + Faster Blaster Pump	1 p



Pilates < €

Soft-Over-Ball 🍻



The Mambo Max Soft-Over-Ball is an inflatable ball used for resistance work in Pilates. It is a wonderful tool for adding focus and energy to your workouts. In addition, when slightly deflated,



it can be used as a soft prop under the head, neck, feet, back... Now available in 3 different colours and sizes, giving a range of 9 different Soft-Over-Balls.



Reference	Description	Unit
04-010105	Mambo Max Pilates Soft-Over-Ball - 17-19 cm Green	1 pc
04-010113	Mambo Max Pilates Soft-Over-Ball - 17-19 cm Blue	1 pc
04-010107	Mambo Max Pilates Soft-Over-Ball - 17-19 cm Black	1 pc
04-010115	Mambo Max Pilates Soft-Over-Ball - 21-23 cm Green	1 pc
04-010123	Mambo Max Pilates Soft-Over-Ball - 21-23 cm Blue	1 pc
04-010117	Mambo Max Pilates Soft-Over-Ball - 21-23 cm Black	1 pc
04-010125	Mambo Max Pilates Soft-Over-Ball - 25-27 cm Green	1 pc
04-010103	Mambo Max Pilates Soft-Over-Ball - 25-27 cm Blue	1 pc
04-010127	Mambo Max Pilates Soft-Over-Ball - 25-27 cm Black	1 pc





The Mambo Max Foam Roller has been designed for Pilates training and spine exercises (back aches) and is made of EVA. Ideal for balance and stabilisation exercises. Foam rollers are used for stretch and strength exercises, as well as for massage and relaxation.



Reference	Description	Unit
04-010102	Mambo Max Pilates Foam Roller - EVA 90cm	1 pc
04-010104	Mambo Max Pilates Foam Roller - EVA 45cm	1 pc



Foam covered fiberglass resistive exercise ring with foam padded grips. For exercises focusing on strength, balance, flexibility and posture. The Mambo Max Pilates Ring adds light resistance to pilates exercises and helps target specific muscle groups. Often used to help tone the arms and thighs, and increase the core strength challenge of an exercise.

Ø 38 cm

Reference	Description	Unit
04-010101	Mambo Max Pilates Ring	1 pc



03-UPPER & LOWER BODY

Massage Tools ← €

The Knobble® II massage tool is made of durable polymer with a band of non-latex. The Knobble® II gives your hand more power when you need it. Use either the tip or the base for applying pressure on yourself or others. The Knobble® II doesn't replace your healing hands, it just helps them do more and last longer.







The Jacknobber® II massage tool can be used in many ways to apply deep pressure to myofascial trigger points and areas of muscle tension. The smaller knobs provide more intense pressure while the larger knobs provide a more diffuse pressure. Use one knob at a time to press deeply into taut bands with embedded trigger points; use two knobs at a time on either side of the spine.

The Index Knobber® II massage tool is a simple device designed to let you apply deep pressure to muscles within easy reach of your hand. It is the ideal instrument for clinic or home use allowing for firm, sustained pressure on targeted muscle. It's unique design allows it to be used in several hand positions and by the right or left hands comfortably and interchangeably.







Reference

04-031101

04-031102



The Palmassager™ is one of the most efficient massage therapy aids available anywhere. Hold it firmly in the palm of the hand and use one, two or three of its round knobs to massage sore muscles. It is used for back, shoulder, arms and legs.

The Backnobber® II massage tool is allows the user to apply as much pain-releasing pressure as long as required to achieve remarkable





effects. A locking mechanism allows it to come apart for convenient travel and storage.

04-031103 Palmassager 04-031104 Backnobber II 1 pc 04-031105 Knobble II 1 pc

Description

Jacknobber II

Index Knobber II

Massage Ball, Roll & Ring <€

1 pc

The Mambo Max Massage Balls, Roll & Ring are used for relaxing muscles and reducing tension by relieving trigger points. The Mambo Max Massage Ball helps to strengthen hands and feet and to improve coordination and flexibility. It can be used for massage as well as to improve dexterity skills.



Reference	Description	Unit
04-030103	Mambo Max Massage Ball 6 cm	1 pc
04-030105	Mambo Max Massage Ball 7 cm	1 pc
04-030107	Mambo Max Massage Ball 8 cm	1 pc
04-030109	Mambo Max Massage Ball 9 cm	1 pc
04-030111	Mambo Max Massage Ball 10 cm	1 pc
04-030121	Mambo Max Massage Roll 16 cm	1 pc
04-030131	Mambo Max Massage Ring	1 pc








Gym Mat





The Mambo Max Gym Mats provide a soft touch and comfortable exercise environment. They are perfect for gym & fitness. Non-slip, well-cushioned surface with a unique closed cell structure. The mats are equipped with plastic rings 40cm apart so they can be suspended while not in use. Each mat comes with a handy storage/carry strap.





With rings

Reference	Description	Unit
05-030001	Mambo Max Gym Mat Wall Hanger for 10 mats	1 pc
05-030303	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) Blue Rings Strap	1 pc
05-030304	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) Black Rings Strap	1 pc
05-030305	Mambo Max Comfort Gym Mat 140 x 60 x 1 (cm) Magenta Rings Strap	1 pc
05-030101	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (cm) Red Rings Strap	1 pc
05-030102	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (cm) Green Rings Strap	1 pc
05-030103	Mambo Max Xtra Comfort Gym Mat 180 x 60 x 1,5 (cm) Blue Rings Strap	1 pc
05-030201	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Red Rings Strap	1 pc
05-030202	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Green Rings Strap	1 pc
05-030203	Mambo Max Xtra Comfort Gym Mat 180 x 100 x 1,5 (cm) Blue Rings Strap	1 pc

07-THERMAL

Aerobic Step < €



Mambo Max

Step aerobics is a form of aerobics that is practised worldwide by millions of people. The reason for its popularity is simple, it offers a great cardio workout without the need for extensive space or equipment.

While the aim of aerobics is to increase the body's consumption of oxygen, step aerobics builds on that goal by offering a more intense workout. The Mambo Max Aerobic Step provides this extra boost by incorporating the action of stepping onto the platform to intensify the workout. It is low impact,

Needle Pump

Ball Stacker

with velcro

Ball Stabiliser

The Ball Stabiliser facilitates balance. Its notches enable combined AB


Gym Ball / MoVeS Band exercises.

closure

and doesn't stress the joints as much as running or

jogging.





Reference	Description	Unit
05-070101	Mambo Max Aerobic Step	1 pc



Electric Compressor Pump





Reference	Description	Unit
05-020101	Mambo Max Ball Stabiliser	1 pc
05-020102	Mambo Max Ball Stacker	3 pcs
05-020201	Mambo Max Faster Blaster Pump 2x120cc	1 pc
05-020202	Mambo Max Hand Pump 2x2000cc	1 pc
05-020203	Mambo Max Needle Pump 100cc	1 pc
05-020204	Mambo Max Foot Pump	1 pc
05-020205	Electric Compressor Pump	1 pc





Fabric Elastic Band < €





The Mambo Max Fabric Elastic Band consists of 10 consecutive loops. This great toning and work-out tool is suitable for early training up to professional athletes. Three different resistances are available, while the 10 consecutive loops allow you to spice things up. Physical training can be intensified by using a lower amount of consecutive loops resulting in a higher resistance.



The Mambo Max Fabric Elastic Band is durable and extremely versatile. It can be used for stretching, toning or even as an alternative for dumbbells as a weight training accessory. The progressive resistance offers you a broad range of exercises.



COLOUR	RESISTANCE	REFERENCE	UNIT
Light	7 kg	05-120107	1 pc
Medium	10 kg	05-120110	1 pc
Heavy	15 kg	05-120115	1 pc





Half sphere, inflatable exercise dome. This versatile product enables the user to perform exercises on both sides of the product, allowing it to be used for a full body work-out as well as for balance exercises. Basically it's Mambo Max' economical answer to the BOSU® system, but with extra tubing and handles!



\$ 21-25cm 175kg

Ø 59cm







02-HAND

03-UPPER & LOWER BODY

Balance Board ← €



The Mambo Max Balance Board is an effective tool for balance training and fitness exercises. Stand on the rounded platform and start balancing. It is also a great tool to strengthen your ankles and thigh muscles. For a more severe balance training, add the included screw-on top to your Mambo Max Balance Board.

Ø 42cm \$ 5-8cm 100kg

g	

Reference	Description	Unit
05-040201	Mambo Max Balance Board	1 pc

Wooden Balance Boards <€

The Mambo Max Wooden Balance Board is your perfect partner if you want to improve your balance on an ecologic and high-end product. The Wooden Balance board is available in 3 versions: Uni-Directional, Standard & Sport. Their different form and balance point make it possible to have 3 products with a very different character and level of difficulty.

Children and elderly could experience a balance board with pivot point as too challenging. The Uni-directional Balance Board offers more confidence. This enables the user to improve body stability without the fear of falling.

A Uni-directional Balance Board offers a more predictable training as it can only pivot in one direction. When turning the Balance board upside-down, the arches can be used as push-up grips. Once the user is ready for the next level of training, the Balance Board Rehab or Sport are ideal next steps in the process.

The Balance Board Standard is lower than the Balance Board Sport. In addition with its smoother pivot point, it reacts completely different than the Sport. It's a great middle step between the Balance Board Uni-Directional and Balance Board Sport.

The most challenging Balance Board is our Sport version, which, as the name specifies, makes it a real sport. Finding your balance and avoiding the edges touching the ground is hard and requires you to use all your muscles in your body. It's a great and fun way to exercise a fullbody workout.

Reference	Description	Unit
05-040001	Mambo Max Wooden Balance Board - Standard	1 pc
05-040002	Mambo Max Wooden Balance Board - Sport	1 pc
05-040003	Mambo Max Wooden Balance Board - Uni-directional	1 pc







Standard: 40cm x 40cm / \$6cm Max inclination: 12,5°







Balance Board Freeman <€

The Mambo Max Balance Board Freeman is an effective tool for balance training. Stand on the rounded platform and start balancing. This is a great way to strengthen your ankles and thigh muscles. The handgrips on the bottom are used for upper body balancing and workouts.







	Reference	Description	Unit
05-040202 Mambo Max Balance Board Freeman 1 pc	05-040202	Mambo Max Balance Board Freeman	1 pc





ANTI-BURST

350kg Yes, 350 kg!!!

The Mambo Max AB (Anti-Burst) Gym Ball increases flexibility, improves strength, fitness and posture. Exercising with this ball also enhances your balance and coordination. All AB Gym Balls are Anti-Burst: if the ball gets punctured by a sharp object the ball won't burst, thus reducing

This extremely durable Mambo Max AB Gym Ball comes in a Mambo Max cardboard retail pack, including AB Gym Ball, Faster Blaster Pump, Exercise Chart, Plug and Plug Remover.



The Mambo Max AB Gym Ball is also available in the value pack: a polybag containing only the ball and plug.



Ecofriendly pvc. Heavy metal content in compliance with EN 71 Standards. Phthalates content meets EU-directive 2005/84/EC.

Ball Size	Colour	Pody Hoight	Mambo Max AB Gym Ball	Mambo Max AB Gym Ball	Unit
Dali Size	Coloui	Body Height	Box + Faster Blaster Pump	Polybag	
45cm	Yellow	< 125 cm	05-010102	n/a	1 pc
55cm	Red	125-150 cm	05-010103	05-010203	1 pc
65cm	Green	150-175 cm	05-010104	05-010204	1 pc
75cm	Blue	175-185 cm	05-010105	05-010205	1 pc
85cm	Black	185-200 cm	05-010106	n/a	1 pc
95cm	Silver	> 200 cm	05-010107	n/a	1 pc



AB Peanut Ball ←



The Mambo Max AB (Anti-Burst) Peanut ball helps users who struggle with the instability of a regular gym ball. The peanut shape offers extra sideway stability, which gives more confidence to start working out with Gym Balls.



The Mambo Max AB Peanut Ball offers a more predictable therapy or fitness training as it can only roll in one direction. This certainly is a huge benefit for children, elderly or people who suffered an injury. In addition, Peanut Balls are very popular amongst pregnant women and are very helpful devices during la-

Similar to our regular Mambo Max AB Gym Ball, the Peanut Ball is

Anti-Burst: if the ball gets punctured by a sharp object the ball won't burst, thus reducing risk of

The Mambo Max AB Peanut Ball comes in a Mambo Max cardboard retail pack, including Faster Blaster Pump and Plug.

Peanut Ball Size	Colour	Mambo Max AB Peanut Ball + Faster Blaster Pump	Unit
50 x 100cm	Red	05-011103	1 pc

Body Twister ←



The Mambo Max Body Twister is your ideal partner to slim your thighs, waist and hips. While twisting your body, the disc below your feet turns and glutes, abdominals and legs are being trained. Continue to twist your body from side to side to tone all muscles in your body. The Mambo Max Body Twister is also an effective product to enhance balance and coordination.

Reference	Description	Unit
05-100101	Mambo Max Body Twister	1 pc





Ekilibre <€



The Mambo Max Balance Plate "Ekilibre" is designed for balance training. Combining balance pad with balance board, it offers the best of both worlds. It improves your sense of equilibrium, coordination, reaction and motor skills. It can be used for lying, kneeling, sitting and stepping movement and it can be incorporated with other exercise equipment.

Beginner	
Commercial	☐ 60cm x 39cm
M. 6.31 1	\$ 9cm
	100kg
	■ C 100kg



Intermediate	Ø 39cm
ALL LAND AND DESCRIPTION OF THE PARTY OF THE	
	1

Reference	Description	Unit
05-040501	Mambo Max Ekilibre Beginner *	1 pc
05-040502	Mambo Max Ekilibre Intermediate *	1 pc

*can not be sold in Canada





With the Mambo Max Wooden Work-Out Station an even higher amount of exercises is suddenly available. The cleverly spread rings allow you to attach a resistive exercise band or tubing for various vertical and diagonal exercises.

The Work-Out Station is designed to fit an Aerobic Step or a set of Oval Balance Pads. This mix of various products results in a complete and challenging workout. Combining balance with resistive exercise triggers many different muscles in your body. The Mambo Max Wooden Work-Out Station truly offers a complete and three-dimensional workout.

Outer Dimension: 103 x 40 (cm) Inner Dimension: 96 x 33 (cm)



Reference	Description	Unit
05-110001	Mambo Max Wooden Work-Out Station + Tube with Handles + Carabiner	Set
05-110099	Mambo Max Pair of Tubes with Handles + Carabiner	Pair

Balance Pad Pro (€

Mambe and my body body

This balance pad is a great training tool to improve your balance, stability, flexibility, strength and motor skills. The Balance Pad Pro creates more instability than our regular Balance Pad, so you constantly need to correct body posture, making the excercises just that bit tougher. This new balance pad is made out of PU-material instead of EVA-foam, which makes it more durable and less prone to small damages.



Dimensions: 40 x 50 x 6cm

Reference	Description		Unit
05-040304	Mambo Max Balance Pad Pro Rectangular Blue	NEW	1 pc



Balance Pad ←

This EVA-foam balance pad provides the comfortable exercise environment needed to perform reinforcing balance exercises. Its soft and thus unstable surface challenges the user to work the foot-, ankle-, and leg muscles.







Available in three versions:


- 1 rectangular right angle shaped pad, 50 x 40 x 6 (cm), blue
- 1 rectangular chamfered shaped pad, 50 x 40 x 6 (cm), anthracite
- 2 oval shaped pads, 36 x 22 x 6 (cm) sold per pair



Description	Unit
Mambo Max Balance Pad Rectangular Blue	1 pc
Mambo Max Balance Pad Rectangular Anthracite	1 pc
Mambo Max Balance Pad Oval Blue	pair
	Mambo Max Balance Pad Rectangular Blue Mambo Max Balance Pad Rectangular Anthracite

Ob ex

Balance Trainer ∈ €



The Mambo Max Balance Trainer is designed for balance training. Use it to improve your sense of equilibrium, coordination, reaction and motor skills. It is also ideal for muscle exercise. The Mambo Max Balance Trainer can be adjusted by can be used as inflating and deflating the air through the valve. It can be used for Sitting Cushion lying, kneeling, sitting and stepping movement and it can be incorporated with other exercise equipment to spice up your rehabilitation.



Reference	Description	Unit
05-040102	Mambo Max Balance Trainer 45	1 pc
05-040103	Mambo Max Balance Trainer 60	1 pc



FAIRS



Therapie Hamburg Hamburg - Germany 22 - 23 October 2021

Therapie Hamburg is the platform in northern Germany for physiotherapists, occupational therapists and sports therapists

as well as speech therapists, masseurs and lifeguards. Innovations, products and industry trends are shown here.



MEDICA Düsseldorf - Germany

MEDICA has been the largest medical trade fair in the world and the event on the agenda of medical experts for over 40 years. In 2019, the fair welcomed 5,100 exhibitors, nearly

15-18 November 2021

130,000 visitors from 144 countries and had a satisfaction rating for closed deals of 96%.



ARAB HEALTH Dubai - UAE

24 - 27 January 2022

Arab Health is the largest gathering of healthcare and trade professionals in the MENA region. With over 55,000 local and international visitors, 4,250 exhibitors and participants from over 150 countries, Arab Health is truly where the healthcare world comes to do business.



7 - 10 April 2022 **FIBO Global Fitness**

is an expo for sports,

nutrition, and healthy lifestyle. It brings global top decision makers of the industry together looking for the latest innovations and trends for their gyms, wellness facilities, doctors and physiotherapy practises, etc...

Swim Noodle <€



Noodles can be used as resistance equipment for muscle strengthening or as floating device. They are made of ultra-flexible foam that can bend into any shape. Ideal for aquagym, for fitness exercises and a perfect way to strengthen muscles in the water, especially the abdominals. Also great for relaxation and fun in the pool.

Ø 7cm ↔ 150cm

	Unit	Reference
05-050101 Mambo Max Swim Noodle	1 pc	05-050101







Coated Dumbbell < €

The Mambo Max Dumbbells are cast iron dumbbells, coated with vinyl or neoprene. The vinyl coating is easy washable or can be treated with a water based disinfectant.

The neoprene coating has an added resin to give them a scrubby feel and superior grip.

Available in 7 weights with matching colour code. The Dumbbell Rack holds 6 Dumbbells (3 pairs).







Description	Weight	Vinyl	Neoprene	Unit
Mambo Max Dumbbell	0,5kg	06-010102	06-010202	pair
Mambo Max Dumbbell	1,0kg	06-010103	06-010203	pair
Mambo Max Dumbbell	1,5kg	06-010104	06-010204	pair
Mambo Max Dumbbell	2,0kg	06-010105	06-010205	pair
Mambo Max Dumbbell	3,0kg	06-010107	06-010207	pair
Mambo Max Dumbbell	4,0kg	06-010108	06-010208	pair
Mambo Max Dumbbell	5,0kg	06-010109	06-010209	pair
Mambo Max Dumbbell rack (empty)	n/a	06-00	00101	1 pc





Mambo Max



SoftMeds <€

Mambo Max SoftMeds are colour coded hand weighted balls that are small and flexibleenough to grasp with one hand, yet big enough to hold with both hands. The soft, pliable vinyl cover makes it easy to hold and is a safe and effective way to add weight training to your workout.

Mambo Max SoftMeds are fully inflated and ready to use, straight out-of-the box. When needed, air pressure can be adjusted with a needle pump.






Ø 12cm

Reference	Weight	Description	Unit
06-040102	0,5kg	Mambo Max SoftMeds Tan	1 pc
06-040103	1,0kg	Mambo Max SoftMeds Yellow	1 pc
06-040104	1,5kg	Mambo Max SoftMeds Red	1 pc
06-040105	2,0kg	Mambo Max SoftMeds Green	1 pc





Wrist & Ankle Weight <€



Reference	Description	Weight	Unit
06-020102	Mambo Max Wrist & Ankle Weights	0,5kg	pair
06-020103	Mambo Max Wrist & Ankle Weights	1kg	pair
06-020104	Mambo Max Wrist & Ankle Weights	1,5kg	pair
06-020105	Mambo Max Wrist & Ankle Weights	2kg	pair
06-020107	Mambo Max Wrist & Ankle Weights	3kg	pair
06-020108	Mambo Max Wrist & Ankle Weights	4kg	pair
06-020109	Mambo Max Wrist & Ankle Weights	5kg	pair
	-		



To quickly build stamina and muscle, simply add Mambo Max Wrist and Ankle Weights to your workout. They are filled with iron sand and will effectively shape and tone your body and can be worn anywhere. Wear them around the office, out shopping or at the gym. Eventually, your body will adjust to running, walking and living with this extra weight.



Medicine Ball < €



A medicine ball (also known as an exercise ball, a med ball, or a fitness ball) is often used for strength training as it plays an important role in the field of sports medicine. Medicine ball training is one of the oldest forms of strength and conditioning training – the first reference to wrestlers training with sand filled bladders appears in Persia nearly 3000 years ago. In ancient Greece the physician Hippocrates had his patients throw them back and forth for injury prevention.

> 28cm x 57cm x 183cm $(I \times w \times h)$

Reference	Description	Ø	Unit
06-030103	Mambo Max Medicine Ball - 1kg	20cm	1 pc
06-030104	Mambo Max Medicine Ball - 1,5kg	20cm	1 pc
06-030105	Mambo Max Medicine Ball - 2kg	20cm	1 pc
06-030107	Mambo Max Medicine Ball - 3kg	22cm	1 pc
06-030108	Mambo Max Medicine Ball - 4kg	22cm	1 pc
06-030109	Mambo Max Medicine Ball - 5kg	27cm	1 pc
06-000301	Mambo Max Medicine Ball rack (empty)		1 pc

Mambo Max

Thumb Lock Wrist Weight <€





The Mambo Max Thumb Lock Wrist Weight enables you to fix the weight to your thumb, ensuring a fixed position during your entire workout. Being it boxing, running or Zumba, the weights will remain in place.

Reference	Description	Weight	Unit
06-020202	Mambo Max Thumb Lock Wrist Weights	0,5kg	pair


Kettlebell < €

The Mambo Max Kettlebells are an important part of any fitness or CrossFit training. Made from high quality iron surrounded with a colourful vinyl coating, they will not damage any surfaces

> when being stored. Starting from 2 kg up to 24 kg, we also offer you a complete set, including a premium quality Kettlebell Rack, which can hold up to 150kg per shelf (total of 300kg).

> > Available in: 2 - 4 - 6 - 8 - 10 - 12 - 16 - 20 - 24kg

120cm x 50cm x 70cm (l x w x h)







02-HAND

the digital thermometer.

Section 2 Para Pro - Paraffin Bath C€

The Para Pro is a professional paraffin wax heating bath which is made to be used on 220/230 Volt power sockets and which can be used 24/7. In combination with our premium paraffin pearls it only needs between 2 and 3 hours to completely melt 2,5kg of paraffin. The temperature is adjustable from 40-60°C and is easy to control due to





- Green fully refined paraffin wax
- Melting point: 46°C
- Tea Tree scented




When ordering the Para Pro, you receive 2,5kg of white paraffin pearls or 6 blocks (2,72kg) of Tea Tree Green Paraffin free of charge.

Para Pro Specifications:

- Outer dimensions: 39cm x 25cm x 18cm (I x w x h)
- Inner dimensions: 32,5cm x 18cm x 13cm (I x w x h)
- Weight: 3,1kg
- 220 V: No transformer needed, which will lengthen the durability.
- Low power consumption once the paraffin has melted
- 24-hours a day fully operational
- Fast melting between 2 3 hours!
- Adjustable temperature (40 60°C)
- Digital thermometer

Paraffin Pearls:

- White fully refined paraffin wax
- Congealing point: 53°C
- Melting point: +/- 54°C
- Meets the purity requirements of BGA1, Ph.EUR.2 and FDA3
- 1. Bundesgesundheitsamt recommendation No. XXV B II;
- 2. European Pharmacopoeia, current edition;



Reference	Description	Unit
07-100101	MoVeS Para Pro Bath + 2,5kg Paraffin Pearls - White EU-plug	1 pc
07-100101 UK	MoVeS Para Pro Bath + 2,5kg Paraffin Pearls - White UK-plug	1 pc
07-100206	MoVeS Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg) - Tea Tree Green EU-plug	1 pc
07-100206 UK	MoVeS Para Pro Bath + Set of 6 Paraffin Blocks 46°C (2,72 kg) - Tea Tree Green UK-plug	1 pc
07-110102	MoVeS Paraffin Pearls - 2,5kg - White	1 pc
07-110125	MoVeS Paraffin Pearls - 25kg - White	1 pc
07-110206	MoVeS Paraffin Blocks 46° Set of 6 x 453g (2,72 kg) - Tea Tree Green	1 pc
07-109801	MoVeS Insulating Gloves for Para Pro Bath Blue	1 pair
07-109802	MoVeS Disposable Gloves for Para Pro Bath	10 pairs

Comes with UK or EU plug!



PearlPack MD (€



PearlPacks are reusable Hot/Cold Packs filled with small pearls. Hot and Cold Therapy can be applied to relieve pain and inflammations associated with muscle strain and sprain. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The small pearls absorb heat or cold, which they can retain for 20 minutes. During these 20 minutes, these packs will relieve a part of the pain.

Not only do they look good, the products are very convenient in use. PearlPacks won't feel wet nor drip after freezing or heating them, avoiding a mess. All materials used are non-toxic and odourless.

Besides relieving pain from an injury or muscle strain, hot/cold packs are often used to reduce swelling after an operation. PearlPacks truly are multifunctional products.

PearlPack Knee/Elbow



This offers a solution for Hot/Cold therapy on the knee or elbow joints. The large opening ensures that your knee or elbow remains free and can continue to move freely, while the Pearl technology sends direct relief to the pain. The flexible velcro strap keeps the pack firmly in place.





PearlPack Eye Mask



This mask adapts to the patient's face and has a fabric strap. Hot treatment can be applied when eyes feel dry or if there's an infection of the eyelid. Cold treatment will help when suffering from headaches, swollen eyes or sinusitis.





PearlPack Back Wrap



This large pack divided in 3 segments is equipped with velcro strap. It fits perfectly around the patient's waist for an optimal pain relief on the back. The pearls can be pushed from one segment to the other if one area of the back needs less or more thermal treatment.



07-040201

PearlPack Sports Pack + strap



The Pearl Pack Sport pack is updated and now comes with a velcro strap.





Reference	Description		Dimension	Unit
07-040102	MoVeS PearlPack - Sport Pack + strap	UPDATED	19cm x 12,5cm	1 pc
07-040201	MoVeS PearlPack - Back Wrap		43cm x 17,5cm	1 pc
07-041101	MoVeS PearlPack - Eye Mask		22,5cm x 7,5cm	1 pc
07-040301	MoVes PearlPack - Knee/Elbow		37cm x 25,6cm	1 pc

04-WELLNESS

60







MoVeS Roll-On packs combine hot/cold treatment with practicality. The tubular shape and flexible material allow an easy roll-on application, even after being heated or frozen. The various sizes will fit perfectly over finger, wrist, ankle or knee. This design ensures that the whole area will benefit from the hot or cold treatment.

From the same factory as the original Torex® products. Reusable, non-toxic and latex free. Available in 5 sizes.





					Access to the second se	
	Reference	Description	Inside Ø	Circumference	Uses:	Unit
	07-050101	MoVeS Roll-On Pack Finger	2,5cm	One size fits all	Finger	1 pc
	07-050102	MoVeS Roll-On Pack Small	10,0cm	Fits up to 25cm	Hand, Wrist, Arm, Elbow	1 pc
	07-050103	MoVeS Roll-On Pack Medium	12,5cm	Fits 25cm - 38cm	Hand, Wrist, Arm, Elbow, Ankle	1 pc
	07-050104	MoVeS Roll-On Pack Large	15,0cm	Fits 38cm - 53cm	Elbow, Ankle, Foot, Calf, Knee	1 pc
-	07-050105	MoVeS Roll-On Pack Extra Large	17,5cm	Fits 53cm - 58cm	Knee, Thigh	1 pc



Torex® Pro Ankle Boot A high performance ankle hot/cold and compression brace designed by PRO Athletic Trainers.



Reference Description 07-030201 Torex Hot/Cold Pro Ankle Boot 1 pc



Hot/Cold Pack MD (€

The MoVeS Hot/Cold Packs are usable for hot as well as cold therapy. This dual applicability makes them a strong alternative to ordinary heating pillows and ice packs. The versatile, reusable packs deliver hot and cold therapy effectively and economically. The use of the MoVeS Hot/Cold Pack is unlimited and its material is non-toxic. Even below -20°C it remains soft. The MoVeS Hot/Cold Pack can be bent in any angle, allowing maximum coverage of joint and muscle. All Hot/Cold Packs can be heated in the microwave and can be frozen in the freezer. Temperature holds for approximately 30 minutes.

MoVeS Hot/Cold Pack | Classic - Gel

The classic (and economical) transparent Hot/Cold Pack. The special gel is non-toxic and biodegradable.

Sold per 10 pieces

VIoVeS

OT / COI

PACK

or home use only

MoVeS Hot/Cold Pack | Standard - Gel

Our standard Hot/Cold Pack is a nylon covered gelpack. The special gel is non-toxic and biodegradable.

MoVeS Hot/Cold Pack | Soft Touch - Clay

The Soft Touch is improved with a more user friendly surface. The frost-free surface allows usage in direct contact with the skin. Filled with clay, a 100% natural product.

Sold per 10 pieces





Protection sleeve

100% polyester Protection Sleeve for hot/cold packs. Machine washable. Available in 3 sizes.



07-010206: XXL Hot/Cold Pack

This is ideal for pain relief on shoulders and back and has a nylon cover.



07-010311: Tri-Sectional **Hot/Cold Pack**



Dimension		N	NoVeS Hot/Cold Pag	cks	Protection Sleeve	Lloit
Diffi	lension	Classic (gel)	Standard (gel)	Soft Touch (clay)	100% polyester	Unit
15cm x 25cm	Small	07-010102	07-010202	07-010302	07-020102	10 pcs
20cm x 30cm	Medium	07-010103	07-010203	07-010303	07-020103	10 pcs
25cm x 35cm	Large	07-010104	07-010204	07-010304	07-020104	10 pcs
33cm x 47cm	XXL	n/a	07-010206	n/a	n/a	1 pc
40cm x 20cm	Tri-Sectional	n/a	07-010211	07-010311	n/a	10 pcs

04-WELLNESS



Infrared Thermometer MD < €

Ear & Forehead Thermometer

The infrared thermometer is used for measuring body temperature by two-in-one functions: Ear and Forehead temperature. It is suitable for medical unit and home use.

- Fast and easy temperature measurement by advanced infrared technology
- Range between 34.0°C~43.0°C Accuracy; <35.0°C and >42.0°C = ±0.3 °C 35.0°C~42.0°C = ±0.2°C
- Tri-colour backlit;
 - <37.4°C = Green
 - 37.4°C~38.0°C = Yellow
 - >38.0°C = Red
- Memory mode
- Suitable for both infants and adults
- Easy operating and cleaning
- Automatic power off after 60 seconds
- 2 AAA alkaline batteries









Shotgun Thermometer

Measures your body temperature without direct contact with the skin. Can measure between 5-15 cm from the skin.

- Detailed 34 memories with time/date
- Response time: 1 second
- Effective distance: 5-15cm
- User guidance by infrared beam
- Fever alarm







Reference	Description	Unit
08-070301	Infrared Thermometer Ear & Forehead auto-sensing 2 x AAA	1 pc
08-070303	Infrared Thermometer Shotgun Forehead auto-sensing 2 x AAA	1 pc



Blood Pressure Monitor MD (€

The blood pressure monitor is used for measuring blood pressure and pulse rate. It is suitable for people older than 12 years.



08-070501 Blood Pressure Monitor + Cuff

1 pc

- Large LCD display: 72mm x 43,5mm
- 2 users, 99 memories per user
- Irregular heart beating sign
- Indicating blood pressure level in 4 colours
- Average of last 3 records
- Low power detection
- Over pressure protection
- Automatic switchoff



Į.	Available early 20



Reference Description

Fingertip Pulse Oximeter MD (€

A Pulse Oximeter is a very important and common device to check patient blood-oxygen saturation (SpO₂) level and pulse rate. As a small, compact, simple, reliable and durable physiological monitoring device, Fingertip Pulse Oximeter greatly enhances patient care.



08-070201

Junior Oximeter

Rechargeable oximeter with OLED-screen, specifically designed for kids from 3-12 years old.

- Rechargeable via USB-cable
- Dual colour OLED displays SpO₂, PR, pulse bar and
- 4-direction & 6-mode display for convenient readings
- Set-up alarm range of SpO₂ and pulse rate
- Menu-function setting (beep sounds, etc.)
- Lithium battery; automatically power off
- Gravity function, automatic rotating reading option



Features:

- Dual colour OLED displays SpO₂, PR, pulse bar and
- SpO₃; measurement range between 70%~99% Accuracy; ±2% on the stage of 80%~99%
- PR; measurement range between 30BPM~240BPM Accuracy; ±1BPM or ±1% (the larger one)
- Low power indicator
- 6 display modes with adjustable brightness
- Automatically power off
- AAA alkaline batteries



08-070202

Reference	Description	Unit
08-070201	Fingertip Pulse Oximeter Adult OLED 2 x AAA	1 pc
08-070202	Fingertip Pulse Oximeter Junior 3-12y old OLED Rechargeable	1 pc

62

\sim Hydraulic 3-piece hand evaluation kit $\subset \in$





Portable carrying case contains all instruments necessary to evaluate the strength and range-of-motion of the hand

Set includes:

- hydraulic hand dynamometer
- hydraulic pinch gauge
- stainless steel finger goniometer 15cm 180° per 5°



Reference	Description	Unit
08-010113	Saehan Hydraulic Hand Evaluation Kit 3 pieces	Set

Hydraulic 7-piece hand evaluation kit $\in \mathbb{C}$





Portable carrying case contains all instruments necessary to evaluate the strength and range-of-motion of the hand

Set includes:

- hydraulic hand dynamometer
- hydraulic pinch gauge
- stainless steel finger goniometer 15cm 180° per 5°
- 2-point discriminator with 3rd point
- Wartenberg pinwheel
- finger gauge
- functional finger motion gauge



08-010117 Saehan Hydraulic Hand Evaluation Kit 7 nieces Set	Reference	Description	Unit
00-010117 Sachan Hydraune Hand Evaluation Rit 7 pieces Set	08-010117	Saehan Hydraulic Hand Evaluation Kit 7 pieces	Set



08-010101 Saehan Hydraulic Hand Dynamometer

The Hydraulic Hand Dynamometer gives accurate grip strength readings without the patient being able to "feel" the handle move. The adjustable handle has 5 positions to accommodate to any hand size. The indicator remains at the patient's maximum reading until it is reset. Has a scale in kilograms. Comes in a hard carrying case.

1 pc

Range: 0 - 90kg (200lb)

Quality

Quality Hard

Software



§SAEHAN`

Digital hand dynamometer with	
Data Transfer Software- DHD-1	$\boldsymbol{\zeta}$



calibration system and up to 40 test data can be stored in memory. It is powered by a rechargeable Li-Ion polymer battery and data will not be lost when changing the battery. The Data Transfer

transfer the test results into .txt-files, which can be easily accessed

using Notepad or MS-Word for editing and printing.

This Digital Hand Dynamometer hard carrying case includes a Digital Hand Dynamometer, the DTS-software and a USB cable.

Range: 0 - 90kg (200lb)

Reference	Description	Unit
08-010201	Saehan Digital Hand Dynamometer - DHD-1	1 pc

Digital hand dynamometer with G-STAR™ Software - DHD-3

DHD-3 Digital Hand Dynamometer (SH1003) prosame features as the DHD-1 when used a stand-alone device, but adds automated Rapid Exchange Grip Testing. Using Auto Mode, up to 10 strength tests for each hand as well as the time allowed to switch hands can be configured to automate the testing procedure.



The Grip Strength Testing and Research (G-STAR) is an easy to use software tool which automates grip strength testing and endurance testing in both clinical and research environments. All patient information and their test results are stored in an Access database.

Range: 0 - 90kg (200lb) / Minimum requirements: Windows XP + USB port

Reference	Description	Unit
08-010202	Saehan Digital Hand Dynamometer - DHD-3 (incl. G-STAR software)	1 pc



SAEHAN

SAEHAN

Smedley Hand Dynamometer < €





The Smedley spring-type Dynamometer has been successfully used for many years to measure grip strength. The ergonomic handle is adjustable to

Unit is calibrated in kilograms and pounds. Stored in a soft carrying case.



Range: 0 - 100kg (220lb)

SAEHAN
e 11 e

Reference	Description	Unit
08-010401	Saehan Smedley Hand Dynamometer - Ergonomic Handle	1 pc

Push-Pull Dynamometer MD (€



A simple, easy-to-use, ergonomically designed instrument that objectively measures push, pull and lift forces for manual muscle testing, functional capacity evaluation and job task evaluation at a remarkably affordable price. Because the instrument is lightweight, small and portable, you can perform precise, objective evaluations in your practice or at the patient's location. The easy-toattach dual grip handle can be used when measuring larger forces.

Lift platform and chain to perform lifting evaluation not included.



Reference	Description	Unit
08-010501	Baseline Hydraulic Push-Pull Dynamometer	1 pc
08-010502	Dual Grip for Push-Pull Dynamometer	1 pc





Body Fat Scale MD (€

Baseline body fat scale is made with tempered glass. This scale measures weight, body fat and water % and memorises data up to 8 persons. Measurements range from 4% - 45% body fat; 37.8% - 66% body water and 330lb/150kg weight capacity. The required 1 x 3V battery is included.


Resolution:

- 0.1% bodyfat

- 0.1% water

08-070101

Dimensions: 35cm x 35cm x 6cm Product weight: 2kg





Pinch Gauge <€

08-020102: Mechanical Pinch Gauge

Used for hand-strength testing by measuring tip, key and palmar pinch force. Made of durable painted stainless steel. Mechanical properties will not degrade over time. It can serve as a diagnostic tool or for monitoring the progress of rehabilitation.



Range: 0 - 20kg (45lb)

Reference	Description	Unit
08-020101	Saehan Hydraulic Pinch Gauge	1 pc
08-020102	Saehan Mechanical Pinch Gauge	1 pc

08-020101: Hydraulic Pinch Gauge

Uses the improved hydraulic system of the hand dynamometer to assure convenience, product reliability, measurement accuracy and repeatability. This results in a more accurate test result for all pinch tests.

\$SAEHAN



Squeeze Dynamometer < €





The Squeeze (bulb) Dynamometer is the inexpensive device to reliably measure hand and finger strength. The dynamometer has an indicator dial that remains at the maximum reading until reset. Unit is calibrated in pounds and kilograms.

Range: 0 - 70kg (150lb)

Reference	Description	Unit
08-010301	Saehan Squeeze Dynamometer 70kg	1 pc

Back-Leg-Chest Dynamometer < €

Measure strength of back, leg and chest. Base provides secure footing. Chain length is adjustable to accommodate for height differences or to vary the point of force application. Shows pounds and kilograms. Pointer remains at maximum until reset. Comes fully assembled with 1,5m chain. Range: 0 - 300kg (660lb)

Reference	Description	Unit
08-010601	Saehan Back-Leg-Chest Dynamometer	1 pc



SAEHAN

03-UPPER & LOWER BODY


Fat Caliper < €



Professional light-weight aluminium caliper. Measurement arms with floating tips to assure accurate and reproducible results. 65mm scale is large and easy to read. Corrective scale for sex and age included in multilingual manual.

- Professional, very accurate unit with spring loaded levers and low-pressure bearings for constant pressure.
- Scale permits reading up to 65mm
- Aluminium casted to insure years of reliable use.
- Includes carrying case
- CE marked and FDA approved

SAEHAN





The FatTrack® FITNESS TK Body Transformation Toolkit enables you to measure (and track) Body Fat %, Lean Body Mass, your workouts, macronutrient information and more by yourself - Easily, accurately and in the privacy of your own home!

Each FatTrack FITNESS TK Includes:

- Updated FatTrack® Fitness Digital Body Fat Caliper
- User Manual with instructions
- **Body Tracker Software**
- MyoTape Body Tape Measure
- FREE copy of "Improve Your Measuring Skill" with practical pointers for getting the best measurement readings and valuable information on ideal weight and muscle mass.

08-040102

Reference	Description	Unit
08-040101	Saehan Skinfold Fat Caliper Aluminium - Green	1 pc
08-040104	Baseline Skinfold Fat Caliper Plastic - Black	1 pc
08-040102	FatTrack Skinfold Fat Caliper Digital	1 pc

Measurement Tape MD < €









Reference	Description	Unit
08-080101	Baseline Measurement Tape with Gulick Attachment - 150cm	1 pc
08-080102	Baseline Measurement Tape with Gulick Attachment - 180cm	1 pc
08-080103	Baseline Measurement Tape with Hands-free Attachment - 150cm	1 pc
08-080201	Baseline Finger Circumference Gauge - 15cm	1 pc

Tuning Fork

These aluminium Tuning Forks are developed to measure sensitivity to vibration which can be an indicator of nerve damage.

It is also used for foot exam by diabetic people to check for signs of neuropathy.

Reference	Description	Unit
08-050301	Tuning Fork 128cps	1 pc
08-050302	Tuning Fork 256cps	1 pc
	08-050301	Reference Description 08-050301 Tuning Fork 128cps 08-050302 Tuning Fork 256cps



Inclinometer

MD (€

Baseline's Bubble and Digital Inclinometers are used to measure range-of-motion (ROM).

Some neck and back measurement protocols require the simultaneous use of 2 inclinometers.







- 10			
	-	-	
	- 10000	18	
	-	,, ex	
	- 3.50		The same of the sa

Reference	Description	Unit
08-060101	Baseline Bubble Inclinometer	1 pc
08-060201	Baseline Digital Inclinometer	1 pc

Scoliometer

Scoliometers are inclinometers designed to measure trunk asymmetry or Angle of Truncal Rotation (ATR) used in scoliosis screening programmes. They are small, non-invasive devices which are placed over the spine while a person is in a forward bending position. The scoliometer reading indicates degree of trunk asymmetry (rotation). Research shows that most people have a small degree of spinal rotation. However, depending on the measurement, spinal rotation is a possible indicator that scoliosis may be present.

Baseline Plastic Scoliometer



Pedi-Scoliometer

MD (E



Reference	Description	Unit
08-060102	Pedi-Scoliometer	1 pc
08-060104	Baseline Plastic Scoliometer	1 pc

Goniometers

Baseline MD (€




30cm Plastic Goniometer 08-030401



30cm Plastic Goniometer - HiRes™ 08-030402




Absolute + Axis® Attachment 08-030411



Reference	Description	Unit
08-030401	Baseline Plastic Absolute+Axis Goniometer - 30cm 360°	1 pc
08-030402	Baseline Plastic Absolute+Axis Goniometer - HiRes™ - 30cm 360°	1 pc
08-030411	Baseline Plastic Absolute+Axis - Attachment Only for 30cm Goniometer	1 pc
08-030601	Baseline Plastic Finger Goniometer - 9cm Flexion to Hyper-Extension	1 pc


Lafayette ∈ €



Reference	Description	Unit
08-030701	Lafayette Gollehon Extendable Goniometer with Magnifying Lens - 24cm - 66cm 180°	1 pc

Goniometer Kit ←





The Plastic set includes:

- Plastic Goniometer 15 cm | 0° to 360° per 1°
- Plastic Goniometer 20 cm | 0° to 360° per 1°
- Plastic Goniometer 30 cm | 0° to 360° per 1°
- Rulong Goniometer 18 cm | 0° to 360° per 2°
- Pocket Goniometer 15 cm | 180° x2
- Plastic Flexion-Hyperextension goniometer



Quality Soft Pouch Case

The Stainless Steel set includes:

- Robinson pocket goniometer 15cm 180°
- Stainless steel goniometer 18cm 180°
- Stainless steel goniometer 35cm 180° x2
- Stainless steel goniometer 35cm 360°
- Stainless steel finger goniometer 15cm
- X-ray goniometer 20cm 180°

Reference	Description	Unit
08-030211	Saehan Goniometer Set Stainless Steel 6 pieces	Set
08-030121	Saehan Goniometer Set Plastic 6 pieces	Set

Plastic Goniometer < €

Clear plastic permits observation of joint's axis of motion, and its range of motion. The 360° head has three scales calibrated to be used with ISOM (International Standard Of Measurements, STFR) system.













Reference	Description	Unit
08-030101	Saehan Plastic Goniometer 15cm 0° to 360° per 1°	1 pc
08-030102	Saehan Plastic Goniometer 20cm 0° to 360° per 1°	1 pc
08-030103	Saehan Plastic Goniometer 30cm 0° to 360° per 1°	1 pc
08-030111	Saehan Rulong Goniometer 18cm 0° to 360° per 2°	1 pc
08-030112	Saehan Pocket Goniometer 15cm 180° x2 per 5°	1 pc

Stainless Steel Goniometer < €





20cm / 1° increments 2 opposing 180° scales 08-030202

35cm / 2 opposing 180° scales / 1° increments Arm tension is controlled by thumb knob 08-030203



Reference	Description	Unit
08-030201	Saehan Stainless Steel Robinson Pocket Goniometer 15cm 180°	1 pc
08-030202	Saehan Stainless Steel Goniometer 20cm 180° x 2	1 pc
08-030203	Saehan Stainless Steel Goniometer 35cm 180° x 2	1 pc
08-030204	Saehan Stainless Steel Goniometer 35cm 360°	1 pc

Finger Goniometer < €

Finger and Flexion/Hyper-Extension Goniometers measure range-of-motion of finger joints. Head has two opposing 180° scales marked in 5° increments. A linear scale in centimetres is on arm.



SAEHAN



70

Monofilament



Monofilament evaluators are used to measure cutaneous sensory perception threshold.

08-050401

08-050421

08-050425

Reference	Description	Unit	
08-050401	Baseline Tactile Monofilament - LEAP Program - Disposable - 5.07 - 10 g		
08-050421	08-050421 Baseline Tactile Monofilament - 2-piece Set - Normal - 2.83, 3.61 g		
08-050425	Baseline Tactile Monofilament - 2-piece Set - Protective Sensation - 4.56, 5.07 g	Set	

Discriminator

MD (€









Reference	Description	Unit
08-050203	Touch-Test 2-point Discriminator	1 pc
08-050204	Baseline 2-point Discriminator	1 pc
08-050205	Baseline 3-point Discriminator	1 pc
08-050206	Baseline 2-point Discrim-A-Gon, 2 Disk Set	1 pc

Neurological Reflex Instruments MD < €

08-050102: Babinski Hammer

Pointed bottom can be used as





08-050103: Taylor Hammer

08-050101: Buck Hammer

Comes with a needle and a brush.

Use to elicit responses during examination

Two heads for cutaneous sensation evaluation.



08-050104: Dejerine Hammer This hammer is used for testing deep tendon reflexes

in a diagnostic evaluation. It features a double-sided head composed of two rounded rubber surfaces for increased patient comfort.



08-050201: Wartenberg Pinwheel

Elicit cutaneous sensory and pain perception responses.



Reference	Description	Unit
08-050101	MoVeS Buck Hammer	1 pc
08-050102	MoVeS Babinski Hammer	1 pc
08-050103	MoVeS Taylor Hammer	1 pc
08-050104	MoVeS Dejerine Hammer	1 pc
08-050105	MoVeS Troemner Hammer	1 pc

MoVeS Wartenberg Pinwheel

08-050105: Troemner Hammer

This hammer is shaped like a two-headed mallet. The larger mallet is used to elicit tendon stretch reflexes, and the smaller mallet is used to elicit percussion myo-



Sit & Reach Unit


MD (€



Sit & Reach Unit

08-060103

When using the Sit & Reach unit you sit with your feet against the "sit and reach box" and flex forward moving the slide bar. Flexion is read directly from the calibrated bar.



72



TENS device

MD (E



The Neuropex TENS is accepted to be far more accurate than a standard TENS machine.

- Highly user friendly with 3 programs and modes:
 - B (burst)
 - C (continuous)
 - M (Pulse Rate Modulation)
- · Regulated pulse width control for adjustable intensity and
- Clear plastic window allows controls to be viewed during operation
- Easy to use and very compact
- Adjustable intensity and stimulation strength
- Protective case for easy electrode storage and transportation
- Complete set of electrodes, in-built timer, complimentary spare







Neuropex Long Lasting Self-Adhesive Electrodes

- Resealable pouch for safe storage
- 4 Electrodes per pack
- 50mm x 50mm
- Ideal for use with all Neuropex TENS and Interferential devices

Reference	Description	Unit
08-070401	2-Channel Professional TENS unit	1 pc
08-070411	Replacement TENS Electrode - 4 piece set	set



We have nice products, but you would like your own brand or your customer's name on the band, loop, superloop...? We'll do it! Together with our manufacturers from all the products in our range, we can create your desired product with the look you want. We can replace our logo with yours, change the colour, even the packaging can be tweaked to your needs. Private Labels will allow you to design the product the way you want to match it to your company's or customer's identity.

PRIVATE LABEL ///



are putting your brand in the spotlights. This is perfect if you need a homogeneous range of products or want to increase your brand awareness. Or maybe you're looking for a one time opportunity: a promotional product during your company's special occasion, a gift for visitors or an event you're supporting?

Please contact us for more information about pricing and volume.

Product description	Minimum Order Quantity
Resistive Exercise Band	5.000m per colour
Resistive Exercise Loop	500 pcs per colour
Soft-Over-Ball	3.000 pcs per colour
Hollow Foam Roller	300 pcs
AB Gym Ball	500 pcs
NBR Comfort Gym Mat NBR Xtra Comfort Gym Mat	300 pcs per colour
Pearl Packs	500 pcs
Plastic Goniometer	2.000 pcs
Various other products	Contact us for more info



MVS IN MOTION CALIBRATION SERVICE

Measurement Instruments are calibrated at manufacturing with a full one year guarantee, before sending them to users all over the world. Ideally, these instruments should be recalibrated every year.

In order to make this happen in an easy and inexpensive way, Hydraulic Hand Dynamometers (Saehan, Jamar, Baseline...) and Hydraulic Pinch Gauges are recalibrated by our experts at the MVS In Motion Centre in Tisselt - Belgium. This operation consists of re-establishing the exact relationship between the force applied to a dynamometer and the measurement it shows. Small repairs can be done as well.



Simply send your request to our customer service department for a return authorisation. After each recalibration, you will receive an updated calibration certificate, certifying accurate and repeatable test results.

If you have many calibrations to do, you might consider buying your own calibration unit.



Reference	Description	Unit
08-019901	Calibration Unit	1 pc
08-019900	Calibration & Repair	-

MVS IN MOTION REPAIR CENTRE



Since machines are designed and assembled by humans, they are not flawless. That is why MVS In Motion installed an inhouse repair centre, to swiftly and correctly solve any problems that may occur with our products.

No matter what the cause, be it faulty fabrication or misuse, our technician will shine his light on the problem.

Send an email to our customer service department to inform us on the type of problem and the serial number/production date of the device. After our authorisation, please send over the patient, accompanied by your name and client number and a brief description of the symptoms to info@mvs-in-motion.com. We'll do our best to have him cured in no time.

3-in-1 Foam Roller	43	Flex-Ion	25	Ring Tubes	20
Accessories	24	Foam Hand Grip	32	Roll-On Packs	60
Adjustable Hand Grip	32	Foot Pump	48	Saehan Cycle	34
Aerobic Step	48	Goniometers	70-71	Safety Handles	24
Ankleciser	21	Gym Ball	51	Scoliometer	69
Babinski Hammer	73	Gym Mat	47	Shoulder Rope Pulleys	22,23
Back Wrap	59	Hand Pump	48	Shoulder Tube Pulley	22
Back-Leg-Chest Dynamometer	67	Handmaster Plus	27	Sit & Reach Unit	73
Balance Boards	50,51	Hollow Foam Roller	43	Sitting Cushions	44
Balance Pads	53	Hot/Cold Packs	61	Smedley Hand Dynamometer	66
Balance Trainers	51,52	Hydraulic Hand Dynamometer	65	Soft Handles	24
Ball Accessories	48	Hydraulic Hand Evaluation Kits	64	SoftMeds	55
Ball Stabiliser	48	Incline Boards	39	Soft-Over-Ball	45
Ball Stacker	48	Inclinometer	69	Spiky Foam Roller	43
Band	8-11,17	Kettlebells	57	Sport Pack + Strap	59
Band & Tubing Scissors	24	Kinesiology Tape	41	Squeeze Balls & Eggs	31
Band Latex-Free	13	Knee Rope Pulley	22	Squeeze Dynamometer	67
Bar	20	Knee/Elbow Wrap	59	Stepping Stones	42
Blood Pressure Monitor	63	Loop	15,19	Superloop	19
Body Fat Scale	66	Massage Ball, Roll & Ring	46	Swim Noodle	54
Body Twister	52	Massage Stick	43	Taylor Hammer	73
Buck Hammer	73	Massage Tools	46	TENS Device	74
Comfort Putty	33	Measurement Tape	68	Theraflex Putty	33
Constant Force X-Tend	27	Medicine Ball	57	Thermometers	62
Cuff-Ring Tube Loop	20	Monofilament	72	Thumb Lock Wrist Weight	56
Curamotion	38	Multi-Grip Exerciser	21	Torex Pro Ankle Boot	60
Dejerine Hammer	73	Needle Pump	48	Triggerpoint Rollers	43
Digi-Extend	26	Neurological Instruments	73	Troemner Hammer	73
Digi-Flex	26	O-Ring Tube Loop	20	Tube	14,18
Digital Hand Dynamometers	65	Oximeters	63	Universal Assist	24
Discriminators	72	Oxycycles	39-37	Varigrip Sport	28
Door Anchor	23	Paraffin Bath	58	Varigrip Therapy	28
Dumbbells	55	Peanut Ball	51	Varigrip Uno	28
Dynadome	49	PearlPacks	59	Wall Mount	24
Ekilibres	52	Pedal Exercisers	39	Wartenberg Pinwheel	73
Electric Compressor Pump	48	Pilates Foam Roller	45	Wheeleo	35
Ergonomic Foam Roller	43	Pilates Ring	45	Wide Loop	15
Eye Mask	59	Pinch Gauge	67	Work-Out Station	52
Fabric Elastic Band	49	Power-Webs	29-30	Wrist & Ankle Weight	56
Faster Blaster Pump	48	Prohands Series	26	Wrist Exerciser	31
Fat Caliper	68	Pumps	48	Yoga Block	42
Figure 8-Ring Tube Loop	20	Push-Pull Dynamometer	66	Yoga Mat	42
Fit Stretch	34	Putties	33	Yoga Strap	42

76

TERMS & CONDITIONS

- 1. Terms: Prepayment, unless otherwise specified. Buyer should bear all bank charges.
- Minimum purchase value: An administration fee of €10,00 will be added on all orders below the minimum purchase value of €250,00.
- 3. Freight Terms: Prices are ex-works Tisselt, Belgium. All shipping/delivery costs are at charge of the buyer, unless otherwise specified. Visible damage to goods delivered, not duly noted on the consignment note of the carrier, cannot be compensated in any form.
- 4. Prices: All prices for merchandise purchased are based on seller's price list in effect on the time of the order. Prices can be subject to change without notice.
- 5. Claims: Unless otherwise provided for herein, claims of any kind in goods sold by seller shall be made by buyer to seller, in writing, within a period of seven days after such goods are delivered to buyer.
- 6. Credit and Exchange: Written authorisation must be obtained before returning merchandise for credit or exchange! Prior to return, buyer shall furnish seller an itemized inventory of the products to be returned, and any other information seller may request, such as product lot number, manufacturing date, expiration date, ... In no event shall seller accept products for return more than thirty (30) days after delivery by seller to buyer and products which are damaged, used or not in their original and unopened packaging. Seller will charge a handling fee of ten percent (10%) of the returned goods.
- 7. Finance Charge: A finance charge of 1% per month (12 % per annum) will be added to all amounts past due. A penalty of 15% on value, with a minimum of €40,00, for late payment will be added. Terms of payment are from date seller's invoice and not contingent upon delivery.
- 8. Governing Law: All disputes arising out of all sales by seller to buyer shall be governed by the laws of Belgium and place of jurisdiction is Antwerp, section Mechelen.
- 9. Warranties: Seller warrants that it has good title to the goods purchased by buyer and that its goods will conform to buyer's specifications, if any, except as otherwise provided herein. All other warranties, whether express or implied, and including, but not limited to, the warranty of merchantability and fitness for a particular purpose are excluded. In all events, buyer's sole and exclusive remedy shall be limited to the amount of the purchase price of any product on which a claim is made as its exclusive remedy for breach of the warranty expressed herein and for negligence and otherwise. Seller reserves the right to repair or replace warranted products in its descretion. Under no circumstances shall seller be liable to buyer or any other person for consequential or incidental damages.
- 10. Property: Seller has the exclusive property on the delivered products until buyer has paid the full amount due. This also implies that if buyer does not pay in time, seller has the right to, without prenotice, claim back the delivered goods.







The choice of physiotherapists for 50 years, with know-how and expertise at your service.

www.ferrox.it









Master Distributor For























VISIT OUR OTHER WEBSITES FOR

- PRODUCT INFORMATION
- EXERCISE CHARTS AND VIDEOS
- SHOP LOCATOR

WWW.MOVES-YOU.COM WWW.MOVES-BAND.COM WWW.MAMBOMAX.COM