

ALMAG-02

A MULTI-FUNCTION PULSED MAGNETIC THERAPY DEVICE FOR TREATING CHRONIC DISEASES



EFFECTIVE

- 2 kinds of magnetic fields: travelling & static
- 3 types of emitters for local & zonal magnetic therapy
- 79 preset exposure programs for various diseases



WHAT DOCTORS GET

- professional device
- possibility of treating people with various diseases
- wider range of therapeutic services in a clinic

THE MAIN EFFECTS:

- improves circulation and metabolism
- accelerates repair processes
- relieves pain
- stops inflammation
- stimulates immunity
- reduces dependency on drugs

APPLICATIONS:

- Locomotor system diseases
- Diseases of circulatory system
- Traumas and injuries
- Nervous system diseases



EASY TO USE

- compact & portable
- intuitive user interface
- step-by-step demonstrative instruction



Licensed as a medical device
Clinically tested

TECHNICAL FEATURES

Magnetic induction	2-45 mT
PEMF frequency	1-100 Hz
Weight, max	11kg
Power supply	220V/50Hz
Power consumption	50 VA



ALMAG-02

PLACEMENT VARIABILITY

