



MoVeS Wrist Exerciser



Adjustable resistance



The MoVeS Wrist Exerciser is an adjustable wrist and forearm trainer which allows the user to exercise with progressive resistance. Its ergonomically designed forearm support keeps the wrist in the right place. Together with the soft handgrip, adjustable to 4 positions, it ensures maximum comfort while rehabilitating your wrist and forearm.

Reference	Description	Unit	CE
02-070201	MoVeS Wrist Exerciser	1 pc	CE




Adjustable Hand Grip



Adjustable resistance

The Mambo Max Adjustable Hand Grip comes in two different versions. Whether it is the 5 to 20 kg light version (orange) or the 10 to 40 kg heavy version (grey), you can always adjust the required tension simply by turning the control knob.



Reference	Description	Unit	CE
02-070102	Mambo Max Adjustable Hand Grip - Light	1 pc	CE
02-070104	Mambo Max Adjustable Hand Grip - Heavy	1 pc	CE



Foam Hand Grip



The Mambo Max Foam Hand Grip assures a firm training grip. The foam handles provide a firm grip for the entire workout session. Take it with you to the office or just use it at home. Wherever you want to use it, the Mambo Max Foam Hand Grip assures a tight grasp.

Gripping power = 20kg



Reference	Description	Unit	CE
02-080101	Mambo Max Foam Hand Grip	1 pair	CE

Wheeleo



Reduced mobility presents 3 challenges:

- Keep the person moving;
- Reduce the risk of falling;
- Loss of autonomy.

Wheeleo, the rolling cane, is designed to offer a solution to these 3 challenges and offers additional benefits:

- Ease of handling, no special training required;
- **Walking** goes a lot **faster** than with a quadripod cane;
- **Safer walking** than with a quadripod cane, as there is constant ground contact;
- **One-handed use** while still maintaining constant ground contact, the other hand can be used to move an object, to open doors, ...
- **Minimal space required** and perfect for indoor use where a rollator is too big.
- **Reinforced stability:** the wheels offer a security system. A lateral push will result in a lateral movement, not a tilting/fall.



Faster and safer walking



Wheeleo is a unique quadripod rolling cane used to improve mobility. It is designed with the aim of offering help without stigmatising the disability of its user.

It will make walking trips smoother and safer which will allow the user to maintain or **improve autonomy**.

The Wheeleo is used with one hand. Thanks to its 4 wheels, you do not have to lift it, which improves stability and **reduces the risk of falling**.

The Wheeleo is super easy to use. Ultra-intuitive, its user will always position it in the right place to guarantee its stability and balance.

An important feature of the Wheeleo is that it **does not require brakes** to prevent its user from falling. Thanks to its rolling design, the Wheeleo will always be placed in the right position for its user.

When losing balance, the user will automatically push the Wheeleo in the right position that is comfortable for him, without lifting the Wheeleo.

Adding brakes would be the same as replacing the Wheeleo with a regular cane because the brakes would block the wheels. As a consequence the Wheeleo would tilt and its user would possibly fall.



Adjustable height

Rolling Cane Wheeleo

↕ 77,5 - 101,5cm
↔ 38 x 42cm
👤 120kg

Reference	Description	Unit	CE
03-060101	Wheeleo NEW	1 pc	CE

